



NeuMask

Ear-loop Procedure Mask

Medgluv Procedure Earloop Face Masks are a level 2 moderate barrier mask appropriate for all general patient care situations and procedures involving low-to-moderate levels of spray, spatter, moisture and airbourne particles.

Soft and strong, ample-length earloops so there is less pulling, and adjustable nosepiece for comfortable fit. Especially useful for laser procedures that generate harmful airborne particles. Comfortable and highly breathable in any settings. Light, comfortable masks for extended wear when necessary. Latex free, ear loops are made with no irritating rubber or plastic touches the skin. Three-ply, glass-free filter media is surrounded by soft, non-irritating inner and outer facing.

Description

- Procedure masks with elastic ear-loops provide complete protection
- Polypropylene material offer a soft, fluid protection barrier.
- A full-width nosepiece guarantees proper fit.
- High Filtration Melt-blown Filter. Easy breathing; Delta-P < 6.0
- Bacterial filtration Efficiency (BFE) > 98% @ 3.0µm.
- Meets and exceeds ASTM F2100-19 Level 2 Standards

Item No: MGM103 (Level 2)
Color: Pink
Packaging: 50 masks/box, 10 boxes/case

