



# Insulin Syringes

## Features

- Latex free, Non-toxic, Non-pyrogenic
- Nearly painless injection with comfort
- Super-secure attachment of needle
- Reduced dead space for less waste
- Bold, clear and large numbers for easy readability
- E.O Gas sterilization
- 5-year shelf life



**We use Thin Wall Needle Tube**  
Wall Thickness: 30G - TW | 31G - ETW | 32G - UTW

Contain	0.3cc	0.5cc	1.0cc
3/16" 5mm	31G	31G	31G
5/16" 8mm	● 30G	● 30G	● 30G
	● 31G	31G	● 31G
	32G	32G	32G
1/2" 13mm	● 29G	29G	29G
	30G	30G	30G



Packaging: 100 syringes/box (blister packaged)