

Air-Tite Aesthetics


THE INVISIBLE NEEDLE™ and other Low Dead Space (LDS) Needles by TSK
available from Air-Tite Aesthetics!

THE INVISIBLE NEEDLE™

YOUR PATIENTS WILL NEVER SEE IT COMING.

NEAR PAINLESS
REDUCED BLEEDING
PRODUCT LOSS CLOSE TO ZERO


LEARN MORE:



THE INVISIBLE NEEDLE™ is TSK's thinnest needle to date, with an outside diameter of just 0.2mm. It is 33% thinner than the popular TSK 30G needles, which may result in a more comfortable injection experience for the patient with possibly less pain, less swelling, and less bruising. THE INVISIBLE NEEDLE™ features the new Low Dead Space needle hub.

The Low Dead Space (LDS) needle's design reduces the excess effective hub space to nearly ZERO, which may result in a potential cost and product savings as it leaves behind so little of the unused precious injectables.

Part #:	Size:
TSK3013LDS	30G x 13mm
TSK3309LDS	33G x 9mm
TSK3313LDS	33G x 13mm
TSK02009LDS (THE INVISIBLE NEEDLE™)	.2mm x 9mm



All Low Dead Space Syringes Are Recommended To Be Used With Luer Lock Low Dead Space Syringes

Air-Tite-Shop.com | Air-TiteAesthetics.com | 800.231.7762