




FRESHLY HATCHED: OUR NEW PNEUMATIC VACCINATOR FOR DAY-OLD CHICKS

Efficacy	adjustable trigger element for safe and consistent positioning of day-old chicks
Accurate	highest dosing accuracy due to reliable double-stroke cylinders and high-quality syringes
Sophisticated	perfectly shaped cover enables best injection angle for healthy day-old chicks
User-friendly	easy change of injection volume with only one syringe due to a patented click system
Smart	easy set-up, maintenance and cleaning procedures for long service life and highest biosecurity standards
Versatile	injection of either one or two different vaccines with only one needle
Resistant	robust and medical grade components allow use with nearly every vaccine

Item No.:

8300048263 – HSW Pneumatic Vaccinator Single

8300048264 – HSW Pneumatic Vaccinator Double



For further information please contact:

Henke-Sass, Wolf GmbH

Keltenstraße 1

78532 Tuttlingen – Germany

Tel. +49 (0)7462 9466 380

E-Mail: michael.speck@henkesasswolf.de

www.henkesasswolf.de