

SET Medikal INJECT

Sterile 3-Part Syringes



COMPONENTS:

- Barrel: PP - Polypropylene
- Plunger: PE - Polyethylene, PP - Polypropylene
- Plunger Gasket: Isoprene Rubber, Latex Free
- Lubricant: Silicone oil
- Sterilization Method: Ethylene Oxide

ATTRIBUTES:

- Made in Turkey
- Human medical use / 510K Registered
- Plastic raw material of glass transparency and medical purity
- Latex-free seal
- Easy and comfortable piston movement
- High impermeability
- Stop point design preventing the easy exit of the piston from the cylinder
- Design allowing single-handed use
- Medical grade blister packing consisting of transparent film

AIR-TITE PART #	SET INJECT SYRINGES	SYRINGES PER BOX	BOXES PER CASE
RL3SET	3mL LUER LOCK	100	24
RL5SET	5mL LUER LOCK	100	26
RL10SET	10mL LUER LOCK	100	12
RL50SET	50/60mL LUER LOCK	60	4

