



TSK MICROCANNULA
STERiGLIDE™



LESS PAIN
LESS BRUISING
LESS SWELLING

MORE PRECISE

50%
Better Gliding

STERiGLIDE™
The cannula that glides.

Air-Tite
Aesthetics



AIR-TITE AESTHETICS PROUDLY DISTRIBUTES
THE AWARD WINNING TSK MICROCANNULA STERiGLIDE™
RECIPIENT OF THE ANTI-AGING AND BEAUTY TROPHY
FOR BEST INJECTION TOOL

WE SUPPLY SPECIALTY HYPODERMIC NEEDLES,
SYRINGES + ACCESSORIES FOR MANY
AESTHETIC PROCEDURES

AIR-TITE AESTHETICS.COM | 800.231.7762



THE 4

TSK MICROCANNULA STERiGLIDE™




KEY REASONS TO USE TSK MICROCANNULA STERiGLIDE™

- 1 LESS PAIN**
Proprietary surface treatment reduces friction within the dermis
- 2 LESS BRUISING**
Patented dome-shaped tip reduces penetration forces by 50%
- 3 LESS SWELLING**
Less edema and trauma due to reduced tissue subcision
- 4 MORE PRECISE**
Clearly visible hub marking for extrusion port location



**“THESE CANNULAS HAVE
REVOLUTIONIZED MY INJECTABLE
PRACTICE AS AN OCULOFACIAL
PLASTIC SURGEON.”**

—Dr. Lopa Gupta
www.drlopagupta.com



THE #1 CHOICE

DOME-SHAPED TIP DESIGN

- Easier cannula introduction
- Reduced risk of bruising

The patented dome-shaped tip of the **STERiGLIDE™** offers less resistance than a traditional cannula tip. This reduces the cannula introduction force needed and the risk of bruising.

NEAREST TO TIP FILLER DELIVERY

- More accurate filler placement
- Reduced product loss

Practitioners can see and feel the end of the cannula tip. Delivery of filler near the end of the tip creates a more accurate sense of filler placement. It also reduces product loss at the distal end of the tip.

PROPRIETARY SURFACE TREATMENT

- Up to 50% better gliding
- Increased patient comfort

STERiGLIDE™ cannulas are made of the highest quality steel with a proprietary surface treatment coating. This reduces friction within the soft tissue, enhancing patient comfort. The surface treatment also improves the handling and maneuverability of the cannula.

CLEAR SIDE-PORT MARKING

- More accurate filler placement
- Better handling

The hard polymer hub of the **STERiGLIDE™** cannula has a clear side-port marking that indicates the position of the cannula opening while serving as a visual guide for more accurate filler delivery.


NEEDLE VS. MICROCANNULA

Originally, the only method for delivering dermal fillers was the sharp-tipped hypodermic needle.


Advancements in the evolution of specialty blunt-tipped microcannulas have changed this procedure forever.

Clinicians have embraced **STERIGLIDE's** dome-shaped tip design over its first generation counterparts.

CANNULA FLEXIBILITY + CONTROL



Filler injection with sharp needle



Filler injection with STERIGLIDE™ microcannula

Sharp-tipped needles are rigid. Microcannulas are a more flexible option that move freely when delivering fillers.

50% Better Gliding

STERIGLIDE™

The cannula that glides.



**“FINALLY A CANNULA THAT
ALLOWS FOR ACCURATE
AND CONTROLLED
GUIDANCE.”**

—Dr. Dalvi Humzah

www.pdsurgery.co.uk

GOLD Standard

Outperforms any other
microcannula available



“THE MOST PRECISE AND
EASIEST TO HANDLE
CANNULA I’VE EVER USED.”

—Dr. Sabine Zenker
www.dr-zenker.de



Anti Aging & Beauty




Expert 2 Expert

TSK introduced its award winning **STERiGLIDE™** cannula in 2012. The **STERiGLIDE™** brand is recognized as the gold standard around the world.



ORDER INFO



**TSK MICROCANNULA
STERiGLIDE™**

CODE #	AWARD WINNING BEST INJECTOR	DIMENSIONS
TSK 2250SG	STERiGLIDE™ Aesthetic Microcannulas	22g x 50mm (2.0")
TSK 2270SG	STERiGLIDE™ Aesthetic Microcannulas	22g x 70mm (2¾")
TSK 2538SG	STERiGLIDE™ Aesthetic Microcannulas	25g x 38mm (1½")
TSK 2550SG	STERiGLIDE™ Aesthetic Microcannulas	25g x 50mm (2.0")
TSK 2725SG	STERiGLIDE™ Aesthetic Microcannulas	27g x 25mm (1.0")
TSK 2738SG	STERiGLIDE™ Aesthetic Microcannulas	27g x 38mm (1½")
TSK 2750SG	STERiGLIDE™ Aesthetic Microcannulas	27g x 50mm (2.0")
TSK 3025SG	STERiGLIDE™ Aesthetic Microcannulas	30g x 25mm (1.0")

Each box contains 20 microcannulas and pilot needles.



TSK 33G x ½" NEEDLE

Newly designed with a stiffer shaft and sharper point, this specialty needle improves patient comfort. # TSK3313



TSK NEEDLE HOLDER

Safely secures the hard-shell plastic capsule of cannulas and needles for safer loading and removal from the syringe. # TSKNH1



HSW SOFT-JECT® 1ml LUER LOCK

Low Dead Space tip to reduce waste of expensive injectables and risk of needle pop-off. # RL1



HSW SOFT-JECT® 1ml LUER SLIP

Low Dead Space syringe tip reduces loss of expensive injectables. # R1



HSW NORM-JECT® 1ml LUER SLIP

Syringes are silicone free with Low Dead Space tip and luer slip fitting. # A1



AWARD WINNING
MICROCANNULA FOR
BEST INJECTION TOOL

STERIGLIDE™
The cannula that glides.

To request samples or place an order visit:

Air-TiteAesthetics.com

Info@Air-TiteAesthetics.com

1.800.231.7762

Air-Tite AESTHETICS

AIR-TITE PRODUCTS CO., INC., 565 CENTRAL DRIVE, VIRGINIA BEACH, VA 23454

T: 800.231.7762 L: 757.340.2501 F: 757.340.2912

