

2-Part Syringe Differences


HENKE-JECT

Air-Tite Part #: A50



PRO-JECT

Air-Tite Part #: MA160



2-Part Syringe Differences

HENKE-JECT

Air-Tite Part #: A50

- FDA Registered
- CE Marked
- Eccentric Tip
- Plunger Safety Stop
- Blister Pack Packaging
- 2mL Graduation Markings
- Numerical Graduations Every 10mL
- Sterilized By Ethylene Oxide
- Ounces Not Indicated
- Grip Wings
- 10 Boxes Per Case
- 50mL (*With Further Graduations To 60mL)
- Less Stiff Operation
- Longer And Thinner Design Dimensions
- Semi-Transparent White Plunger
- Clearer Barrel
- Blue Ink
- Made In Turkey
- Expiration Date Indicated (5 Years)

PRO-JECT

Air-Tite Part #: MAI60

- Not FDA Registered - Lab/Vet Use Only
- Not CE Marked
- Centric Tip
- No Plunger Safety Stop
- Poly Bag Packaging
- 1cc Graduation Markings
- Numerical Graduations Every 5cc
- Sterilized By Gamma Irradiation
- Ounces Indicated
- Flat Wings
- 6 Boxes Per Case
- 60cc
- Stiffer Operation
- Shorter And Wider Design Dimensions
- Completely Opaque White Plunger
- Slightly Cloudy Barrel
- Black Ink
- Made In USA
- Expiration Date Not Noted