

2-PART SYRINGE COMPARISON

HENKE-JECT®

Air-Tite Part #: A50

DISCONTINUED



- Manufactured by Henke-Sass Wolf
- No Rubber, No Silicone Oil
 - Eccentric Tip
 - Plunger Safety Stop
 - Blister Pack Packaging
- Sterilized by Ethylene Oxide
- 2mL Graduation Markings
- Numerical Graduations Every 10mL
 - Grip Wings
 - 30 Pieces/Box
 - 10 Boxes Per Case
- 50mL (*With Further Graduations To 60mL)
 - Polypropylene (PP) Barrel
- High-Density Polyethylene (HDPE) Plunger
 - Semi-Transparent White Plunger
 - FDA Registered
 - CE Marked
- 25lbs 32"x15"x15" Case Dimensions
 - Blue Ink
 - Made In Turkey
- Expiration Indicated (5 Years)
- Passed MEA Testing

AIR-TITE PREMIUM

Air-Tite Part #: AMS50



- Manufactured for Air-Tite Products
- No Rubber, No Silicone Oil
 - Eccentric Tip
 - Plunger Safety Stop
 - Blister Pack Packaging
- Sterilized by Ethylene Oxide
- 1mL Graduation Markings
- Numerical Graduations Every 10mL
 - Grip Wings
 - 30 Pieces/Box
 - 12 Boxes Per Case
- 50mL (*With Further Graduations To 60mL)
 - Polypropylene (PP) Barrel
- High-Density Polyethylene (HDPE) Plunger
 - Semi-Transparent White Plunger
 - Not FDA Registered - Lab/Vet Use Only
 - Not CE Marked
- 34lbs 23"x18"x19" Case Dimensions
 - Blue Ink
 - Made In China
- Expiration Indicated (5 Years)
- Passed MEA Testing