

# 2-PART SYRINGE COMPARISON

## HENKE-JECT®

Air-Tite Part #: AL50

**DISCONTINUED**



- Manufactured by Henke-Sass Wolf
- No Rubber, No Silicone Oil
- Centric Tip
- Plunger Safety Stop
- Blister Pack Packaging
- Sterilized by Ethylene Oxide
- 2mL Graduation Markings
- Numerical Graduations Every 10mL
- Grip Wings
- 30 Pieces/Box
- 10 Boxes Per Case
- 50mL (\*With Further Graduations To 60mL)
- Polypropylene (PP) Barrel
- High-Density Polyethylene (HDPE) Plunger
- Semi-Transparent White Plunger
- FDA Registered
- CE Marked
- 25lbs 32"x15"x15" Case Dimensions
- Blue Ink
- Made In Turkey
- Expiration Indicated (5 Years)
- Passed MEA Testing

## AIR-TITE PREMIUM

Air-Tite Part #: AML50



- Manufactured for Air-Tite Products
- No Rubber, No Silicone Oil
- Centric Tip
- Plunger Safety Stop
- Blister Pack Packaging
- Sterilized by Ethylene Oxide
- 1mL Graduation Markings
- Numerical Graduations Every 10mL
- Grip Wings
- 30 Pieces/Box
- 12 Boxes Per Case
- 50mL (\*With Further Graduations To 60mL)
- Polypropylene (PP) Barrel
- High-Density Polyethylene (HDPE) Plunger
- Semi-Transparent White Plunger
- Not FDA Registered - Lab/Vet Use Only
- Not CE Marked
- 34lbs 23"x18"x19" Case Dimensions
- Blue Ink
- Made In China
- Expiration Indicated (5 Years)
- Passed MEA Testing