



2024 HEALTHCARE AND MEDICAL CATALOG



YOUR SYRINGE AND NEEDLE SPECIALISTS SINCE 1926

Air-Tite Products Co., Inc.

565 Central Drive, Virginia Beach, VA 23454
T: 800.231.7762 F: 757.340.2912 L: 757.340.2501 | Air-Tite.com



WE SUPPLY SPECIALTY HYPODERMIC NEEDLES,
SYRINGES, & ACCESSORIES FOR MANY
AESTHETIC PROCEDURES


AWARD WINNING
MICROCANNULA FOR
BEST INJECTION TOOL

STERiGLIDE™

The cannula that glides.



Air-Tite
Aesthetics



Request samples or place orders at:
www.Air-Tite.com
atinfo@air-tite.com
800.231.7762

See page 32 for more information



**ANTI - AGING &
BEAUTY AWARD**

Each box of **STERiGLIDE™**
is packaged with 20
microcannulas and pilot needles
for entry hole preparation.

SCAN FOR
PRICE LIST



TABLE OF CONTENTS

SYRINGES

2-Part Syringes	1
Henke-Ject Syringes	3
Exel Syringes	5
BD Syringes	7
MinimLL Syringes	8
Omnifix Syringes	9
Oral Syringes	10

SYRINGES WITH NEEDLES

Henke-Ject Syringes with Needles	11
EXEL Syringes with Needles	13
BD Syringes with Needles	15

NEEDLES

EZ FLO Needles	17
TSK Needles	19
Exel International Needles	21
HSW Needles	23
JBP Needles	25
BD Needles	27
Exel Spinal Needles	29
BD Spinal Needles	30

MICROCANNULAS

TSK CSH Cannulas	31
TSK STERiJECT Cannulas	32

IV AND BLOOD COLLECTION

Intravenous Sets	33
Blood Collection	35

SUPPLIES

Sharps Collection, Sharps Mail-Back System	37
Gloves	39
General Supplies	41

2-PART SYRINGES

NORM-JECT®

B | BRAUN

ALL PLASTIC SYRINGES


Norm-Ject® syringes are latex free, contain no rubber or silicone oil. These syringes are the choice for any situation needing an inert, hypo-allergenic syringe. Because of their purity, applications by users have included IVF, embryo-transfer, chromatography and many laboratory procedures. Manufactured only from laboratory grade polypropylene and polyethylene, these unique plastic syringes have a positive safety stop to help prevent accidental spills. Norm-Ject® syringes are individually sterile strip packed and are also available in bulk, unsterile packaging. 1mL - 20mL syringes are Norm-Ject from B Braun Melsungen AG. FDA registered, CE Mark, ISO 13485 certified.



CONVENTIONAL RUBBER TIP SYRINGES

Although in general use for medical injections, rubber tipped disposable syringes have several drawbacks compared to Norm-Ject® all plastic syringes. Since they contain natural or synthetic rubber and silicone oil, they may contaminate any media exposed to the syringe. Adverse effects from rubber curing agents, which contain zinc, sulfur, and other elements, have been observed on everything from chromatographs to in-vitro procedures. Synthetic rubber plungers are also much more reactive to organic solvents and the silicone oil can cause separation of emulsions. In short, your best choice is the Norm-Ject® syringe if you are interested in purity and accuracy.








STERILE 2-PART SYRINGES

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
A1	1mL Tuberculin	Luer Slip	Sterile	100	1800
A3	3mL Centric	Luer Slip	Sterile	100	2500
A5	5mL Centric	Luer Slip	Sterile	100	2000
A10	10mL Eccentric	Luer Slip	Sterile	100	1200
A20	20mL Eccentric	Luer Slip	Sterile	100	800

** A1 Syringe has Low Dead Space Plunger Tip

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
AL3	3mL	Luer Lock	Sterile	100	2500
AL5	5mL	Luer Lock	Sterile	100	2000
AL10	10mL	Luer Lock	Sterile	100	1200
AL20	20mL	Luer Lock	Sterile	100	800



STERILE HENKE-JECT®

Conventional Syringes



STERILE HENKE-JECT®

Henke-Sass, Wolf GmbH Henke-Ject® (formerly Soft-Ject®) syringes are a high quality disposable rubber piston syringe offering substantial savings. A lubricated synthetic rubber tipped piston allows for smooth draw while a positive stop helps to prevent accidental pullout. These syringes are ETO gas sterilized and are latex and pyrogen free. All Henke-Ject® syringes are packaged in individual sterile strip packs. FDA registered, CE Mark, ISO 13485 certified.





STERILE HENKE-JECT®

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
R1**	1mL Tuberculin	Luer Slip	Sterile	100	2000
R2	2.5mL Centric	Luer Slip	Sterile	100	2500
R5	5mL Centric	Luer Slip	Sterile	100	2000
R10	10mL Eccentric	Luer Slip	Sterile	100	1200
R20	20mL Eccentric	Luer Slip	Sterile	100	800
R30HJ	30mL Eccentric	Luer Slip	Sterile	100	600

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
RL1 **	1mL	Luer Lock	Sterile	100	2000
RL3	3mL	Luer Lock	Sterile	100	2400
RL5	5mL	Luer Lock	Sterile	100	2000
RL10	10mL	Luer Lock	Sterile	100	1200
RL20	20mL	Luer Lock	Sterile	100	800
RL30HJ	30mL	Luer Lock	Sterile	100	600

**R1 and RL1 have low dead space plunger tip

EXEL INTERNATIONAL®

DISPOSABLE SYRINGES



EXEL is known for unsurpassed quality at a low price. These syringes are latex free, sterile, non-pyrogenic, and non-toxic. The design includes bold print, clear barrels, and siliconized gaskets for smooth, consistent plunger motion and a positive stop.

1mL syringes utilize a low dead space plunger

5mL syringes graduated to 6mL

10mL syringes graduated to 12mL

20mL syringes graduated to 25mL

30mL syringes graduated to 35mL

50mL syringes graduated to 60mL





EXEL INTERNATIONAL

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
ES1	1mL Tuberculin	Luer Slip	Sterile	100	1000
ES3	3mL Centric	Luer Slip	Sterile	100	1000
ES5	5mL Centric	Luer Slip	Sterile	100	800
ES10	10mL Centric	Luer Slip	Sterile	100	800
ES20	20mL Eccentric	Luer Slip	Sterile	50	500
ES30	30mL Eccentric	Luer Slip	Sterile	50	250
ES60	60mL Eccentric	Luer Slip	Sterile	25	150

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
EL1	1mL	Luer Lock	Sterile	100	1000
EL3	3mL	Luer Lock	Sterile	100	1000
EL5	5mL	Luer Lock	Sterile	100	800
EL10	10mL	Luer Lock	Sterile	100	800
EL20	20mL	Luer Lock	Sterile	50	500
EL30	30mL	Luer Lock	Sterile	50	250
EL60	60mL	Luer Lock	Sterile	25	150

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
EC30	30mL	Catheter Tip	Sterile	50	250
EC60	60mL	Catheter Tip	Sterile	25	150



BD[®] SYRINGES



Premium quality from the world's largest syringe manufacturer, Becton Dickinson. BD syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, and positive plunger stop. Available in both Luer-Lok and Luer Slip tips, BD syringes are made in America.

BD - BECTON DICKENSON

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
BD309659	1mL Tuberculin	Luer Slip	Sterile	200	1600
BD309656	3mL Centric	Luer Slip	Sterile	200	800
BD309647	5mL Centric	Luer Slip	Sterile	125	500
BD302831	20mL Centric	Luer Slip	Sterile	48	192
BD309654	50mL Centric	Luer Slip	Sterile	40	160



Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
BD309628	1mL Polycarbonate	Luer Lock	Sterile	100	800
BD309657	3mL	Luer Lock	Sterile	200	800
BD309646	5mL	Luer Lock	Sterile	125	500
BD302995	10mL	Luer Lock	Sterile	200	400
BD302830	20mL	Luer Lock	Sterile	48	192
BD302832	30mL	Luer Lock	Sterile	56	224
BD309653	50mL	Luer Lock	Sterile	40	160



Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
BD309620	50mL	Catheter Tip	Sterile	40	160

MINIMLL SYRINGES



THE ORIGINAL 0.5mL LOW DEAD SPACE LUER LOCK SYRINGE


- PRECISE DOSING
- MINIMAL PRODUCT LOSS
- EXCHANGE NEEDLES
- REDUCE NEEDLE POP-OFF
- SUPPORT BREAST CANCER FOUNDATIONS

The MinimLL syringe is the newest addition to the Air-Tite Products portfolio of specialty medical supplies. Created due to market feedback from clinicians who saw a need for very specific features, we are excited to introduce this first option of a less-than-1mL, Luer Lock syringe.

Features include:

- Easy to read 0.01 mL graduations for precise dosing of the smallest volumes
- Low Dead Space Plunger Tip allows for minimal residual volume
- Luer Lock Barrel Tip lowers the risk of needle pop-off
- Pink Plunger indicates a percentage of our sales will be donated to breast cancer research and awareness programs

MinimLL Syringes are manufactured with Polypropylene, which has a long history of use in various medical devices and is considered a safe vessel for most uses. The low dead space (LDS) feature minimizes loss of material, and can be combined with low dead space needles for additional product savings. In addition to helping reduce needle pop-off, the luer lock mechanism allows for the exchange of needle (such as using a large gauge needle for draw-up and a fine gauge for injection).



MINIMLL SYRINGES

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
MINI05	0.5mL Low Dead Space	Sterile	100	2000

Air-Tite's line of Pink products contribute a share of the sales to Breast Cancer Research and Awareness Programs. Other products include Pink Gloves and Pink Masks.

OMNIFIX® SYRINGES

B | BRAUN

GERMAN ENGINEERED SYRINGES

B. Braun Omnifix® syringes are high quality German engineered disposable rubber piston luer lock syringes. These syringes have a highly transparent polypropylene barrel, high contrast black graduations for easy reading, polystyrene plungers with an integrated backstop to prevent unintentional plunger withdrawal, and a rubber piston with double sealing ring made of polyisoprene allowing for slow aspiration and injection. These syringes are PVC, latex, and DEHP free. All Omnifix® syringes are packaged in individual sterile strip packs. FDA registered, ISO 9001 certified.

**Low Dead Space (LDS) Plunger allows for minimal residual volume.
Fine 0.01 ml graduation allows for precise dosage of small volumes.**



OMNIFIX® SYRINGES

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
BBL1	1mL Low Dead Space	Sterile	100	2400

ORAL DOSING SYRINGES

B | BRAUN

B BRAUN BULK ALL-PLASTIC ORAL SYRINGES

Catalog Number	Specifications	Tip Type	Sterility	Quantity
BBORAL1	1mL Low Dead Space	Oral Dosing	Unsterile	7000
ADB5SP	5mL	Oral Dosing	Unsterile	3000



COMAR ALL-PLASTIC ORAL SYRINGES

Catalog Number	Specifications	Tip Type	Sterility	Quantity	Packaging
CD1BX	1mL Low Dead Space	Oral Dosing Short	Unsterile	100	Bag
CD1	1mL Low Dead Space	Oral Dosing Short	Unsterile	2500	Case
CD3LBX	3mL	Oral Dosing Long	Unsterile	100	Bag
CD3L	3mL	Oral Dosing Long	Unsterile	1900	Case
CD3S	3mL	Oral Dosing Short	Unsterile	1500	Case
CD5S	5mL	Oral Dosing Short	Unsterile	1000	Case
CD5L-900	5mL	Oral Dosing Long	Unsterile	900	Case
CD10S	10mL	Oral Dosing Short	Unsterile	600	Case
CD10L	10mL	Oral Dosing Long	Unsterile	700	Case



BD ORAL SYRINGES

Catalog Number	Specifications	Tip Type	Sterility	Quantity
BD305218	5mL	Oral Dosing	Unsterile	500
BD305219	10mL	Oral Dosing	Unsterile	500




STERILE HENKE-JECT®

SYRINGES WITH NEEDLES



HSW

Henke-Ject® (Formerly "Soft-Ject®") fixed needle syringes offer a crystal clear cylinder with easy to read scale. The piston, piston-rod, and cylinder consist of polypropylene. The fixed needle syringes are free of pyrogens, PVC, and latex, and are non-toxic. The 0.5cc syringes are in a sterile polybag with 10 pieces each, the 1cc syringes are sterile individually blister packaged.



The By-packs from Henke-Sass, Wolf offer a convenient syringe with needle combination packaged side-by-side in individual strip packs.





HENKE-JECT® FIXED-NEEDLE SYRINGES

Catalog Number	Syringe Size	Tip Type	Sterility	Quantity per Box	Quantity per Case
RU13012	1cc Tuberculin	30g x 1/2"	Sterile	100	2000

HENKE-JECT® SYRINGES / NEEDLE BY-PACKS

Catalog Number	Syringe Size	Tip Type	Sterility	Quantity per Box	Quantity per Case
R22234	2.5mL LS	22g x 3/4"	Sterile	100	2500
R2221	2.5mL LS	22g x 1"	Sterile	100	2500
R22558	2.5mL LS	25g x 5/8"	Sterile	100	2500
R2018	20mL LS	18g x 1-1/2"	Sterile	100	800



EXEL INTERNATIONAL®

SYRINGES WITH NEEDLES



EXEL

EXEL's Comfort Point fixed needle syringes offer a permanently attached needle, reduced dead space and smooth plunger action for greater filling accuracy. The auto-ground needle tip offers greater comfort. These syringes are sterile, non-toxic, and non-pyrogenic. All EXEL fixed needle syringes are packaged 10 per bag and 100 per box.

EXEL syringes with mounted needle offer convenience at an economic price. We carry EXEL's wide selection of both luer slip and luer lock syringes with standard size needles pre-attached. These syringes allow the user to reduce the number of inventory items and there is no need to open multiple packs to assemble a ready-to-use syringe.

EXEL FIXED NEEDLE SYRINGES

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
EU032913	0.3cc Insulin	29g x 1/2"	Sterile	100	500
EU03308	0.3cc Insulin	30g x 5/16"	Sterile	100	500
EU03318	0.3cc Insulin	31g x 5/16"	Sterile	100	500
EU052813	0.5cc Insulin	28g x 1/2"	Sterile	100	500
EU052913	0.5cc Insulin	29g x 1/2"	Sterile	100	500
EU05308	0.5cc Insulin	30g x 5/16"	Sterile	100	500
EU05318	0.5cc Insulin	31g x 5/16"	Sterile	100	500
EU12558	1cc Tuberculin	25g x 5/8"	Sterile	100	500
EU12813	1cc Insulin	28g x 1/2"	Sterile	100	500
EU12913	1cc Insulin	29g x 1/2"	Sterile	100	500
EU1308	1cc Insulin	30g x 5/16"	Sterile	100	500
EU1318	1cc Insulin	31g x 5/16"	Sterile	100	500

EXEL LUER SLIP SYRINGES WITH MOUNTED NEEDLE

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
ES12558	1cc	25g x 5/8"	Sterile	100	1000
ES1251	1cc	25g x 1"	Sterile	100	1000
ES12612	1cc	26g x 1/2"	Sterile	100	1000
ES12712	1cc	27g x 1/2"	Sterile	100	1000
ES3201	3cc	20g x 1"	Sterile	100	1000
ES3211	3cc	21g x 1"	Sterile	100	1000
ES32234	3cc	22g x 3/4"	Sterile	100	1000
ES3221	3cc	22g x 1"	Sterile	100	1000
ES3231	3cc	23g x 1"	Sterile	100	1000
ES32558	3cc	25g x 5/8"	Sterile	100	1000

EXEL LUER LOCK SYRINGES WITH MOUNTED NEEDLE

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
EL318112	3cc	18g x 1-1/2"	Sterile	100	1000
EL3201	3cc	20g x 1"	Sterile	100	1000
EL320112	3cc	20g x 1-1/2"	Sterile	100	1000
EL3211	3cc	21g x 1"	Sterile	100	1000
EL321114	3cc	21g x 1-1/4"	Sterile	100	1000
EL321112	3cc	21g x 1-1/2"	Sterile	100	1000
EL32234	3cc	22g x 3/4"	Sterile	100	1000
EL3221	3cc	22g x 1"	Sterile	100	1000
EL322114	3cc	22g x 1-1/4"	Sterile	100	1000
EL322112	3cc	22g x 1-1/2"	Sterile	100	1000
EL3231	3cc	23g x 1"	Sterile	100	1000
EL323112	3cc	23g x 1-1/2"	Sterile	100	1000
EL32558	3cc	25g x 5/8"	Sterile	100	1000
EL3251	3cc	25g x 1"	Sterile	100	1000
EL325112	3cc	25g x 1-1/2"	Sterile	100	1000
EL327114	3cc	27g x 1-1/4"	Sterile	100	1000
EL5201	5cc	20g x 1"	Sterile	100	800
EL5211	5cc	21g x 1"	Sterile	100	800
EL521112	5cc	21g x 1-1/2"	Sterile	100	800
EL5221	5cc	22g x 1"	Sterile	100	800
EL522112	5cc	22g x 1-1/2"	Sterile	100	800
EL10181	10cc	18g x 1"	Sterile	100	800
EL10201	10cc	20g x 1"	Sterile	100	800
EL10211	10cc	21g x 1"	Sterile	100	800
EL1021112	10cc	21g x 1-1/2"	Sterile	100	800
EL10221	10cc	22g x 1"	Sterile	100	800
EL1022112	10cc	22g x 1-1/2"	Sterile	100	800



BD[®]

SYRINGES WITH NEEDLES



BD

Becton Dickinson fixed needle syringes offer Ultra-Fine™ permanently attached needles. BD insulin syringes are the #1 selling brand in the US.

BD Mounted Needle syringes offer a convenient BD syringe and Precision Glide™ needle - no assembly required. The clear barrels on these syringes feature bold scales and the needle helps improve dosage accuracy.

Air-Tite carries a mixture of BD Safety syringes in the most popular sizes. The Eclipse line allows clinicians to activate the safety mechanism after the injection while the pivoting shield offers low angle injections and efficient safety. The Safety-Lok™ line offers the only syringe of its type and affords clear visual confirmation of the sliding sleeve in the lock position. SafetyGlide™ has a single-hand activated shielding mechanism.

BD FIXED NEEDLE SYRINGES

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box
BD328438	3/10 mL Insulin	31g x 5/16"	Sterile	100
BD328468	1/2 mL Insulin	31g x 5/16"	Sterile	100
BD329424	1 mL Insulin	28g x 1/2"	Sterile	100
BD328418	1 mL Insulin	31g x 5/16"	Sterile	100

BD SYRINGES WITH MOUNTED NEEDLES

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box
BD309623	1mL TB Luer Slip	27g x 1/2"	Sterile	100
BD309637	1mL Luer Lock	20g x 1"	Sterile	100
BD309580	3mL Luer Lock	18g x 1-1/2"	Sterile	100
BD309570	3mL Luer Lock	25g x 5/8"	Sterile	100



BD SAFETY SYRINGES WITH MOUNTED NEEDLES

Catalog Number	Syringe Size	Needle	Sterility	Quantity	Safety Brand
BD305781	3mL	25g x 5/8"	Sterile	300	Eclipse
BD305782	3mL	23g x 1"	Sterile	300	Eclipse



EZ FLO NEEDLES



Smooth Entry • Smooth Delivery

Air-Tite Product's newest addition, the EZ FLO line of hypodermic needles, offers high quality and a broad selection at an affordable price. Ranging from 18 gauge to 30 gauge, these needles offer a standard bevel and regular-wall to combine comfort and durability. Each needle is sterile, individually packaged in a box of 100. Made of surgical stainless steel and a polypropylene hub. 18g to 25g are thin wall. 27g and 30g are regular wall thickness.



EZ FLO BY AIR-TITE

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NEZ18034	18	3/4"	Sterile	100	2000
NEZ18100	18	1"	Sterile	100	2000
NEZ18112	18	1-1/2"	Sterile	100	2000
NEZ18200	18	2"	Sterile	100	2000
NEZ19100	19	1"	Sterile	100	2000
NEZ19112	19	1-1/2"	Sterile	100	2000
NEZ20100	20	1"	Sterile	100	2000
NEZ20112	20	1-1/2"	Sterile	100	2000
NEZ21100	21	1"	Sterile	100	2000
NEZ21112	21	1-1/2"	Sterile	100	2000
NEZ21200	21	2"	Sterile	100	2000
NEZ22034	22	3/4"	Sterile	100	2000
NEZ22100	22	1"	Sterile	100	2000
NEZ22112	22	1-1/2"	Sterile	100	2000
NEZ22200	22	2"	Sterile	100	2000
NEZ22300	22	3"	Sterile	100	2000
NEZ22400	22	4"	Sterile	100	2000
NEZ23034	23	3/4"	Sterile	100	2000
NEZ23100	23	1"	Sterile	100	2000
NEZ23112	23	1-1/2"	Sterile	100	2000
NEZ25058	25	5/8"	Sterile	100	2000
NEZ25100	25	1"	Sterile	100	2000
NEZ25112	25	1-1/2"	Sterile	100	2000
NEZ25200	25	2"	Sterile	100	2000
NEZ27100	27	1"	Sterile	100	2000
NEZ27114	27	1-1/4"	Sterile	100	2000
NEZ27112	27	1-1/2"	Sterile	100	2000
NEZ27200	27	2"	Sterile	100	2000
NEZ30012	30	1/2"	Sterile	100	2000
NEZ30100	30	1"	Sterile	100	2000

Needles

TSK STERIJECT™ NEEDLES



TSK STERIJECT™ Hypodermic Needles are produced with a focus on needle sharpness and quality. All needles are thin-wall or ultra-thin wall. Made in Japan from one of the global leaders for ultra-fine tip needles, TSK, these needles are commonly used for aesthetic procedures, mesotherapy, and ophthalmology among other uses. Available in hard-to-find sizes for the specialist. 100 needles per box, surgical stainless steel, and poly hub.

TSK's Low Dead Space Needles feature a unique hub that is designed to nearly eliminate product loss due to dead space typically seen in a conventional needle hub, which may result in significant cost savings. This hard polymer hub also helps prevent flexing and pop-off, with external hub threading to reduce risk of leaks. Their TSKiD™ large needle barrels have enlarged inner diameters for less extrusion force. Recommended for use with Luer Lock syringes. The INVISIBLE NEEDLE features the low dead space hub on TSK's smallest needle ever! At just 0.2mm in diameter, it is thinner than even a 33G needle.

TSK LOW DEAD SPACE NEEDLES

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box	Quantity per Case
TSK3009LDS	30	9mm	Ultra-Thin	Sterile	100	6000
TSK3013LDS	30	13mm	Ultra-Thin	Sterile	100	6000
TSK3309LDS	33	9mm	Thin	Sterile	100	6000
TSK3313LDS	33	13mm	Thin	Sterile	100	6000
TSK02009LDS (THE INVISIBLE NEEDLE)	0.2mm	9mm	Thin	Sterile	100	6000

TSK STERIJECT™ NEEDLES

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box	Quantity per Case
TSK2513U	25	13mm	Ultra-Thin	Sterile	100	6000
TSK2713U	27	13mm	Ultra-Thin	Sterile	100	6000
TSK2719U	27	19mm	Ultra-Thin	Sterile	100	6000
TSK2738U	27	38mm	Ultra-Thin	Sterile	100	6000
TSK304U	30	4mm	Ultra-Thin	Sterile	100	6000
TSK306U	30	6mm	Ultra-Thin	Sterile	100	6000
TSK3013U	30	13mm	Ultra-Thin	Sterile	100	6000
TSK3113	31	13mm	Thin	Sterile	100	6000
TSK324	32	4mm	Thin	Sterile	100	6000
TSK3206	32	6mm	Thin	Sterile	100	6000
TSK3209C	32	9mm	Thin	Sterile	100	6000
TSK3213	32	13mm	Thin	Sterile	100	6000
TSK3304	33	4mm	Thin	Sterile	100	6000
TSK3313	33	13mm	Thin	Sterile	100	6000



EXEL INTERNATIONAL®

NEEDLES



EXEL International® offers a variety of hypodermic needles with regular wall thickness, made of 304 surgical stainless steel, non-toxic and pyrogen-free. Each box is color coded to help the clinician easily identify the appropriate gauge. Each needle is individually packaged, sterilized, and ready to use.

EXEL INTERNATIONAL® HYPODERMIC NEEDLES

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
EN181	18	1"	Sterile	100	2,000
EN18112	18	1-1/2"	Sterile	100	2,000
EN191	19	1"	Sterile	100	2,000
EN19112	19	1-1/2"	Sterile	100	2,000
EN201	20	1"	Sterile	100	2,000
EN20112	20	1-1/2"	Sterile	100	2,000
EN211	21	1"	Sterile	100	2,000
EN21112	21	1-1/2"	Sterile	100	2,000
EN2234	22	3/4"	Sterile	100	2,000
EN221	22	1"	Sterile	100	2,000
EN22112	22	1-1/2"	Sterile	100	2,000
EN2334	23	3/4"	Sterile	100	2,000
EN231	23	1"	Sterile	100	2,000
EN23112	23	1-1/2"	Sterile	100	2,000
EN2558	25	5/8"	Sterile	100	2,000
EN251	25	1"	Sterile	100	2,000
EN25112	25	1-1/2"	Sterile	100	2,000
EN2612	26	1/2"	Sterile	100	2,000
EN2712	27	1/2"	Sterile	100	2,000
EN27114	27	1-1/4"	Sterile	100	2,000
EN27112	27	1-1/2"	Sterile	100	2,000
EN2834	28	3/4"	Sterile	100	1,000
EN3012	30	1/2"	Sterile	100	2,000
EN301	30	1"	Sterile	100	2,000
EN3212	32	1/2"	Sterile	100	1,000



Needles



SecureTouch needles feature one-handed, finger-stroke activation. No solid surface required. An audible click as well as tactile feedback and visual indicator confirm safety mechanism engagement. Needle remains locked inside activated cover. Fits any standard brand syringe.


EXEL® SECURETOUCH SAFETY NEEDLES

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
ENS18112	18	1-1/2"	Sterile	100	1.000
ENS201	20	1"	Sterile	100	1.000
ENS231	23	1"	Sterile	100	1.000
ENS2558	25	5/8"	Sterile	100	1.000
ENS251	25	1"	Sterile	100	1.000
ENS25112	25	1-1/2"	Sterile	100	1.000
ENS27114	27	1-1/4"	Sterile	100	1.000
ENS27112	27	1-1/2"	Sterile	100	1.000
ENS3012	30	1/2"	Sterile	100	1.000

HENKE-JECT® NEEDLES



Henke-Ject® (formerly Fine-Ject®) hypodermic needles by Henke-Sass, Wolf GmbH are high quality yet economical. Made of surgical stainless steel, these needles are sterile, disposable, individually wrapped, and available in thin wall up to 25 gauge.



HENKE-JECT® HYPODERMIC NEEDLES BY HENKE-SASS, WOLF GMBH

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box	Quantity per Case
NH15112	15	1-1/2"	Thin	Sterile	100	4,000
NH172	17	2"	Thin	Sterile	100	3,000
NH18112	18	1-1/2"	Thin	Sterile	100	4,000
NH182	18	2"	Thin	Sterile	100	3,000
NH201	20	1"	Thin	Sterile	100	4,000
NH20112	20	1-1/2"	Thin	Sterile	100	4,000
NH212	21	2"	Thin	Sterile	100	3,000
NH221	22	1"	Thin	Sterile	100	4,000
NH22112	22	1-1/2"	Thin	Sterile	100	4,000
NH222	22	2"	Thin	Sterile	100	3,000
NH22234	22	2-3/4"	Thin	Sterile	100	3,000
NH224	22	4"	Thin	Sterile	100	2,000
NH2558	25	5/8"	Thin	Sterile	100	4,000
NH251	25	1"	Thin	Sterile	100	4,000
NH252	25	2"	Thin	Sterile	100	3,000
NH272	27	2"	Regular	Sterile	100	3,000
NH3012	30	1/2"	Regular	Sterile	100	4,000



Needles



JBP NANONEEDLES



JBP (S. Korea) Nanoneedles™ are the integration of craftsmanship and technology which has led to an ultra-thin wall needle making injections of high viscosity materials through extra fine needles possible. JBP uses a special process called Micron Sharpness to reduce frictional forces. With gauges ranging from 30 to 34, these needles are ideally suited for those seeking ultra-fine tips and the smallest of injection needles.



JBP NANONEEDLES™


Catalog Number	Gauge	Length	Wall Type	Package Type	Sterility	Quantity per Box	Quantity per Case
JBP3004	30	4mm	Ultra-thin	Capsule	Sterile	100	800
JBP3004B	30	4mm	Ultra-thin	Blister	Sterile	100	800
JBP3008	30	8mm	Ultra-thin	Capsule	Sterile	100	800
JBP3304B	33	4mm	Ultra-thin	Blister	Sterile	100	800
JBP3308	33	8mm	Ultra-thin	Capsule	Sterile	100	800
JBP3308B	33	8mm	Ultra-thin	Blister	Sterile	100	800
JBP3313	33	13mm	Ultra-thin	Capsule	Sterile	100	800
JBP3313B	33	13mm	Ultra-thin	Blister	Sterile	100	800
JBP3404	34	4mm	Ultra-thin	Capsule	Sterile	100	800
JBP3406	34	6mm	Ultra-thin	Capsule	Sterile	100	800
JBP3408	34	8mm	Ultra-thin	Capsule	Sterile	100	800

*Additional sizes available. Please call for more information.



CONVENTIONAL NEEDLE

VS



JBP NANONEEDLE (UTW)

BD[®] NEEDLES



BD Precision Glide™ needles are available in a wide range of gauges, lengths and bevel designs.

BD Safety Glide™ Needles feature BD Activation Assist™ technology making a safe injection easier and more efficient.





BD® PRECISION GLIDE™ NEEDLES

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box
BD305198	16	1-1/2"	Regular	Sterile	100
BD305216	16	1"	Thin	Sterile	100 NOKOR ADMIX
BD305195	18	1"	Regular	Sterile	100
BD305196	18	1-1/2"	Regular	Sterile	100
BD305215	18	1-1/2"	Thin	Sterile	100 NOKOR ADMIX
BD305180	18	1-1/2"	Regular	Sterile	100 (Blunt tip)
BD305167	21	1-1/2"	Regular	Sterile	100
BD305129	21	2"	Regular	Sterile	100
BD305145	23	1"	Regular	Sterile	100
BD305120	23	1-1/4"	Regular	Sterile	100
BD305122	25	5/8"	Regular	Sterile	100
BD305125	25	1"	Regular	Sterile	100
BD305127	25	1-1/2"	Regular	Sterile	100
BD305110	26	3/8"	Regular	Sterile	100
BD305109	27	1/2"	Regular	Sterile	100
BD305106	30	1/2"	Regular	Sterile	100
BD305128	30	1"	Regular	Sterile	100



BD® SAFETY GLIDE NEEDLES

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box
BD305915	21	1"		Sterile	500
BD305902	23	1"		Sterile	500
BD305901	25	5/8"		Sterile	500
BD305916	25	1"		Sterile	500



EXEL INTERNATIONAL®

SPINAL NEEDLES



Quincke bevel spinal needles. Siliconized for smooth, comfortable penetration. Transparent hub. Latex free. Convenient blister packaging.



EXEL® SPINAL NEEDLES


Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NE183112A	18	3-1/2"	Sterile	50	100
NE203112A	20	3-1/2"	Sterile	50	100
NE22234A	22	2-3/4"	Sterile	50	100
NE223112A	22	3-1/2"	Sterile	50	100
NE253112A	25	3-1/2"	Sterile	50	100



BD[®] SPINAL NEEDLES







The key/slot arrangement of stylet and cannula hubs on the BD[®] Quincke spinal needles facilitates proper needle bevel orientation. Fitted stylet designed to reduce tissue coring. Unique translucent window hub featuring contact clarity that helps allow visualization of CSF.



BD[®] SPINAL NEEDLES


Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
BD408360	18	6"	Sterile	10	50
BD405182	20	3-1/2"	Sterile	25	100
BD405161	22	1-1/2"	Sterile	N/A	100
BD405181	22	3-1/2"	Sterile	25	100
BD405148	22	5"	Sterile	10	50
BD405149	22	7"	Sterile	10	50



TSK MICROCANNULAS



TSK STERiJECT™ CSH CANNULA




VALUE PRICED CANNULA

- Traditional Product Delivery
- Ultra Thin Wall



ROUNDED, CLOSED-TIP DESIGN

- Reduced Risk of Trauma
- Easier Handling



Conventional needle

ULTRA THIN WALL

- Less force required for injection
- Smoother and More Controlled Injection




Regular Wall

CODE #	DESCRIPTION	DIMENSIONS
TSK2250SP25	25-PACK STERiJECT CSH Microcannulas	22g x 50mm
TSK2538SP25	25-PACK STERiJECT CSH Microcannulas	25g X 38mm
TSK2550SP25	25-PACK STERiJECT CSH Microcannulas	25g x 50mm
TSK2738SP25	25-PACK STERiJECT CSH Microcannulas	27g x 38mm
TSK2750SP25	25-PACK STERiJECT CSH Microcannulas	27g x 50mm

Introducer pilot needles sold separately.

TSK STERIGLIDE™ CANNULA

AWARD WINNING BEST INJECTION TOOL



CODE #	DIMENSIONS
TSK2238SG	22g x 38mm (1½")
TSK2250SG	22g x 50mm (2.0")
TSK2270SG	22g x 70mm (2¾")
TSK2538SG	25g x 38mm (1½")
TSK2550SG	25g x 50mm (2.0")
TSK2725SG	27g x 25mm (1.0")
TSK2738SG	27g x 38mm (1½")
TSK2750SG	27g x 50mm (2.0")
TSK3025SG	30g x 25mm (1.0")

Each box contains 20 microcannulas and pilot needles.

Microcannulas


PREMIUM CANNULA

Smooth Gliding, Near Tip Product Delivery, Ultra-Thin Wall

Flexible enough to navigate anatomy, rigid enough to place accurately, with improved tactile sensitivity


DOME SHAPED TIP DESIGN

- Improved cannula introduction
- Reduced risk of trauma
- Easier handling




NEAREST-TO-TIP DELIVERY

- More Accurate Placement
- Reduced Product Loss




PROPRIETARY SURFACE TREATMENT

- Less force to advance
- Reduced Patient Discomfort



CLEAR SIDE-PORT MARKING

- Directional accuracy
- Better Handling



INTRAVENOUS SETS



EXEL IV ADMINISTRATION

Catalog Number	Specs	Sterility	Quantity Per Case
EIV7815	78" Tubing, 15 drops/mL, vented / non-vented, Luer Slip	Sterile	50
EIV7815LL	78" Tubing, 15 drops/mL, vented / non-vented, Luer Lock	Sterile	50
EIV7815-2C	78" Tubing, 15 drops/mL, vented / non-vented, Luer Lock/Luer Slip with Dual Clamp	Sterile	50
EIV7860	78" Tubing, 60 drops/mL, vented / non-vented, Luer Slip	Sterile	50
EIV7860LL	78" Tubing, 60 drops/mL, vented / non-vented, Luer Lock	Sterile	50
EIV8420F	0.2 Micron IV filter, 4" tubing, roller clamp and occlusion clamp	Sterile	50
EIV10515	Mullti-site, 2 Y-injection, 105" tubing, 28mL priming volume	Sterile	50


EXEL IV CATHETERS

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EIV14114	14g x 1-1/4"	Sterile	50	200
EIV142	14g x 2"	Sterile	50	200
EIV16114	16g x 1-1/4"	Sterile	50	200
EIV162	16g x 2"	Sterile	50	200
EIV18114	18g x 1-1/4"	Sterile	50	200
EIV18112	18g x 1-1/2"	Sterile	50	200
EIV182	18g x 2"	Sterile	50	200
EIV201	20g x 1"	Sterile	50	200
EIV20114	20g x 1-1/4"	Sterile	50	200
EIV202	20g x 2"	Sterile	50	200
EIV221	22g x 1"	Sterile	50	200
EIV2434	24g x 3/4"	Sterile	50	200



EXEL WINGED INFUSION SETS

Catalog Number	Specs	Sterility	Quantity Per Box	Quantity Per Case
EWI1934	19g x 3/4" w/ 12" Tubing	Sterile	50	500
EWI2034	20g x 3/4" w/ 12" Tubing	Sterile	50	500
EWI2134	21g x 3/4" w/ 12" Tubing	Sterile	50	500
EWI2334	23g x 3/4" w/ 12" Tubing	Sterile	50	500
EWI2534	25g x 3/4" w/ 12" Tubing	Sterile	50	500
EWI2734	27g x 3/4" w/ 12" Tubing	Sterile	50	500



EXEL SAFETY WINGED INFUSION SETS

Catalog Number	Specs	Sterility	Quantity Per Box	Quantity Per Case
EWI2134S	21g x 3/4" Safety w/ 12" Tubing	Sterile	50	500
EWI2334S	23g x 3/4" Safety w/ 12" Tubing	Sterile	50	500
EWI2534S	25g x 3/4" Safety w/ 12" Tubing	Sterile	50	500



EXEL INJECTION PLUGS

Catalog Number	Specs	Sterility	Quantity Per Box	Quantity Per Case
EIP2	Injection Plug with Cap Luer Lock	Sterile	50	500

BLOOD COLLECTION



BD VACUTAINER® BLOOD COLLECTION TUBES

Catalog Number	Specs	Sterility	Quantity per Box
BD366668	3mL Red, Plastic, No Additive, w/ Silicone	Sterile	100
BD367812	4mL Red, Plastic, No Additive, w/ Silicone	Sterile	100
BD367815	6mL Red, Plastic, No Additive, w/ Silicone	Sterile	100
BD366430	10mL Red, Glass, No Additive, w/ Silicone	Sterile	100
BD367985	10mL Red/Grey, Plastic, Plus SST Brand	Sterile	100





EXEL BLOOD COLLECTION SETS WITH 12" TUBING

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EBS2134	21g x 3/4"	Sterile	50	200
EBS2334	23g x 3/4"	Sterile	50	200
EBS2534	25g x 3/4"	Sterile	50	200



EXEL MULTI-TYPE BLOOD COLLECTION NEEDLES

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EBC201	20g x 1"	Sterile	100	1000
EBC20112	20g x 1-1/2"	Sterile	100	1000
EBC211	21g x 1"	Sterile	100	1000
EBC21112	21g x 1-1/2"	Sterile	100	1000
EBC221	22g x 1"	Sterile	100	1000
EBC22112	22g x 1-1/2"	Sterile	100	1000



IV & Blood Collection

EXEL BLOOD COLLECTION ADAPTOR

For use with winged infusion sets

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EWI2BC	Adaptor	Sterile	100	1000



EXEL BLOOD COLLECTION HOLDER WITH MULTI-SAMPLE ADAPTOR

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EBCA	Holder with Adaptor	Sterile	50	200

SHARPS COLLECTION AND MAIL-BACK SYSTEMS



Bemis® SharpSentinel® Containers

Part Number	Specifications	Quantity	Quantity per Case
SC101R	1 Quart (4.25"L x 4.25"W x 6.5"H - fits phlebotomy trays)	Each	100
SC105R	5 Quart Rotating Lid (13"L x 5.5"W x 9"H - 1.25 gallons - wall mount or cabinet)	Each	32
SC108R	8 Quart (11.75"L x 7.75"W x 8.75"H - 2 gallons - table top or wall mount)	Each	30
SC112R	12 Quart (15"L x 6.5"W x 15"H - 3 gallons)	Each	12
SC212R	12 Quart Rotating Lid (15"L x 6.5"W x 15"H - 3 gallons)	Each	12
SC132R	8 Gallon (15.25"L x 11"W x 15.5"H)	Each	10
SC144R	11 Gallon (15.25"L x 11"W x 22"H)	Each	6
SC101H	1 Quart Holder (fits SC101R)	Each	20
SC100B	SharpSentinel® Mounting Bracket (fits 5 quart - SC105R & 8 quart - SC108R)	Each	10
SC100K	SharpSentinel® Bracket Key - fits SC100B	Each	10
SC100W	3 Gallon Wire Bracket (fits 12 quart SC112R & SC212R)	Each	5
SC5QWC	5 Quart Wall Cabinet	Each	1
SC3GWC	3 Gallon Wall Cabinet	Each	1
SCGH1	Glove Box Holder (snaps onto bottom of SC5QWC & SC3GWC)	Each	1

Non-Sharps Medical Waste Mail-Back Systems **May not contain any LOOSE sharps

Part Number	Specifications	Dimensions
MWM1G	1 Gallon System	8.375" diameter x 8" high
SCM5G	5 Gallon System	10.45" diameter x 14.51" high

Sharps Mail-Back System

Each Kit includes all components required for collection and transportation of regulated medical waste.

Each Kit includes:

Step-by-step instructions, waste shipping containers, prepaid and pre-addressed shipping labels, prepaid disposal & tracking documents.



1 quart (QTY: 5)

1 gallon

2 gallon

3 gallon

Part Number	Specification	Dimension
SCM1Q	1 Quart System	Container Dimensions: 3.75" x 3.5" x 8" Each container holds approximately 42 3cc syringes.
SCM1Q-5	1 Quart System (QTY: 5)	Container Dimensions: 3.75" x 3.5" x 8" Each container holds approximately 42 3cc syringes.
SCM1G	1 Gallon System	Container Dimensions: 9" x 5.5" x 8.5" Each container holds approximately 98 3cc syringes.
SCM1G-2	1 Gallon System (QTY: 2)	Container Dimensions: 9" x 5.5" x 8.5" Each container holds approximately 98 3cc syringes.
SCM2G	2 Gallon System	Container Dimensions: 9" x 5.5" x 14.25" Each container holds approximately 245 3cc syringes.
SCM2G-2	2 Gallon System (QTY: 2)	Container Dimensions: 9" x 5.5" x 14.25" Each container holds approximately 245 3cc syringes.
SCM3G	3 Gallon System	Container Dimensions: 9" x 5.5" x 19.75" Each container holds approximately 358 3cc syringes.

Universal Brackets

Universal mounting brackets fit 1 quart, 1 gallon, 2 gallon, and 3 gallon Mail-Back systems (see above). Locking brackets include lock & key.

Part Number	Specification
SCMSB123	Universal Bracket
SCMSB123L	Universal Bracket - Locking

GLOVES



Latex Exam Gloves *100 gloves/box.*

Part Number	Specifications	Quantity	Quantity per Case
LS100S	Powder FREE Small	100	1000
LS100M	Powder FREE Medium	100	1000
LS100L	Powder FREE Large	100	1000
LS100X	Powder FREE XL	100	1000

Vinyl Exam Gloves *100 gloves/box.*

Part Number	Specifications	Quantity	Quantity per Case
VS100S	Powder FREE Small	100	1000
VS100M	Powder FREE Medium	100	1000
VS100L	Powder FREE Large	100	1000
VS100X	Powder FREE XL	100	1000

Latex High Risk Gloves

50 gloves/box. Extra long, textured.

Part Number	Specifications	Quantity	Quantity Per Case
LD300S	Powder FREE Small	50	500
LD300M	Powder FREE Medium	50	500
LD300L	Powder FREE Large	50	500
LD300X	Powder FREE Extra Large	50	500

Nitrile Textured Gloves

100 gloves/box

Part Number	Specifications	Quantity	Quantity Per Case
KT100S	Powder FREE Small	100	1000
KT100M	Powder FREE Medium	100	1000
KT100L	Powder FREE Large	100	1000
KT100X	Powder FREE Extra Large	100	1000

Black Nitrile Gloves

100 gloves/box

Part Number	Specifications	Quantity	Quantity Per Case
KB100S	Powder FREE Small	100	1000
KB100M	Powder FREE Medium	100	1000
KB100L	Powder FREE Large	100	1000
KB100X	Powder FREE Extra Large	100	1000

Pink Nitrile Gloves

200 gloves/box

Part Number	Specifications	Quantity	Quantity Per Case
KP200XS	Powder FREE Extra Small	200	2000
KP200S	Powder FREE Small	200	2000
KP200M	Powder FREE Medium	200	2000
KP200L	Powder FREE Large	200	2000
KP200XL	Powder FREE Extra Large	200	2000

Air-Tite
Products Co., Inc.



Supplies


SUPPLIES

Air-Tite Products Co. offers a wide range of supplies for the clinicians practice. From Sharps Containers to Personal Protective Equipment, from gauze to Instruments - Air-Tite carries these products for your convenience. Please see the current Price List or visit our website at www.air-tite.com for part numbers, pricing, ordering and more information.

HIGHLIGHT PRODUCTS

TSK NEEDLE HOLDER is a high quality, reusable needle stand designed to allow a safe procedure for engaging the syringe to the needle, while also providing a convenient holder during the procedure.

Part # **TSKNH1**



SYRINGE CONNECTORS are high quality female-to-female adapters with low dead space. Compatible with both Luer Slip and Luer Lock syringes. Made in the USA.

Part # **IM57-400**



CURACATOR™ designed and developed by Dr. Janine Hopkins, board-certified dermatologist, to meet the need for hands-free, needle-free, controlled application of products to skin. Made in the USA, the patented CuraCator™ device is a sterile, single-use only applicator that simply attaches to any Luer lock syringe, providing a solution that alleviates safety concerns and risks of product waste or contamination. Dripping products such as platelet-rich plasma, exosomes, or other serums from syringes through a needle is dangerous. Applying products by dripping them onto tissue from an open-ended syringe leads to waste and little control. Wiping product onto tissue with a gloved hand can lead to contamination or adulteration of product. The goal of the CuraCator™ is to become the new standard of care for allowing product (whether liquid, cream, or ointment), to be safely and accurately applied to the skin or other tissue. Made from durable Triton plastic. FDA Class 1 Medical Device Accessory.

Part # **AESCURA**



GENERAL SUPPLIES

Access to nearly the entire range of Miltex surgical instruments including scissors, needle holders, forceps, clamps and many more.



Aesthetic marking pencils are essential for any aesthetician to have in their kit in order to mark areas of interest. White cosmetic wax pencils have been used for marking prior to non-surgical aesthetic procedures and during consultations as the coloring can be easily removed.

Part # **AESPEN**



Gauze is available in sterile, unsterile, woven, and non-woven options.



Supplies


In addition to the supplies mentioned above, Air-Tite Products offers Alcohol and Iodine Prep Pads, Cotton Products (Rolls, Balls, Tip Applicators), Disinfectant wipes, Hand Creams, Tongue depressors, Cohesive bandage wraps, Thermometers, Stethoscopes, Instrument Trays, Lubricating Jelly, and more. Please see the price list for our full offerings, or call our Customer Service team today!



LOW DEAD SPACE (LDS) NEEDLES

TSK


Compared to conventional needles, the LDS hub design minimizes product loss which may lead to significant cost savings when injecting high cost drugs.



LDS Hub

Regular Hub

The LOW DEAD SPACE needle hub is designed to reduce the effective material left behind as close to zero as possible.



THE INVISIBLE NEEDLE™ features the LDS needle hub on TSK's thinnest needle ever!



www.air-tite.com



atinfo@air-tite.com

800.231.7762

Welcome to the **Air-Tite** family!

We supply the world with high quality disposable + reusable hypodermic products for Medical, Laboratory, and Veterinary use offering personalized service with a 100% satisfaction guarantee. Our sales and service team are always available to answer your questions and guide your selection with product that improve your business.

In addition to our own line of proven products, we represent and distribute other manufacturers such as:

Henke-Sass, Wolf GmbH (Germany), B. Braun (Germany),
Becton-Dickinson (USA), EXEL (USA), TSK laboratory (Japan),
Poulsen & Graf (Germany), plus many more.

TO REQUEST A CATALOG OR SAMPLES, OR TO PLACE
AN ORDER VISIT

www.Air-Tite.com
Email atinfo@Air-Tite.com
or call 1.800.231.7762
office: 757.340.2501 fax: 757.340.2912
We guarantee our quality 100%

HOURS: MON - FRI 8:30 AM - 6:00 PM Eastern Time

Orders placed by 2:00PM Eastern Time are packaged and shipped the same day

Orders for delivery in continental US are shipped UPS Ground, prepaid and added to your invoice. Expedited delivery available upon request.

Orders outside the continental US shipped Air Freight prepaid and added to your invoice.

Export sales welcomed. Please call for details.

HEALTHCARE ACCOUNTS REPRESENTATIVE

Sarah Villarreal

SarahV@air-tite.com

Air-Tite Human Health Distribution Partners

Accuderm
Cardinal Health
Cobalt Medical Supply, Inc.
Delasco
Henry Schein
McKesson
Medical Spa Supply
Merit Pharmaceutical
Rally, Inc.
The Torrance Company
& More!

ICE-STICK

AESTHETIC COOLING DEVICE

Aesthetic
Group

This cooling device from Aesthetic Group in France provides a great non-drip solution in pre- and post-procedure for your patients. The ICE-STICK core can be removed from the plastic casing and placed in your freezer until needed. Used frequently in aesthetic medicine practices, this ICE-STICK can be cleaned and reused.

The ICE-STICK is composed of an ice-cold cartridge with refreshing aluminum head to cool areas of the face and body. Helpful prior to injection for pain reduction and to reduce the risk of hematomas.

PART # ICESTICK



Air-Tite
Products Co., Inc.



Scan For Website

