

2024 | CATALOG


YOUR SYRINGE AND NEEDLE SPECIALISTS SINCE 1926

Air-Tite.com

T: 800.231.7762 F: 757.340.2912 L: 757.340.2501

565 Central Drive, Virginia Beach, VA 23454



ANIMAL HEALTH 

Highest Quality • Economical Price • Low Dead Space

DESIGNED FOR VETERINARIANS

1mL Luer Lock Syringes

with and without needles

1ml Luer-Lock Syringe without needle #ML1
Box of 100

1ml Luer-Lock Syringe with 25g x 5/8" needle #ML12558
Box of 100

Extra secure needle connection

Strong regular wall 25g x 5/8" needle

Air-Tite
Products

ONE-STOP SHOPPING

In addition to our needle and syringe products, we also carry gloves, gowns, blood collection tubes, scalpels, tape, gauze, vials, bottles and much more.



50mL All Plastic 2-Part Syringes

Now Accepting Orders!



Luer Lock #AML50
Luer Slip #AMS50

Bulk Unsterile Luer Lock #AMB50LL
Bulk Unsterile Luer Slip #AMB50LS

RETURN POLICY

For over 95 years, Air-Tite has served its customers with the highest quality products at the lowest prices. If, however, you need to return an item, simply contact Customer Service at 1-800-231-7762 for return authorization. For full credit on the purchase price of stocked items, return only full, unopened and undamaged boxes with at least 18 months remaining before product expiration. Non-stocked/special-order items are subject to the return policy of the manufacturer, which may include a restocking fee or other restrictions.

TABLE OF CONTENTS

SYRINGES

Air-Tite Syringes - Sterile	1
2-Part All-Plastic Syringes - Sterile	3
Henke-Ject® 3-Part Syringes - Sterile	5
EXEL International® Syringes - Sterile	7
Bulk Air-Tite Syringes - Unsterile	9
Bulk 2-Part All-Plastic Syringes - Unsterile	11
Bulk EXEL International® Syringes - Unsterile	13
Oral Dose Syringes - Unsterile	15
Repeater/Multi-Dose Syringes	17

SYRINGES WITH NEEDLES

Air-Tite Syringes with Needles	19
HSW Syringes with Needles	21
BD® Syringes with Needles	22
EXEL International® Syringes with Needles	23

NEEDLES

Air-Tite Needles	25
EXEL International® Needles	27
EZ FLO Needles	29
Henke-Ject® Needles	31
TSK STERJECT® Needles	33
Metal Hub Needles	35
Spinal Needles	36

BLOOD COLLECTION

Blood Collecting Needles	37
Blood Collection Tubes, Sets, and Adapters	38

MISCELLANEOUS SUPPLIES

I.V. Sets and Infusion Sets	39
Catheters	40
Pneumatic Vaccinator	41
EPIG® Needleless Vaccinator	42
Syringe Caps and Dispensing Tips	43
Sharps Containers and Mailback Systems	44
Surgical Blades, Handles, Scaples, Sutures	45
Scissors	46
Clamps and Needle Holders	47
Forceps, Shears, Spay & Neuter Pack	48
Clippers, Instrument Trays, Sterilization	49
Pads, Gauze, and Cotton	50
Tape and Bandage Wraps	51
Gloves	52
Wipes, Drapes, and Covers	53
Medicine Vials, Biopsy Punches, Stethoscope	54

**SCAN FOR
PRICE LIST**




Air-Tite Products is happy to assist customers by completing quality surveys and other documentation. We reserve the right to charge a nominal fee for extensive requests.

AIR-TITE BRAND STERILE SYRINGES

Air-Tite
Products

Air-Tite Brand Veterinary and Laboratory syringes are the crossroads where quality meets economy. A basic 3-part syringe with consistent quality is ideal for applications that simply require a low cost transfer vessel. Made to industry standard specifications, these syringes are latex free, non-pyrogenic, and non-toxic. The design includes quality print and siliconized gaskets for a consistently smooth plunger motion and a positive plunger stop.

NOT FDA REGISTERED FOR HUMAN USE





Air-Tite
Products Co., Inc.

STERILE AIR-TITE BRAND

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
MS1	1mL Tuberculin	Luer Slip	Sterile	100	1000
MS3	3mL Centric	Luer Slip	Sterile	100	1000
MS5	5mL Centric	Luer Slip	Sterile	100	1000
MS10	10mL Centric	Luer Slip	Sterile	100	1000
MS20	20mL Eccentric	Luer Slip	Sterile	50	500
MS30	30mL Eccentric	Luer Slip	Sterile	50	500
MS60	60mL Eccentric	Luer Slip	Sterile	25	250



Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
ML1	1mL	Luer Lock	Sterile	100	1000
ML3	3mL	Luer Lock	Sterile	100	1000
ML5	5mL	Luer Lock	Sterile	100	1000
ML10	10mL	Luer Lock	Sterile	100	1000
ML20	20mL	Luer Lock	Sterile	50	500
ML30	30mL	Luer Lock	Sterile	50	500
ML60	60mL	Luer Lock	Sterile	25	250



2-PART SYRINGES NORM-JECT® / AIR-TITE PREMIUM

B | BRAUN

Air-Tite
Products

ALL PLASTIC SYRINGES



2-Part syringes are latex free, contain no rubber or silicone oil. These syringes are the choice for any situation needing an inert, hypoallergenic syringe. Because of their purity, applications by users have included IVF, embryo-transfer, chromatography, and many laboratory procedures. Manufactured only from laboratory grade polypropylene and polyethylene, these unique plastic syringes have a positive safety stop to help prevent accidental spills. These syringes are individually sterile strip packed and are also available in bulk, unsterile packaging. Nearly all 1mL - 20mL syringes are Norm-Ject® brand from B Braun Melsungen AG. 1mL Luer Lock and all 50mL syringes are Air-Tite Premium brand.



CONVENTIONAL RUBBER TIP SYRINGES

Although in general use for medical injections, rubber tipped disposable syringes have several drawbacks compared to the 2-part all plastic syringes. Since they contain natural or synthetic rubber and silicone oil, they may contaminate any media exposed to the syringe. Adverse effects from rubber curing agents, which contain zinc, sulfur and other elements, have been observed on everything from chromatographs to in-vitro procedures. Synthetic rubber plungers are also much more reactive to organic solvents and the silicone oil can cause separation of emulsions. In short, your best choice is a 2-Part syringe if you are interested in purity and accuracy.






STERILE 2-PART SYRINGES

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
A1	1mL Tuberculin	Luer Slip	Sterile	100	1800
A3	3mL Centric	Luer Slip	Sterile	100	2500
A5	5mL Centric	Luer Slip	Sterile	100	2000
A10	10mL Eccentric	Luer Slip	Sterile	100	1200
A20	20mL Eccentric	Luer Slip	Sterile	100	800
AMS50	50mL Eccentric	Luer Slip	Sterile	30	360

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
AML1	1mL	Luer Lock	Sterile	100	2000
AL3	3mL	Luer Lock	Sterile	100	2500
AL5	5mL	Luer Lock	Sterile	100	2000
AL10	10mL	Luer Lock	Sterile	100	1200
AL20	20mL	Luer Lock	Sterile	100	800
AML50	50mL	Luer Lock	Sterile	30	360

** A1 & AML1 Syringes have Low Dead Space Plunger Tip



HENKE-JECT® STERILE CONVENTIONAL SYRINGES



STERILE HENKE-JECT® 3-PART SYRINGES

Henke-Sass, Wolf GmbH Henke-Ject® (formerly Soft-Ject®) syringes are a high quality disposable rubber piston syringe offering substantial savings. A lubricated synthetic rubber tipped piston allows for smooth draw while a positive stop helps to prevent accidental pullout. These syringes are ETO gas sterilized and are latex and pyrogen free. All Henke-Ject® syringes are packaged in individual sterile strip packs. FDA registered, CE Mark, ISO 13485 certified.



STERILE HENKE-JECT®

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
R1**	1mL Tuberculin	Luer Slip	Sterile	100	2000
R2	2.5mL Centric	Luer Slip	Sterile	100	2500
R5	5mL Centric	Luer Slip	Sterile	100	2000
R10	10mL Eccentric	Luer Slip	Sterile	100	1200
R20	20mL Eccentric	Luer Slip	Sterile	100	800
R30HJ	30mL Eccentric	Luer Slip	Sterile	100	600

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
RL1 **	1mL	Luer Lock	Sterile	100	2000
RL3	3mL	Luer Lock	Sterile	100	2400
RL5	5mL	Luer Lock	Sterile	100	2000
RL10	10mL	Luer Lock	Sterile	100	1200
RL20	20mL	Luer Lock	Sterile	100	800
RL30HJ	30mL	Luer Lock	Sterile	100	600

**R1 and RL1 have low dead space plunger tip



EXEL INTERNATIONAL® STERILE CONVENTIONAL SYRINGES

EXELINT®

EXEL is known for unsurpassed quality at a low price. These syringes are latex free, sterile, non-pyrogenic, and non-toxic. The design includes bold print, clear barrels, and siliconized gaskets for smooth, consistent plunger motion and a positive stop.

1mL syringes utilize a low dead space plunger

5mL syringes graduated to 6mL

10mL syringes graduated to 12mL

20mL syringes graduated to 25mL

30mL syringes graduated to 35mL

50mL syringes graduated to 60mL



EXEL INTERNATIONAL®


Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
ES1	1mL Tuberculin	Luer Slip	Sterile	100	1000
ES3	3mL Centric	Luer Slip	Sterile	100	1000
ES5	5mL (6mL) Centric	Luer Slip	Sterile	100	800
ES10	10mL (12mL) Centric	Luer Slip	Sterile	100	800
ES20	20mL (25mL) Eccentric	Luer Slip	Sterile	50	500
ES30	30mL (35mL) Eccentric	Luer Slip	Sterile	50	250
ES60	50mL (60mL) Eccentric	Luer Slip	Sterile	25	150

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
EL1	1mL	Luer Lock	Sterile	100	1000
EL3	3mL	Luer Lock	Sterile	100	1000
EL5	5mL (6mL)	Luer Lock	Sterile	100	800
EL10	10mL (12mL)	Luer Lock	Sterile	100	800
EL20	20mL (25mL)	Luer Lock	Sterile	50	500
EL30	30mL (35mL)	Luer Lock	Sterile	50	250
EL60	50mL (60mL)	Luer Lock	Sterile	25	150

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
EC30	30mL	Catheter Tip	Sterile	50	250
EC60	60mL	Catheter Tip	Sterile	25	150



Air-Tite
Products Co., Inc.





AIR-TITE BRAND BULK UNSTERILE SYRINGES

Air-Tite
Products

Air-Tite Brand Veterinary and Laboratory bulk, unsterile syringes are our most economical non-sterile option. A general purpose, 3-part syringe for applications that simply require a low cost transfer vessel. Made to our exacting specifications, these syringes are latex free, non-pyrogenic, and non-toxic. The design includes quality print and siliconized gaskets for a consistently smooth plunger motion and a positive plunger stop.

NOT FDA REGISTERED FOR HUMAN USE





UNSTERILE AIR-TITE BRAND

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Bag	Quantity per Case
MS1B	1mL Tuberculin	Luer Slip	Unsterile	450	4500
MSB3	3mL Centric	Luer Slip	Unsterile	350	3500
MSB5	5mL Centric	Luer Slip	Unsterile	200	2000
MSB10	10mL Centric	Luer Slip	Unsterile	130	1300
MSB20	20mL Eccentric	Luer Slip	Unsterile	50	500
MSB30	30mL Eccentric	Luer Slip	Unsterile	50	500

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Bag	Quantity per Case
MLB1	1mL	Luer Lock	Unsterile	450	4500
MLB3	3mL	Luer Lock	Unsterile	350	3500
MLB5	5mL	Luer Lock	Unsterile	200	2000
MLB10	10mL	Luer Lock	Unsterile	130	1300
MLB20	20mL	Luer Lock	Unsterile	50	500
MLB30	30mL	Luer Lock	Unsterile	50	500
MLB60	60mL	Luer Lock	Unsterile	25	250



2-PART SYRINGES UNSTERILE BULK / MINI-BULK



The unsterile version of the Henke-Ject® 2-Part, All-Plastic syringes are manufactured to the high standards expected from Henke-Sass, Wolf GmbH. These syringes are very popular in a wide range of animal health and laboratory applications due to their inert components. Manufactured from polypropylene and polyethylene only, these syringes do not contain rubber or silicone oil. Available in both "Mini-Bulk" laboratory convenience packages, as well as large bulk cases. 1mL Luer Lock and all 50mL syringes are Air-Tite Premium brand.

UNSTERILE SYRINGES ARE MADE AT THE SAME
FACTORY AS STERILE, BUT ARE NOT FDA REGISTERED



UNSTERILE "MINI-BULK" 2-PART SYRINGES


Catalog Number	Specifications	Tip Type	Sterility	Quantity per Bag	Quantity per Case
ABC1	1mL	Luer Slip	Unsterile	100	7000
ABC3LS	3mL Centric	Luer Slip	Unsterile	100	6000
ABC5LS	5mL Centric	Luer Slip	Unsterile	100	3200
ABC10LS	10mL Centric	Luer Slip	Unsterile	100	1900
ABC20LS	20mL Centric	Luer Slip	Unsterile	100	1000
AMB50LS	50mL Eccentric	Luer Slip	Unsterile	30	360

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Bag	Quantity per Case
AMB1LL	1mL	Luer Lock	Unsterile	100	4200
ABC3LL	3mL	Luer Lock	Unsterile	100	6000
ABC5LL	5mL	Luer Lock	Unsterile	100	3200
ABC10LL	10mL	Luer Lock	Unsterile	100	1900
ABC20LL	20mL	Luer Lock	Unsterile	50	1000
AMB50LL	50mL	Luer Lock	Unsterile	30	360

UNSTERILE LARGE BULK 2-PART SYRINGES

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Case
AB1	1mL	Luer Slip	Unsterile	7000
AB3LS	3mL Centric	Luer Slip	Unsterile	6300
AB5LS	5mL Centric	Luer Slip	Unsterile	3600
AB10LS	10mL Eccentric	Luer Slip	Unsterile	2000
AB20LS	20mL Eccentric	Luer Slip	Unsterile	1000

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Case
AB3LL	3mL	Luer Lock	Unsterile	6300
AB5LL	5mL	Luer Lock	Unsterile	3600
AB10LL	10mL	Luer Lock	Unsterile	2000
AB20LL	20mL	Luer Lock	Unsterile	1000



EXEL INTERNATIONAL® BULK UNSTERILE SYRINGES



EXEL International's bulk, unsterile syringes offer convenience for those who do not need a sterile product for their application. Offering a cost savings in addition to the convenience, EXEL bulk, unsterile syringes may be just the solution you are seeking for your process.

UNSTERILE SYRINGES ARE MADE AT THE SAME FACTORY AS STERILE, BUT ARE NOT FDA REGISTERED



UNSTERILE EXEL SYRINGES

Air-Tite
Products Co., Inc.

Unsterile Syringes

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Case
EB1LS-4	1mL	Luer Slip	Unsterile	4000
EB3LS	3mL Centric	Luer Slip	Unsterile	1500
EB5LS	5mL Centric	Luer Slip	Unsterile	2000
EB10LS	10mL Centric	Luer Slip	Unsterile	1500
EB20LS	20mL Centric	Luer Slip	Unsterile	1000
EB30LS	30mL Eccentric	Luer Slip	Unsterile	400
EB60LS	60mL Eccentric	Luer Slip	Unsterile	400

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Case
EB1LL	1mL	Luer Lock	Unsterile	4000
EB3LL	3mL	Luer Lock	Unsterile	1500
EB5LL	5mL	Luer Lock	Unsterile	2000
EB10LL	10mL	Luer Lock	Unsterile	1500
EB20LL	20mL	Luer Lock	Unsterile	1000
EB30LL	30mL	Luer Lock	Unsterile	400
EB60LL	60mL	Luer Lock	Unsterile	400

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Case
EB30C	30mL	Catheter	Unsterile	400
EB50C	50mL	Catheter	Unsterile	400



B BRAUN EXADORAL ORAL DOSE SYRINGES®

B | BRAUN

These 2-part oral syringes by B Braun are designed for upside down withdrawal of liquid medication. The easy handling, syringe-type design allows for exact, safe and simple withdrawal and administration of oral liquid drugs. Made from polypropylene and polyethylene with a centric oral cone, these syringes have a highly transparent barrel with contrasting black graduations and purple plungers for easy identification. The safe plunger backstop helps avoid accidental plunger removal and the silicone-oil free production process leads to an inert solution.






B BRAUN EXADORAL® ORAL DOSE SYRINGES

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Case
BBORAL1	1mL	Oral	Unsterile	7000
BBORAL2	2(3)mL	Oral	Unsterile	6000
BBORAL5	5mL	Oral	Unsterile	3000
BBORAL10	10mL	Oral	Unsterile	2000
BBORAL20	20mL	Oral	Unsterile	1000

Oral Syringes



REPEATER SYRINGES



Eco-Matic®

Durable, self-filling syringe for injections made from high strength plastic, featuring a strong metal Luer Lock needle fitting. Available in a tubing version which uses a draw off, or a bottle mount version for secure attachment of the medication bottle (50-100ml) directly on the syringe. Ergonomic, compact V-grip design for ease of use, especially in high volumes. Accurate, continuously adjustable dose setting. Removable barrel for cleaning.

SIZE	PART#
0.3mL Bottle Mount Type	PS51B
1mL Tubing Type	PS53
1mL Bottle Mount Type	PS53B
2mL Tubing Type	PS52
2mL Bottle Mount Type	PS52B
5mL Tubing Type	PS55
5mL Bottle Mount Type	PS55B



Multi-Matic®

Repeater syringe for multiple dosing of livestock. Durable, lightweight pressure-cast aluminum handle. Durable medical grade polycarbonate barrel with metal Luer Lock needle fitting molded to the end. Biocompatible. Fully autoclavable (134°C / 273°F). Chemically resistant to most medications. UV-resistant and break-resistant. Exceptionally accurate dosing (Tolerance < 1%). Easy dose setting via 0.5ml (25 ml) or 1ml (50 ml) dosage wheel. Ergonomic design for ease of use, especially for high volume uses

SIZE	PART#
25mL	PS20
50mL	PS50



Vet-Matic®

High quality, automatic self-filling tube type syringe for injections. Metal/glass version. Durable piston rod. Easy to assemble/dismantle for servicing. Accurate, continuously adjustable dose setting. Metal Luer Lock needle fitting.

SIZE	PART#
10mL	PS100



Uni-Matic®

Easy-to-handle self-filling syringe which is ideal for poultry, swine, and small animal injections. Available in a tubing version which uses a draw off, or a bottle mount version for secure attachment of the medication bottle (50-100ml) directly on the syringe. Glass barrel for easy disassembly and cleaning. Quality stainless steel valve system for trouble-free dosing. Pin lock system: Bottle holder or tubing attachment can be positioned as required. Durable and precise steel plunger (no O-ring). Innovative, ergonomic design provides fatigue-free function. Medication bottles can be replaced quickly and easily with new bottle adapter. Exceptionally accurate dosing. Bottle adapter interchangeable with Ferro-Matic®

SIZE	PART#
1mL Tubing Type	PS91
2mL Tubing Type	PS92
2mL Bottle Mount Type	PS92B
5mL Tubing Type	PS95



Ferro-Matic®

High quality, automatic self-filling tube type syringe for injections. Metal/glass version. Durable piston rod. Easy to assemble/dismantle for servicing. Accurate, continuously adjustable dose setting.

SIZE	PART#
5mL	PS98



Repeater Syringe Replacement Parts Available - Contact Us For Additional Information.

REPEATER & MULTI-DOSE SYRINGES



Drench-Matic®

Automatic self-filling syringe for the administration of solutions and suspensions for oral and pour-on use. Pressure cast aluminum handle. High strength plastic barrel. Easy to use and clean. Dosing accuracy (Tolerance < 1% of dose). Dosage is continuously adjustable. Drenching pressure is smoothly adjustable. Automatic lubrication due to twin O-ring design.

SIZE	PART#
10mL	PS28
30mL	PS30



Henke-Ject TBC®

Dosage syringe for tuberculin testing. Fully visible glass barrel for visual verification. Ergonomic light weight metal, plastic design. 2ml barrel with 0.1ml increment trigger. Safe and easy needle change out via Luer Lock connection.

SIZE	PART#
2mL	PSD2.0TBC



Ideal Mega-Shot™ Vaccination System

Ideal's MegaShot Pistol-Grip Syringe repetitively injects doses of 1 to 5 cc quickly and easily without sacrificing accuracy, as its rear dosage dial locks securely to insure a precise dose. Its frame and plunger disassemble easily for cleaning and sterilization. MegaShot features a rubber plunger and virtually indestructible Lexan® barrel.

SIZE	PART#
50mL	ID1000



Ideal Pro-Shot™ Vaccination System

50ml pistol grip syringe. Easily adjustable 1ml to 5ml. Repetitive dosing. Durable polypropylene barrel.


SIZE	PART#
50mL Plastic Handle	ID1008
50mL Aluminum Handle	ID1009



Ideal Nylon Syringes

Heavy Duty field use syringe. Easy to disassemble and clean. Popular farm use syringe that is durable and great for all animals. Adjustable dosing nut allows for accurate dosage in 1ml increments.

SIZE	PART#
10mL (Blue)	ID9810
20mL (Red)	ID9812



Metal Syringes

Heavy Duty field use syringe. Easy to disassemble and clean. Popular farm use syringe that is durable and great for all animals. Adjustable dosing nut allows for accurate dosage in 1ml increments.

SIZE	PART#
50mL	PS56




Repeater Syringes

AIR-TITE BRAND SYRINGES WITH NEEDLES

Air-Tite
Products

Air-Tite's line of Veterinary and Laboratory syringes with attached needle are quality products at budget-achieving costs. Ideally suited for general applications requiring a non-human-medical solution. Custom solutions are possible with this line - contact us for more details.

NOT FDA REGISTERED FOR HUMAN USE



AIR-TITE LUER SLIP SYRINGES WITH MOUNTED NEEDLE

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
MS12558	1mL Tb	25g x 5/8"	Sterile	100	1000
MS12712	1mL	27g x 1/2"	Sterile	100	1000
MS3181	3mL	18g x 1"	Sterile	100	1000
MS3201	3mL	20g x 1"	Sterile	100	1000
MS32234	3mL	22g x 3/4"	Sterile	100	1000
MS3221	3mL	22g x 1"	Sterile	100	1000
MS32558	3mL	25g x 5/8"	Sterile	100	1000
MS5181	5mL	18g x 1"	Sterile	100	1000



AIR-TITE LUER LOCK SYRINGES WITH MOUNTED NEEDLE

Catalog Number	Specifications	Needle	Sterility	Quantity per Box	Quantity per Case
ML12558	1mL	25g x 5/8"	Sterile	100	1000
ML3181	3mL	18g x 1"	Sterile	100	1000
ML3201	3mL	20g x 1"	Sterile	100	1000
ML320112	3mL	20g x 1-1/2"	Sterile	100	1000
ML32234	3mL	22g x 3/4"	Sterile	100	1000
ML3221	3mL	22g x 1"	Sterile	100	1000
ML32558	3mL	25g x 5/8"	Sterile	100	1000
ML5181	5mL	18g x 1"	Sterile	100	1000
ML520112	5mL	20g x 1-1/2"	Sterile	100	1000
ML1018112	10mL	18g x 1-1/2"	Sterile	100	1000
ML2018112	20mL	18g x 1-1/2"	Sterile	50	500



Syringes with Needles

HENKE-JECT® SYRINGES WITH NEEDLES



Henke-Ject® (Formerly "Soft-Ject®") fixed needle syringes offer a crystal clear cylinder with easy to read scale. The piston, piston-rod, and cylinder consist of polypropylene. The fixed needle syringes are free of pyrogens, PVC, and latex, and are non-toxic. The 1cc syringes are sterile, individually blister packaged.



HENKE-JECT FIXED NEEDLE SYRINGES

Catalog Number	Specifications	Needle	Sterility	Quantity per Box	Quantity per Case
RU13012	1cc Tuberculin	30g x 1/2"	Sterile	100	2000

The By-packs from Henke-Sass, Wolf offer a convenient syringe with needle combination packaged side-by-side in individual strip packs.



HENKE-JECT SYRINGES / NEEDLE BY-PACKS

Catalog Number	Specifications	Needle	Sterility	Quantity per Box	Quantity per Case
R22234	2.5mL LS	22g x 3/4"	Sterile	100	2500
R2221	2.5mL LS	22g x 1"	Sterile	100	2500
R22558	2.5mL LS	25g x 5/8"	Sterile	100	2500
R2018	20mL LS	18g x 1-1/2"	Sterile	100	800



BD SYRINGES WITH NEEDLES




BD INSULIN SYRINGES WITH FIXED NEEDLES

Catalog Number	Specifications	Needle	Sterility	Quantity per Box	Quantity per Case
BD328438	3/10mL	31g x 5/16"	Sterile	100	500
BD328468	1/2mL	31g x 5/16"	Sterile	100	500
BD329424	1mL	28g x 1/2"	Sterile	100	500
BD328418	1mL	31g x 5/16"	Sterile	100	500



BD SYRINGES WITH MOUNTED NEEDLES

Catalog Number	Specifications	Needle	Sterility	Quantity per Box	Quantity per Case
BD309623	1mL TB	27g x 1/2"	Sterile	100	800
BD309637	1mL LL	20g x 1"	Sterile	100	800
BD309580	3mL LL	18g x 1-1/2"	Sterile	100	800
BD309570	3mL LL	25g x 5/8"	Sterile	100	800



Syringes with Needles

EXEL INTERNATIONAL® SYRINGES WITH NEEDLES



EXEL's Comfort Point fixed needle syringes offer a permanently attached needle, zero dead space, and smooth plunger action for greater filling accuracy. The auto-ground needle tip offers greater comfort. These syringes are sterile, non-toxic, and non-pyrogenic. All EXEL fixed needle syringes are packaged 10 per bag and 100 per box.

EXEL syringes with mounted needle offer convenience at an economic price. A wide selection of luer slip and luer lock syringes with standard size needles pre-attached. These syringes allow the user to reduce the number of inventory items and there is no need to open multiple packs to assemble a ready-to-use syringe.

EXEL VET FIXED NEEDLE INSULIN SYRINGES (U-40 GRADUATIONS)

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
EUV052913	0.5cc Insulin	29g x 1/2"	Sterile	100	500
EUV12913	1cc Insulin	29g x 1/2"	Sterile	100	500

EXEL FIXED NEEDLE INSULIN SYRINGES (U-100 GRADUATIONS)


Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
EU032913	0.3cc Insulin	29g x 1/2"	Sterile	100	500
EU03308	0.3cc Insulin	30g x 5/16"	Sterile	100	500
EU03318	0.3cc Insulin	31g x 5/16"	Sterile	100	500
EU052813	0.5cc Insulin	28g x 1/2"	Sterile	100	500
EU052913	0.5cc Insulin	29g x 1/2"	Sterile	100	500
EU05308	0.5cc Insulin	30g x 5/16"	Sterile	100	500
EU05318	0.5cc Insulin	31g x 5/16"	Sterile	100	500
EU1211	1cc Tuberculin	21g x 5/8"	Sterile	100	500
EU12558	1cc Tuberculin	25g x 5/8"	Sterile	100	500
EU12813	1cc Insulin	28g x 1/2"	Sterile	100	500
EU12913	1cc Insulin	29g x 1/2"	Sterile	100	500
EU1308	1cc Insulin	30g x 5/16"	Sterile	100	500
EU1318	1cc Insulin	31g x 5/16"	Sterile	100	500

EXEL LUER SLIP SYRINGES WITH MOUNTED NEEDLE

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
ES12558	1cc	25g x 5/8"	Sterile	100	1000
ES1251	1cc	25g x 1"	Sterile	100	1000
ES12612	1cc	26g x 1/2"	Sterile	100	1000
ES12712	1cc	27g x 1/2"	Sterile	100	1000
ES3201	3cc	20g x 1"	Sterile	100	1000
ES3211	3cc	21g x 1"	Sterile	100	1000
ES32234	3cc	22g x 3/4"	Sterile	100	1000
ES3221	3cc	22g x 1"	Sterile	100	1000
ES3231	3cc	23g x 1"	Sterile	100	1000
ES32558	3cc	25g x 5/8"	Sterile	100	1000

EXEL LUER LOCK SYRINGES WITH MOUNTED NEEDLE

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
EL318112	3cc	18g x 1-1/2"	Sterile	100	1000
EL3201	3cc	20g x 1"	Sterile	100	1000
EL320112	3cc	20g x 1-1/2"	Sterile	100	1000
EL3211	3cc	21g x 1"	Sterile	100	1000
EL321112	3cc	21g x 1-1/2"	Sterile	100	1000
EL32234	3cc	22g x 3/4"	Sterile	100	1000
EL3221	3cc	22g x 1"	Sterile	100	1000
EL322114	3cc	22g x 1-1/4"	Sterile	100	1000
EL322112	3cc	22g x 1-1/2"	Sterile	100	1000
EL3231	3cc	23g x 1"	Sterile	100	1000
EL323112	3cc	23g x 1-1/2"	Sterile	100	1000
EL32558	3cc	25g x 5/8"	Sterile	100	1000
EL3251	3cc	25g x 1"	Sterile	100	1000
EL325112	3cc	25g x 1-1/2"	Sterile	100	1000
EL327114	3cc	27g x 1-1/4"	Sterile	100	1000
EL5201	5cc	20g x 1"	Sterile	100	800
EL5211	5cc	21g x 1"	Sterile	100	800
EL521112	5cc	21g x 1-1/2"	Sterile	100	800
EL5221	5cc	22g x 1"	Sterile	100	800
EL522112	5cc	22g x 1-1/2"	Sterile	100	800
EL10181	10cc	18g x 1"	Sterile	100	800
EL10201	10cc	20g x 1"	Sterile	100	800
EL10211	10cc	21g x 1"	Sterile	100	800
EL1021112	10cc	21g x 1-1/2"	Sterile	100	800
EL10221	10cc	22g x 1"	Sterile	100	800
EL1022112	10cc	22g x 1-1/2"	Sterile	100	800



Syringes with Needles

AIR-TITE BRAND VETERINARY NEEDLES





Air-Tite
Products

Air-Tite Products' line of veterinary hypodermic needles offers one of the broadest selection of needles on the market. "Hard-to-find" sizes are available for nearly any non-human application, and custom sizes may be available upon request. Manufactured with a polypropylene hub and stainless steel to our specification, these regular wall needles will stand up to most applications. If you don't see the size you are seeking, please contact our customer service department for details on unique solutions.



AIR-TITE BRAND VETERINARY NEEDLES

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
N14112	14	1-1/2"	Sterile	100	10,000
N142	14	2"	Sterile	100	10,000
N14314	14	3-1/4"	Sterile	100	5,000
N1634	16	3/4"	Sterile	100	1,000
N161	16	1"	Sterile	100	1,000
N16112	16	1-1/2"	Sterile	100	1,000
N163D	16	3"	Sterile	100	5,000
N164	16	4"	Sterile	100	1,000
N165	16	5"	Sterile	25	500
N1814	18	1/4"	Sterile	100	4,000
N1812	18	1/2"	Sterile	100	1,000
N1834	18	3/4"	Sterile	100	10,000
N181	18	1"	Sterile	100	1,000
N18112	18	1-1/2"	Sterile	100	1,000
N182	18	2"	Sterile	100	1,000
N183	18	3"	Sterile	100	5,000
N19112NS*	19	1-1/2"	Sterile	100	5,000
N193D	19	3"	Sterile	100	5,000
N2014	20	1/4"	Sterile	100	4,000
N2038	20	3/8"	Sterile	100	1,000
N2012	20	1/2"	Sterile	100	1,000
N2034	20	3/4"	Sterile	100	10,000
N201	20	1"	Sterile	100	1,000
N20112	20	1-1/2"	Sterile	100	1,000
N211	21	1"	Sterile	100	1,000
N21112	21	1-1/2"	Sterile	100	1,000
N2238	22	3/8"	Sterile	100	1,000
N2234	22	3/4"	Sterile	100	1,000
N221	22	1"	Sterile	100	1,000
N22112	22	1-1/2"	Sterile	100	1,000
N223	22	3"	Sterile	100	1,000
N22312	22	3-1/2"	Sterile	100	1,000
N224	22	4"	Sterile	100	1,000
N231	23	1"	Sterile	100	10,000
N2558	25	5/8"	Sterile	100	1,000
N251	25	1"	Sterile	100	10,000
N25112	25	1-1/2"	Sterile	100	10,000
N2612	26	1/2"	Sterile	100	1,000
N2712	27	1/2"	Sterile	100	10,000
N276	27	6mm	Sterile	100	5,000



*N19112NS is non-siliconized

EXEL INTERNATIONAL® HYPODERMIC NEEDLES



EXEL International® offers a variety of hypodermic needles with regular wall thickness, made of 304 surgical stainless steel, non-toxic and pyrogen-free. Each box is color coded to help the clinician easily identify the appropriate gauge. Each needle is individually packaged, sterilized, and ready to use.

EXEL INTERNATIONAL® HYPODERMIC NEEDLES

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
EN181	18	1"	Sterile	100	2,000
EN18112	18	1-1/2"	Sterile	100	2,000
EN191	19	1"	Sterile	100	2,000
EN19112	19	1-1/2"	Sterile	100	2,000
EN201	20	1"	Sterile	100	2,000
EN20112	20	1-1/2"	Sterile	100	2,000
EN211	21	1"	Sterile	100	2,000
EN21112	21	1-1/2"	Sterile	100	2,000
EN2234	22	3/4"	Sterile	100	2,000
EN221	22	1"	Sterile	100	2,000
EN22112	22	1-1/2"	Sterile	100	2,000
EN2334	23	3/4"	Sterile	100	2,000
EN231	23	1"	Sterile	100	2,000
EN23112	23	1-1/2"	Sterile	100	2,000
EN2558	25	5/8"	Sterile	100	2,000
EN251	25	1"	Sterile	100	2,000
EN25112	25	1-1/2"	Sterile	100	2,000
EN2612	26	1/2"	Sterile	100	2,000
EN2712	27	1/2"	Sterile	100	2,000
EN27114	27	1-1/4"	Sterile	100	2,000
EN27112	27	1-1/2"	Sterile	100	2,000
EN2834	28	3/4"	Sterile	100	1,000
EN3012	30	1/2"	Sterile	100	2,000
EN301	30	1"	Sterile	100	2,000
EN3212	32	1/2"	Sterile	100	1,000




SecureTouch needles feature one-handed, finger-stroke activation. No solid surface required. An audible click as well as tactile feedback and visual indicator confirm safety mechanism engagement. Needle remains locked inside activated cover. Fits any standard brand syringe.



EXEL® SECURETOUCH SAFETY NEEDLES

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
ENS18112	18	1-1/2"	Sterile	100	1,000
ENS201	20	1"	Sterile	100	1,000
ENS231	23	1"	Sterile	100	1,000
ENS2558	25	5/8"	Sterile	100	1,000
ENS251	25	1"	Sterile	100	1,000
ENS25112	25	1-1/2"	Sterile	100	1,000
ENS27114	27	1-1/4"	Sterile	100	1,000
ENS27112	27	1-1/2"	Sterile	100	1,000
ENS3012	30	1/2"	Sterile	100	1,000



Needles

EZ FLO HYPODERMIC NEEDLES



Smooth Entry • Smooth Delivery

Air-Tite Product's newest addition, the EZ FLO line of hypodermic needles, offers high quality and a broad selection at an affordable price. Ranging from 18 gauge to 30 gauge, these needles offer a standard bevel and regular-wall to combine comfort and durability. Each needle is sterile, individually packaged in a box of 100. Made of surgical stainless steel and a polypropylene hub. 18g to 25g are thin wall. 27g and 30g are regular wall thickness.





EZ FLO BY AIR-TITE

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NEZ18034	18	3/4"	Sterile	100	2000
NEZ18100	18	1"	Sterile	100	2000
NEZ18112	18	1-1/2"	Sterile	100	2000
NEZ18200	18	2"	Sterile	100	2000
NEZ19100	19	1"	Sterile	100	2000
NEZ19112	19	1-1/2"	Sterile	100	2000
NEZ20100	20	1"	Sterile	100	2000
NEZ20112	20	1-1/2"	Sterile	100	2000
NEZ21100	21	1"	Sterile	100	2000
NEZ21112	21	1-1/2"	Sterile	100	2000
NEZ21200	21	2"	Sterile	100	2000
NEZ22034	22	3/4"	Sterile	100	2000
NEZ22100	22	1"	Sterile	100	2000
NEZ22112	22	1-1/2"	Sterile	100	2000
NEZ22200	22	2"	Sterile	100	2000
NEZ22300	22	3"	Sterile	100	2000
NEZ22400	22	4"	Sterile	100	2000
NEZ23034	23	3/4"	Sterile	100	2000
NEZ23100	23	1"	Sterile	100	2000
NEZ23112	23	1-1/2"	Sterile	100	2000
NEZ25058	25	5/8"	Sterile	100	2000
NEZ25100	25	1"	Sterile	100	2000
NEZ25112	25	1-1/2"	Sterile	100	2000
NEZ25200	25	2"	Sterile	100	2000
NEZ27100	27	1"	Sterile	100	2000
NEZ27114	27	1-1/4"	Sterile	100	2000
NEZ27112	27	1-1/2"	Sterile	100	2000
NEZ27200	27	2"	Sterile	100	2000
NEZ30012	30	1/2"	Sterile	100	2000
NEZ30100	30	1"	Sterile	100	2000

Needles

HENKE-JECT® HYPODERMIC NEEDLES




Henke-Ject® (formerly Fine-Ject®) hypodermic needles by Henke-Sass, Wolf GmbH are high quality yet economical. Made of surgical stainless steel, these needles are sterile, disposable, individually wrapped, and available in thin wall up to 25 gauge.





HENKE-JECT® HYPODERMIC NEEDLES BY HENKE-SASS, WOLF GMBH

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box	Quantity per Case
NH15112	15	1-1/2"	Thin	Sterile	100	4,000
NH172	17	2"	Thin	Sterile	100	3,000
NH18112	18	1-1/2"	Thin	Sterile	100	4,000
NH182	18	2"	Thin	Sterile	100	3,000
NH201	20	1"	Thin	Sterile	100	4,000
NH20112	20	1-1/2"	Thin	Sterile	100	4,000
NH212	21	2"	Thin	Sterile	100	3,000
NH221	22	1"	Thin	Sterile	100	4,000
NH22112	22	1-1/2"	Thin	Sterile	100	4,000
NH222	22	2"	Thin	Sterile	100	3,000
NH22234	22	2-3/4"	Thin	Sterile	100	3,000
NH224	22	4"	Thin	Sterile	100	2,000
NH2558	25	5/8"	Thin	Sterile	100	4,000
NH251	25	1"	Thin	Sterile	100	4,000
NH252	25	2"	Thin	Sterile	100	3,000
NH272	27	2"	Regular	Sterile	100	3,000
NH3012	30	1/2"	Regular	Sterile	100	4,000



Needles

TSK STERiJECT™ HYPODERMIC NEEDLES

TSK

TSK STERiJECT™ hypodermic needles are produced with a focus on needle sharpness and quality. All needles are thin-wall or ultra-thin wall. Made in Japan from one of the global leaders for ultra-fine tip needles, TSK, these needles are commonly used for aesthetic procedures, mesotherapy, and ophthalmology among other uses. Available in hard-to-find sizes for the specialist. 100 needles per box, surgical stainless steel, and poly hub.





TSK STERIJECT™ NEEDLES

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box	Quantity per Case
TSK2513U	25	13mm	Ultra-Thin	Sterile	100	6000
TSK2713U	27	13mm	Ultra-Thin	Sterile	100	6000
TSK2719U	27	19mm	Ultra-Thin	Sterile	100	6000
TSK2738U	27	38mm	Ultra-Thin	Sterile	100	6000
TSK304U	30	4mm	Ultra-Thin	Sterile	100	6000
TSK306U	30	6mm	Ultra-Thin	Sterile	100	6000
TSK3013U	30	13mm	Ultra-Thin	Sterile	100	6000
TSK3113	31	13mm	Thin	Sterile	100	6000
TSK324	32	4mm	Thin	Sterile	100	6000
TSK3206	32	6mm	Thin	Sterile	100	6000
TSK3209C	32	9mm	Thin	Sterile	100	6000
TSK3213	32	13mm	Thin	Sterile	100	6000
TSK3304	33	4mm	Thin	Sterile	100	6000
TSK3313	33	13mm	Thin	Sterile	100	6000

LOW DEAD SPACE NEEDLES

TSK3009LDS	30	9mm	Ultra-Thin	Sterile	100	6000
TSK3013LDS	30	13mm	Ultra-Thin	Sterile	100	6000
TSK3309LDS	33	9mm	Thin	Sterile	100	6000
TSK3313LDS	33	13mm	Thin	Sterile	100	6000
TSK02009LDS	0.2mm	9mm	Thin	Sterile	100	6000



Needles

METAL HUB NEEDLES

Ideal Aluminum Hub Needles

Hypodermic Veterinary Needles Regular Wall

Compares to Monoject®. Superior strength and durability. Strong, regular wall #304 stainless tubing. Sterile plastic capsule.



Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NA14112	14	1-1/2"	Sterile	100	1000
NA1658	16	5/8"	Sterile	100	1000
NA161	16	1"	Sterile	100	1000
NA16112	16	1-1/2"	Sterile	100	1000
NA1858	18	5/8"	Sterile	100	1000
NA1834	18	3/4"	Sterile	100	1000
NA181	18	1"	Sterile	100	1000
NA18112	18	1-1/2"	Sterile	100	1000



HSW ECO Metal Hub Needles

Hypodermic Veterinary Needles Regular Wall

Strong reusable needles offered in unsterile dozens, perfect for repetitive use in Hogs, Cattle and Poultry. Packaged in a plastic 12 pack.

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NHM1838	18	3/8"	Unsterile	12	120
NHM1858	18	5/8"	Unsterile	12	120
NHM181	18	1"	Unsterile	12	120
NHM2038	20	3/8"	Unsterile	12	120
NHM2058	20	5/8"	Unsterile	12	120
NHM2218	22	1/8"	Unsterile	12	120
NHM2238	22	3/8"	Unsterile	12	120



Ideal Udder Infusion Cannulas

Ideal udder infusion cannulas allow you to apply medication (such as Mastitis ointment) directly into the teat with ease and without the risk of contamination. Developed with the health and safety of both the animal and the operator in mind, these cannulas are essential to every dairy operation. Made of nickel-finished brass with a square needle hub, these udder infusion cannulas are available in two sizes and are built to last.

Catalog Number	Length	Sterility	Quantity per Box
ID1367C	2-1/2"	Unsterile	3
ID1369C	3-3/4"	Unsterile	3



SPINAL NEEDLES



EXEL Spinal Needles

Quincke bevel spinal needles. Siliconized for smooth, comfortable penetration. Transparent hub. Latex free. Convenient blister packaging.


Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NE183112A	18	3-1/2"	Sterile	50	100
NE203112A	20	3-1/2"	Sterile	50	100
NE22234A	22	2-3/4"	Sterile	50	100
NE223112A	22	3-1/2"	Sterile	50	100
NE253112A	25	3-1/2"	Sterile	50	100



BD Spinal Needles

The key/slot arrangement of stylet and cannula hubs on the BD® Quincke spinal needles facilitates proper needle bevel orientation. Fitted stylet designed to reduce tissue coring. Unique translucent window hub featuring contact clarity that helps allow visualization of CSF.

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
BD408360	18	6"	Sterile	10	50
BD405182	20	3-1/2"	Sterile	25	100
BD405161	22	1-1/2"	Sterile	N/A	100
BD405181	22	3-1/2"	Sterile	25	100
BD405148	22	5"	Sterile	10	50
BD405149	22	7"	Sterile	10	50



Needles

BLOOD COLLECTING NEEDLES


Air-Tite Products

Air-Tite Vet Premium Blood Collecting Needles

Regular Wall

Produced in China from the highest quality regular wall surgical stainless steel tubing available, they meet rigorous quality control specifications and have been an industry standard for over 25 years. These double ended needles are for use with blood collecting needle holders and tubes. Packaged in individual sterile packs.

Catalog Number	Type	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NBC181	Single Use	18	1"	Sterile	100	5000
NBC181M	Multi-Type	18	1"	Sterile	100	5000
NBC18112	Single Use	18	1-1/2"	Sterile	100	5000
NBC18112M	Multi-Type	18	1-1/2"	Sterile	100	5000
NBC201	Multi-Type	20	1"	Sterile	100	5000
NBC20112	Multi-Type	20	1-1/2"	Sterile	100	5000



EXEL^{INT}®

EXEL Blood Collecting Needles

Multi-sample needles permit several samples to be taken with a single puncture. Sharp & smooth edges help make penetration painless and provides easy connection to rubber stoppers. These products are non-toxic and do not contain natural rubber latex.

Catalog Number	Type	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
EBC201	Multi-Type	20	1"	Sterile	100	1000
EBC20112	Multi-Type	20	1-1/2"	Sterile	100	1000
EBC211	Multi-Type	21	1"	Sterile	100	1000
EBC21112	Multi-Type	21	1-1/2"	Sterile	100	1000
EBC221	Multi-Type	22	1"	Sterile	100	1000
EBC22112	Multi-Type	22	1-1/2"	Sterile	100	1000



BLOOD COLLECTION TUBES, SETS, & ADAPTERS



Generic Vet Blood Collection Tubes

Catalog Number	Specs	Color	Additive	Quantity per Pack	Quantity per Case
MV3R	3mL Glass	Red Top	None (with Silicone)	100	1200
MV3P	3mL Plastic	Red Top	None (with Silicone)	100	1200
MV5PSSST	5mL Plastic	Gold Top	Gel & Clot Activator	100	1200
MV6P	6mL Plastic	Red Top	None (with Silicone)	100	1200



BD Vacutainer®

Catalog Number	Specs	Color	Additive	Quantity per Pack	Quantity per Case
BD366668	3mL Plastic	Red Top	Clot Activator (with Silicone)	100	1000
BD367812	4mL Plastic	Red Top	Clot Activator (with Silicone)	100	1000
BD367815	6mL Plastic	Red Top	Clot Activator (with Silicone)	100	1000
BD366430*	10mL Glass	Red Top	None (with Silicone)	N/A	1000
BD367985*	10mL Plastic	Red/Gray Top	Gel & Clot Activator	N/A	1000

*Special Order - Case Quantities Only. No Returns.



Blood Collection Holders & Adapters


Standard Holders for use with standard blood collection tubes and double sided needles. Unsterile.
Catalog Number: NBT1 Sold Individually

EXEL Multi-Sample Blood Collection Luer Adapter Connect this device to your butterfly or blood collection needle and insert the blood collection tube to draw. Single use, sterile.

Catalog Number: EBCA 50/box

EXEL Blood Collection Adapter for use with winged infusion sets.


Catalog Number: EW12BC 100/box



EXEL Blood Collection Sets

12" Tubing

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
EBS2134	21	3/4"	Sterile	50	200
EBS2334	23	3/4"	Sterile	50	200
EBS2534	25	3/4"	Sterile	50	200



IV SETS & WINGED INFUSION SETS



Intravenous Administration Sets

Non-toxic, non-pyrogenic fluid path. Vented/non-vented with roller clamp and occlusion clamp with Y-injection site.

Catalog Number	Drops/mL	Specs	Tubing	Sterility	Quantity per Box
EIV7815	15	Luer Slip	78"	Sterile	50
EIV7815LL	15	Luer Lock	78"	Sterile	50
EIV7815-2C	15	Lock/Slip With Dual Clamp	78"	Sterile	50
EIV10515	15	Multi-Site*	105"	Sterile	50
EIV8420F	20	Filtered**	84"	Sterile	50
EIV7860	60	Luer Slip	78"	Sterile	50
EIV7860LL	60	Luer Lock	78"	Sterile	50



*Two Y-injection sites. Priming volume 28mL/approximate inner diameter is 3.0mm and outer diameter is 4.1mm
** .2 Micron IV Filter

Winged Infusion Sets

These EXEL scalp vein sets have flexible wings to hold the needle in position and prevent inadvertent pullout. Thin stainless steel siliconized cannula for easy insertion. Sterilized with ethylene oxide and ready for immediate use. Latex free, individually packaged and color coded for easy identification. All sets have 12" tubing.

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
EWI1934	19	3/4"	Sterile	50	500
EWI2034	20	3/4"	Sterile	50	500
EWI2134	21	3/4"	Sterile	50	500
EWI2334	23	3/4"	Sterile	50	500
EWI2534	25	3/4"	Sterile	50	500
EWI2734	27	3/4"	Sterile	50	500



CATHETERS



Intravenous Catheter

Functional design case permits one-handed opening and no rolling. Color coded casing cap allows for easy identification. Translucent catheter hub and flashback chamber allows for easy detection of blood flashback. Precision finished radio opaque Teflon catheter assures stable flow and eliminates catheter kink during venipuncture

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
EIV14114	14	1-1/4"	Sterile	50	200
EIV142	14	2"	Sterile	50	200
EIV16114	16	1-1/4"	Sterile	50	200
EIV162	16	2"	Sterile	50	200
EIV18114	18	1-1/4"	Sterile	50	200
EIV182	18	2"	Sterile	50	200
EIV201	20	1"	Sterile	50	200
EIV20114	20	1-1/4"	Sterile	50	200
EIV202	20	2"	Sterile	50	200
EIV221	22	1"	Sterile	50	200
EIV2434	24	3/4"	Sterile	50	200



PNEUMATIC VACCINATOR



FRESHLY HATCHED: OUR NEW PNEUMATIC VACCINATOR FOR DAY-OLD CHICKS




- | | |
|----------------------|--|
| Efficacy | Adjustable trigger element for safe and consistent positioning of day-old chicks |
| Accurate | Highest dosing accuracy due to reliable double-stroke cylinders and high-quality syringes |
| Sophisticated | Perfectly shaped cover enables best injection angle for healthy day-old chicks |
| User-friendly | Easy change of injection volume with only one syringe due to a patented click system |
| Smart | Easy set-up, maintenance and cleaning procedures for long service life and highest biosecurity standards |
| Versatile | Injection of either one or two different vaccines with only one needle |
| Resistant | Robust and medical grade components allow use with nearly every vaccine |



PART #
VPVAC1
VPVAC2

SPECS
Single Vaccinator
Double Vaccinator



EPIG

NEEDLELESS VACCINATOR



THE FIRST ALL-IN-ONE HANDHELD NEEDLE-FREE DEVICE FOR INTRAMUSCULAR VACCINATION

QUALITY
MADE IN
GERMANY

- Efficacy** Optimum vaccination results thanks to direct feedback on each injection from multicolored LED bars on each side of the device
- Accurate** Highest dosing accuracy due to an integrated sensor control detecting and reporting air inside the liquid circuit
- Simple** No external power sources required - insert battery and go
- User-friendly** Easy set-up, maintenance and cleaning procedures for long service life and highest biosecurity standards
- Smart** Integrated Bluetooth enables the evaluation of vaccination data after use (EPIG® App required, available for iOS and Android)
- Resistant** Robust medical grade components ensure dust and water tightness (splash water resistant acc. IP54)
- Versatile** Dosage adjustable to 1ml and 2ml, vial sizes between 20 ml and 500 ml, compatible with most common intramuscular vaccines
- Sophisticated** Increasing vaccination quality and piglet welfare while reducing cross-contamination



Special Order Product
Call 1-800-231-7762 for more information.



To Order: (800) 231-7762 | atinfo@air-tite.com

MISCELLANEOUS SUPPLIES

Plastic Dispensing Tips

For medium viscosity fluid dispensing. Fits all standard luer slip and luer lock syringes. Non-sterile, clear, 100/bag.

Catalog Number	Specs
NDT-STR	Straight Tip
NDT-CRV	Curved Tip
NDT-CAP	Caps



Syringe Caps

Polypropylene syringe caps. Fit all luer slip and luer lock syringes.

Catalog Number	Color	Specs	Sterility	Quantity per Pack
BUC	Black	Translucent	Unsterile	1000
BUC/ADB	Blue	Opaque	Unsterile	1000
BUC/ADR	Red	Opaque	Unsterile	1000
BUCB	Light Blue	Translucent	Unsterile	1000
BUCC-V2	Clear	Translucent	Unsterile	1000
BUCP	Pink	Translucent	Unsterile	1000
BUCWP-100	White	Opaque, Premium Design	Unsterile	100
BUCWS	White	Opaque	Unsterile	1000
BUCWS-100	White	Opaque	Unsterile	100
BUCY	Yellow	Translucent	Unsterile	1000
STC	Blue	Translucent	Sterile	100



BUCWP-100



STC



BEMIS® SHARPS CONTAINERS

All sharps container units are designed for safe use, effectively blocking hand access from discarded sharps. Units snap together, permanently locking with tight seams between the bases and covers. Walls are puncture-resistant.

Catalog Number	Volume	Dimensions	Type	Specs
SC101R	1 Quart	4.25"L x 4.25"W x 6.5"H	Phlebotomy	Fits Phlebotomy trays
SC105R	5 Quart (1.25 Gallons)	13"L x 5.5"W x 9"H	Sentinel	Rotating Lid, Wall Mount or Cabinet
SC108R	8 Quart (2 Gallons)	11.75"L x 7.75"W x 8.75"H	SafeView	Table Top or Wall Mount
SC112R	12 Quart (3 Gallons)	15"L x 6.5"W x 15"H	SafeView	Standard Lid
SC212R	12 Quart (3 Gallons)	15"L x 6.5"W x 15"H		Rotating Lid
SC132R	8 Gallon	15.25"L x 11"W x 22"H		
SC144R	11 Gallon	15.25"L x 11"W x 22"H		



BEMIS® SHARPS CONTAINER ACCESSORIES

Catalog Number	Description	Fits
SC101H	1 Quart Holder	SC101R
SC100B	Mounting Bracket	SC105R, SC108R
SC100K	Bracket Key	SC100B
SC100W	3 Gallon Wire Bracket	SC112R, SC212R
SC5QWC	5 Quart Wall Cabinet	SC105R
SC3GWC	3 Gallon Wall Cabinet	SC212R
SCGH1	Glove Box Holder	Snaps onto SC5QWC & SC3GWC



UNIVERSAL BRACKETS (FOR MAIL-BACKS)

Universal mounting brackets fit 1 quart, 1 gallon, 2 gallon, and 3 gallon Mail-Back Systems (See below). Locking brackets include lock & key.

Catalog Number	Description
SCMSB123	Universal Bracket
SCMSB123L	Universal Bracket - Locking



SHARPS MAIL-BACK SYSTEM

Each Kit includes all components required for collection and transportation of sharps and regulated medical waste. Step-by-step instructions, Waste shipping containers, Pre-paid and pre-addressed shipping labels, Pre-paid disposal, Tracking documents. Approved by the U.S. Postal Service.

Catalog Number	Volume	Dimensions	Quantity	Specs
SCM1Q	1 Quart System	3.75"L x 3.5"W x 8"H	1 Container	Holds approximately 42 3cc syringes
SCM1Q-5	1 Quart System	3.75"L x 3.5"W x 8"H	5 Containers	Holds approximately 42 3cc syringes
SCM1G	1 Gallon System	9"L x 5.5"W x 8.5"H	1 Container	Holds approximately 98 3cc syringes
SCM1G-2	1 Gallon System	9"L x 5.5"W x 8.5"H	2 Containers	Holds approximately 98 3cc syringes
SCM2G	2 Gallon System	9"L x 5.5"W x 14.25"H	1 Container	Holds approximately 245 3cc syringes
SCM2G-2	2 Gallon System	9"L x 5.5"W x 14.25"H	2 Containers	Holds approximately 245 3cc syringes
SCM3G	3 Gallon System	9"L x 5.5"W x 19.75"H	1 Container	Holds approximately 358 3cc syringes



NON-SHARPS MEDICAL WASTE MAIL-BACK SYSTEM

Each mail-back system includes all components required for collection and transportation of hazardous medical waste. Price includes shipping to your location, container, packing material, shipping to the disposal site, prepaid disposal and tracking documents. IMPORTANT: Must NOT contain any LOOSE sharps. Sharps must be contained inside another FDA approved sharps container before being placed in one of these wide-mouth containers.

Catalog Number	Volume	Dimensions	Quantity
MWM1G	1 Gallon System	8.375" x Diameter 8"H	1 Container
SCM5G	5 Gallon System	10.45" x Diameter 14.51"H	1 Container



Mail-Back Systems: Not responsible for any applicable registration or licensing fees. Returns are subject to a 25% restocking fee.


To Order: (800) 231-7762 | atinfo@air-tite.com

MISCELLANEOUS SUPPLIES

Medi-Cut®

Sheffield stainless steel surgical blades with finely honed cutting edge guaranteed for sharpness and exacting performance. Sterile, individually foil packed.

100 blades per box

SPECS	PART#	10	11	12	15	20	22	24
#10	SB5010							
#11	SB5011							
#12	SB5012							
#15	SB5015							
#20	SB5020							
#22	SB5022							
#24	SB5024							

Miltex®

Surgical Blade Handles

SPECS	PART#
#3 (fits 10, 11, 12, 15)	MH4-7
#4 (fits 20, 21, 22, 23, 24)	MH4-8

Economy Handles

Surgical Blade Handles

SPECS	PART#
#3 (fits 10, 11, 12, 15)	SBH03
#4 (fits 20, 21, 22, 23)	SB4H

EXEL International®

Sterile Surgical Scalpels

Rigid plastic handles with stainless steel blades. 10/box

SPECS	PART#
#10	SDE10
#11	SDE11
#15	SDE15

Suture Removal Kits By Dynarex®

- Plastic tray provides sterile field
- Removable Tyvek lid
- Contents...Littauer scissors, 4" metal forceps, gauze sponge

PART#
PPG4521



Staple Remover

Disposable Skin Staple Remover

SPECS	PART#
Wide/Regular	NEW88008



Braided Synthetic Suture

Absorbable (Compares to Vicryl®)

PART#	SIZE	LENGTH	NEEDLE
SMV-J452	3-0	30"	FS-1
SMV-J453	2-0	30"	FS-1
SMV-J332	3-0	27"	CT-2
SMV-J333	2-0	30"	CT-2
SMV-J334	0	27"	CT-2
SMV-J451	4-0	27"	FS-1
SMV-J316	3-0	27"	SH

Monofilament Synthetic Suture

Absorbable (Compares to PDS II-Violet®)

PART#	SIZE	LENGTH	NEEDLE
SMV-Z397	4-0	30"	FS-2
SMV-Z398	3-0	30"	FS-2
SMV-Z451	2-0	30"	FS-1
SMV-Z452	3-0	30"	FS-1
SMV-Z467	0	30"	CP-1
SMV-Z333	2-0	30"	CT-2

Sterile Chromic Gut Suture

Absorbable

PART#	SIZE	LENGTH	NEEDLE
SMV-636	3-0	30"	FS-2
SMV-637	2-0	30"	FS-1
SMV-810	3-0	30"	CT-1

Nylon Monofilament Suture

Non Absorbable (Compares to Ethilon®)

PART#	SIZE	LENGTH	NEEDLE
SMV-662	4-0	30"	FS-2
SMV-663	3-0	30"	FS-1
SMV-664	2-0	30"	FS
SMV-666	1-0	30"	FS

Disposable Skin Stapler

Available in wide and regular widths.

35 staples per unit.

SPECS	PART#
35 wide (7.22mm x 4.3mm)	NEW88006
35 regular (5.4mm x 3.6mm)	NEW88005



Miltex® MeisterHand

MeisterHand® instruments are expertly crafted from surgical grade stainless steel German forgings and meet the very highest operating room quality standards and specifications. This instrument line is advantageously priced while offering great quality and value. We are a full line Miltex® dealer, so if you cannot find what you are looking for, contact us!

Economy Surgical Instruments

Economical instruments manufactured to meet high standards, guaranteed to be free from defects. Stainless steel, made in Pakistan.




Mayo Scissors

SPECS

- Economy 5-1/2" Straight
- Miltex®MH 5-1/2" Straight
- Economy 5-1/2" Curved
- Miltex®MH 5-1/2" Curved
- Economy 6-3/4" Straight
- Miltex®MH 6-3/4" Straight
- Economy 6-3/4" Curved
- Miltex®MH 6-3/4" Curved

PART#

- PRO24180
- MH5-120
- PRO24190
- MH5-122
- PRO24184
- MH5-124
- PRO24194
- MH5-126




Iris Scissors

SPECS

- Miltex®MH 4-1/2" Straight

PART#

- MH5-304



Lister Bandage Scissors

SPECS

- Economy 5-1/2"
- Miltex®MH 5-1/2"
- Economy 7-1/2"
- Miltex®MH 7-1/4"

PART#

- PRO38104
- MH5-514
- PRO38106
- MH5-516



Littauer Stitch Scissors

SPECS

- Economy 5"
- Miltex®MH 5-1/2"

PART#

- PRO24360
- MH9-104



Spencer Stitch Scissors

SPECS

- Miltex®MH 3-1/2"

PART#

- MH9-100




Metzenbaum Scissors

SPECS

- Economy 5-1/2" Straight
- Miltex®MH 5-1/2" Straight
- Economy 5-1/2" Curved
- Miltex®MH 5-1/2" Curved
- Economy 7" Straight
- Miltex®MH 7" Straight
- Economy 7" Curved
- Miltex®MH 7" Curved

PART#

- PRO24230
- MH5-179
- PRO24232
- MH5-180
- PRO24240
- MH5-181
- PRO24250
- MH5-182




Operating Scissors

SPECS

- Economy 4-1/2" Straight/Blunt
- Miltex®Vantage 4-1/2" Straight/Blunt

PART#

- PRO24099
- MH95-12V



MISCELLANEOUS SUPPLIES

TOWEL CLAMPS

Backhaus Towel Clamps

SPECS

Economy 3-1/2"
Miltex® MH 3-1/2"
Economy 5-1/2"
Miltex® MH 5-1/2"

PART#

PRO28912
MH7-504
PRO28914
MH7-506



Lorna Towel Clamps

SPECS

3-1/2"

PART#

PRO28915



NEEDLE HOLDERS


Olsen Hegar Needle Holders

SPECS

Economy 5-1/2"
Miltex® MH 5-1/2"
Economy 6-1/2"
Miltex® MH 6-1/2"
Economy 7-1/2"
Miltex® MH 7-1/4"

PART#

PRO36260
MH8-15
PRO36262
MH8-16
PRO36264
MH8-17




Mayo Hegar Needle Holders

SPECS

Economy 5-1/2"
Miltex® MH 5"
Economy 6" Satin Finish
Miltex® MH 6"
Economy 7" Satin Finish
Miltex® MH 7"

PART#

PRO36130
MH8-42
PRO36131S
MH8-44
PRO36132S
MH8-46



DRESSING & TISSUE FORCEPS

Thumb Dressing Forceps

SPECS

Economy 6" Serrated
Miltex® MH 6" Serrated

PART#

PRO26123
MH6-10



Adson Dressing Forceps

SPECS

4-3/4" Serrated

PART#

PRO26300



Thumb Tissue Forceps

SPECS

Economy 1x2 teeth 5-1/2"
Miltex® MH 1x2 teeth 5-1/2"
Economy 2x3 teeth 5-1/2"

PART#

PRO26142
MH6-44
PRO26152



Brown-Adson Tissue Forceps

SPECS

Miltex® MH 7x7 4-3/4"

PART#

MH6-124



Foerster Sponge Forceps



SPECS

	<u>PART#</u>
Economy 7" Straight	PRO30310
Economy 9-1/2" Straight	PRO30330
Miltex®MH 9-1/2" Straight	MH7-602

Rochester Pean Forceps



SPECS

	<u>PART#</u>
Economy 7-1/4" Curved	PRO28263
Miltex®MH 7-1/4" Curved	MH7-140
Economy 7-1/4" Straight	PRO28253
Miltex®MH 7-1/4" Straight	MH7-122
Economy 8" Curved	PRO28254

Hartman Alligator Forceps



SPECS

	<u>PART#</u>
Economy 3-1/2" Serrated	PRO50290
Economy 5-1/2" Serrated	PRO50291
Economy 6-1/2" Serrated	PRO50292

Spay & Neuter Pack

Products are quality German-made Hi-Tech Brand and carry a limited lifetime guarantee.

PART#

PPGSP1

Each kit contains:

- 1 - Mayo Scissors 6-3/4" Straight
- 1 - #3 Scalpel Handle
- 2 - Halstead Mosquito Forceps 5" Straight
- 2 - Halstead Mosquito Forceps 5" Curved
- 2 - Kelly Forceps 5-1/2" Straight
- 2 - Kelly Forceps 5-1/2" Curved
- 3 - Backhaus Towel Clamps 3-1/2"
- 1 - Olsen Hegar Needle Holder 6-1/2"
- 1 - Adson Brown Forceps 4-3/4" 9x9
- 1 - Allis Forceps 6" 4x5" (Canine Only)
- 1 - Snook Hook



Hartman Mosquito Forceps



SPECS

	<u>PART#</u>
Economy 3-1/2" Straight	PRO28120
Miltex®MH 3-1/2" Straight	MH7-24
Economy 3-1/2" Curved	PRO28122
Miltex®MH 3-1/2" Curved	MH7-26

Medicut Utility Shears

General use utility shears. Great for bandage removal, etc.

SPECS

	<u>PART#</u>
Economy 5-1/4"	PRO300
Economy 7-1/4"	PRO302



Kelly Hemostatic Forceps



SPECS

	<u>PART#</u>
Economy 5-1/2" Curved	PRO28162
Miltex®MH 5-1/2" Curved	MH7-38
Economy 5-1/2" Straight	PRO28160
Miltex®MH 5-1/2" Straight	MH7-36


Miltex® Vantage White Toenail Scissors



SPECS

	<u>PART#</u>
Miltex®Vantage - 4-1/2"	MH1718SSV

Rochester-Carmalt Forceps



SPECS

	<u>PART#</u>
Economy 6-1/4" Curved	PRO07-138

MISCELLANEOUS SUPPLIES

Heavy Duty Nippers

SPECS
5-1/4"

PART #
PRO881900



Nail Clippers Guillotine Style

SPECS
5-1/2"

PART #
PRO1911



Ideal V-Type Ear Notcher

Stainless Steel

SPECS
Small Notch

PART #
ID7005



Miltex® Spray Lube

Surgical Instrument Lubricant

Pre-mixed non silicone-base, anti-corrosive and non-greasy formula. Steam permeable and water soluble. Helps prevent spotting, staining and rusting. 8 oz. bottle

SPECS
Miltex® Spray Lube

PART #
M3-700



Miltex® Cleaner

Surgical Instrument Cleaner

Neutral pH & Phosphate-free formula. Rinses free, leaves no residue, and never stains instruments. Made for manual cleaning. 8 oz. bottle

SPECS
Miltex® Instrument Cleaner

PART #
M3-720



Miltex® Stain Removing Powder

Powder formula removes stains, tarnish and restores luster to surgical instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, etc. Ideal for use on stainless steel, chrome, nickel, brass, copper and silver.

SPECS
Miltex® Stain Remover - 3oz.

PART #
M3-740



Instrument Trays

Stainless steel trays suitable for autoclaving instrument packs.

SPECS
10" x 6-1/2" x 3/4" (standard)
19-1/8" x 12-5/8" x 3/4" (Mayo)

PART #
PPG49-365
PPG49-363



Miltex® Medicine Cup

Stainless Steel 2oz Cup

SPECS
Miltex® 2oz.

PART #
M3-914



Sterilization Pouches

Self-Seal Sterilization Pouches with Bio Indicator
200/box

SPECS - usable area	PART #
1-1/2" x 3-1/8" (2-1/4" x 4")	SP25
2-3/4" x 4-1/2" (3-1/2" x 5-1/4")	SP36
2-3/4" x 7-7/8" (3-1/2" x 10")	SP310
4-1/2" x 8-3/8" (5-1/4" x 10")	SP511
6-1/2" x 11-7/8" (7-1/2" x 13")	SP513



Cloth Instrument Wraps

- High quality, 2-ply, safety-surged
- Made of durable, colorfast, washable surgical sheeting

SPECS	PART #
12" x 12"	PPGW1212
18" x 18"	PPGW1818
24" x 24"	PPGW2424
30" x 30"	PPGW3030
36" x 36"	PPGW3636



C.S.R. Wraps

A superior non-woven wrapping fabric for gas and steam sterilization.

- Folds and drapes like muslin - no spring back.
- Strong - tear resistant - folds easily
- Repellent to water, blood and saline
- Effective against both air and water borne bacteria

SPECS	PART #	QTY
12" x 12"	PPG4412	1000
15" x 15"	PPG4415	1000
24" x 24"	PPG4424	500
36" x 36"	PPG4436	250



Disposable UnderPads

Highly absorbent, soft fluff fill. Waterproof, non-skid blue backing. Sealed on all 4 sides to prevent leakage. 50 pads/box

SPECS

23" x 36" Pads

PART #

PPG1343



Alcohol Prep Pads

Medium Size Individually Wrapped
200/box

PART #

PPA1103



Povidone Iodine Pads

Medium Size Individually Wrapped
100/box

PART #

PPA1108



Non-Adherent Pads, Sterile

- 1 pad/peel-open pouch, 100 pads/box
- Easily removable - will not stick to wound.
- Lets air in - speeds healing
- Highly absorbent
- Provides cushioning protection

SPECS

2" x 3"

PART #

PPG3423

3" x 4"

PPG3434



Sterile Gauze Sponges

8 Ply - 2 per pouch

SPECS

2" x 2"

PART #

PPG3322

4" x 4"

PPG3342



Cotton Rolls

1 lb non-sterile roll made from 100% cotton. Ideal for facial, nail, skin care applications, swabs, padding, leg wraps, etc. The highest quality available for applications where absorbency, flexibility, cleanliness, and softness are key requirements. 1 lb. bag.

PART #

PPG2661



Unsterile Gauze Sponges

12 Ply - 200 per bag

SPECS

2" x 2"

PART #

PPG3223

3" x 3"

PPG3233

4" x 4"

PPG3243



Cotton Balls

Medium Size - General Use
2000/box

PART #

PPG3170



Unsterile Stretch Gauze

White, 5 yards/roll, 12 rolls/bag.

PART #

PPG3102

PPG3103

PPG3104




Cotton Tipped Applicators

6" length, non-sterile with wooden sticks
100 per bag, 10 bags per box

PART #

PPG4302



X-Ray Detectable

Sterile Gauze Sponges

Multi-purpose gauze sponges, ideal for use in surgical procedures. Sterile, highly-absorbent, 16 ply. Fine mesh 100% gauze, virtually lint free. Latex free. Packed in 10-count plastic trays, 128 trays/case.

SPECS

4" x 4"

PART #

PPG3351



MISCELLANEOUS SUPPLIES

Co-Flex® Cohesive Bandage

5 yards per roll

Multicolored lightweight cohesive bandage; adheres to itself without clips or fasteners. Great for wound dressing, stabilizing catheters or situations requiring light compression. Multicolored packs contain Neon Pink, Blue, Purple, Light Blue, Neon Green, & Red.

PART #	SIZE	BOX
PPG3200 (Multi)	2"	36 rolls
PPG3200TN (Tan)	2"	36 rolls
PPG3200RD (Red)	2"	36 rolls
PPG3200GR (Green)	2"	36 rolls
PPG3300 (Multi)	3"	24 rolls
PPG3400 (Multi)	4"	18 rolls
PPG3400T (Tan)	4"	18 rolls



Sensi Wrap Self Adherent Bandage

5 yards per roll

A lightweight compression bandage which sticks to itself, but not to other materials or skin. Easily torn without scissors. Will not slip. Individually poly bagged. Tan or assorted color pack.

PART #	SIZE	BOX
PPG32SW (Color Pack)	2"	36 rolls
PPG32SWT (Tan)	2"	36 rolls
PPG33SW (Color Pack)	3"	24 rolls
PPG33SWT (Tan)	3"	24 rolls
PPG34SW (Color Pack)	4"	18 rolls
PPG34SWT (Tan)	4"	18 rolls



Paper Tape

10 yards per roll

Compares to 3M Micropore tape. Excellent breathability and comfort, lightweight and hypoallergenic. Holds when wet and won't harm delicate wounds.




SPECS	PART #	BOX
Paper Tape 1"	PPG3552	12 rolls
Paper Tape 2"	PPG3553	6 rolls
Paper Tape 3"	PPG3554	4 rolls

Cloth Tape

10 yards per roll

Compares to J & J Dermicell tape. Hypoallergenic cloth tape tears easily without shredding.

SPECS	PART #	BOX
Cloth Tape 1"	PPG3562	12 rolls
Cloth Tape 2"	PPG3563	6 rolls



Transparent Tape

10 yards per roll

Compares to 3M Transpore tape. Perforated for easy tearing - porous for greater breathability. Conforms & stretches to accommodate body movement.

SPECS	PART #	BOX
Transparent Tape 1"	PPG3572	12 rolls
Transparent Tape 2"	PPG3573	6 rolls



Cotton Elastic Adhesive Tape

Features a strong rubber-based adhesive to permit flexibility of joints without the loosening of the tape. Ideal for dynamic/pressure dressings requiring elasticity. Molds to body contours. Even compression and support.

SPECS	PART #	BOX
1" x 5 yards	PPGEA1	12 rolls
2" x 5 yards	PPGEA2	6 rolls
3" x 5 yards	PPGEA3	4 rolls



Autoclave Tape

Use with C.S.R. wraps and Cloth Instrument wraps. Indicates when sterilization is complete.

SPECS	PART #
1" x 60 yd.	PPGAUTO1



Latex Exam Gloves

Dynarex SafeTouch™ brand

3mil, 9" beaded cuff, ambidextrous.

Powder Free

SIZE	PART #	QUANTITY PER BOX
S	LS100S	100 Gloves
M	LS100M	100 Gloves
L	LS100L	100 Gloves
XL	LS100X	100 Gloves



Vinyl Exam Gloves

Dynarex SafeTouch™ brand

3mil, 9" beaded cuff, ambidextrous.

Powder Free

SIZE	PART #	QUANTITY PER BOX
S	VS100S	100 Gloves
M	VS100M	100 Gloves
L	VS100L	100 Gloves
XL	VS100X	100 Gloves



Latex Sterile

Procedure Gloves

DermAssist™ brand

5mil, 10" cuff, ambidextrous.

Sterilized by gamma radiation and packaged by the pair.

Powder Free

SIZE	PART #	QUANTITY PER BOX
S	LD200S	100 Gloves
M	LD200M	100 Gloves
L	LD200L	100 Gloves
XL	LD200X	100 Gloves



Latex High Risk Gloves

EP Blue® brand

15mil, 12" beaded cuff, ambidextrous.

Extra long, textured.

Powder Free

SIZE	PART #	QUANTITY PER BOX
S	LD300S	50 Gloves
M	LD300M	50 Gloves
L	LD300L	50 Gloves
XL	LD300X	50 Gloves



Glove Dispenser

Holds one box horizontally or 2 boxes vertically. Mounting hardware included.

SPECS
Wire 10-3/4" x 3-1/8"

PART #
PPGWGR



Nitrile Textured Gloves

Tillotson True Advantage brand

5mil, 9" beaded cuff, ambidextrous.

Made from advanced synthetic rubber formulation that is free of skin irritating proteins, silicone and plasticizers.

Powder Free

SIZE	PART #	QUANTITY PER BOX
S	KT100S	100 Gloves
M	KT100M	100 Gloves
L	KT100L	100 Gloves
XL	KT100X	100 Gloves



Black Nitrile Gloves

Dynarex brand

4mil, 9.5" beaded cuff, ambidextrous.

Ideal for those sensitive to rubber. Extraordinary strength, puncture resistant, a more natural feel than latex or vinyl. Textured fingers for improved grip. Commonly used by tattoo artists, police, DEA, Homeland Security, & postal workers.

Powder Free

SIZE	PART #	QUANTITY PER BOX
S	KB100S	100 Gloves
M	KB100M	100 Gloves
L	KB100L	100 Gloves
XL	KB100X	100 Gloves



Pink Nitrile Gloves

Soft, comfortable, puncture and tear resistant. Both Air-Tite and the manufacturer donate a portion of pink glove sales towards breast cancer research and awareness programs.

Powder Free

SIZE	PART #	QUANTITY PER BOX
XS	KP200XS	200 Gloves
S	KP200S	200 Gloves
M	KP200M	200 Gloves
L	KP200L	200 Gloves
XL	KP200XL	200 Gloves



Shoulder Length Gloves

Economical OB/A.I. glove designed for use by veterinarians and A.I. technicians. The Breeder-Sleeve® is 36" long and 1.0 mil thick. The red gloves are in a convenient dispenser box containing 100 gloves/box. **Latex Free.**

PART# KSP336



MISCELLANEOUS SUPPLIES

DisCide® Hand Sanitizing Wipes

New and improved formula which kills 99.9% of bacteria while moisturizing your skin with Aloe and vitamin E. Fresh wintergreen fragrance leaves skin tingling. Dries in seconds to ensure staff compliance with hand asepsis guidelines.



<u>SPECS</u>	<u>PART #</u>	<u>QTY</u>
6" x 8"	DI3500	160/can

DisCide® Disinfecting Towelettes

Impregnated with a broad spectrum disinfectant, DisCide towelettes are convenient, effective, and user friendly for hard surface disinfecting. Adenovirus, Herpes, HIV, TB are killed in 1 minute. Wiping with DisCide leaves a pleasant fresh herbal scent.



<u>SPECS</u>	<u>PART #</u>	<u>QTY</u>
6.0" x 6.75"	Di1	160/can
10" x 10"	Di2	60/can

Oasis Sterile Surgical Drapes

Sterile surgical drapes made of SMS non-woven material. Consists of two outer layers of SMS and a middle layer of polypropylene. For Veterinary use only.

<u>SPECS</u>	<u>PART #</u>
15" X 26" (2 per pack)	PPGDR1526
18" x 26"	PPGDR1826
40" x 42"	PPGDR4042



Disposable Boots

Disposable boots designed for protective & sanitary use in farm & industrial operations. Elasticized band sealed at the top for easy application & better protection. 4.0 mil thick, 15" tall, 25 pairs/bag. Size M/L fits up to men's size 12.

<u>SPECS</u>	<u>PART #</u>
Size M/L	PPGBC400



Shoe Covers

50 pairs per box
 SPECS PART #
 One size fits all PPG2131



Bouffant Caps

100 per box
 SPECS PART #
 One size fits all PTCP240



Deluxe Dual Head Stethoscope

Chrome binaural silver chest piece with mushroom eartips. Overall length is 30". No color choice.

SPECS	PART #
22" Y-tubing	PPGSTDH



Digital Thermometer

Simple and easy to use with digital readings in 60 seconds. Accurate to 2/10 of a degree. For rectal or oral use. Reads in °C or °F.

SPECS	PART #
Digital thermometer	PPGDIG



Lubricating Jelly

An ideal water soluble lubricant for general needs.

SPECS	PART #
4oz. with flip top	PPGLJ



Miltex® Disposable Biopsy Punches

Seamless, stainless steel, razor-sharp cutting edge. Ribbed handle provides better comfort and control. Convenient, disposable single-use. Sterile, 5 year shelf life.

SPECS	PART #
1mm	MH33-31AA
1.5mm	MH33-31A
2mm	MH33-31
2.5mm	MH33-31B
3mm	MH33-32
3.5mm	MH33-33
4mm	MH33-34
5mm	MH33-35
6mm	MH33-36
8mm	MH33-37



Medicine Vials

Push & turn amber child resistant vials with 5% extra snap caps.

SPECS	PART #	PIECES
6 Dram (22ml)	AV6	600
8 Dram (30ml)	AV8	410
13 Dram (48ml)	AV13	320
20 Dram (73ml)	AV20	360
30 Dram (110ml)	AV30	240
40 Dram (148ml)	AV40	180
60 Dram (220ml)	AV60	115



Green Veterinary Vials

Distinguishing translucent green color. Push & turn caps have veterinary logo. 5% extra snap caps.

SPECS	PART #	PIECES
6 Dram (22ml) *	AV6GNL	600
8 Dram (30ml)	AV8G	410
13 Dram (48ml)	AV13G	320
30 Dram (110ml)	AV30G	240
40 Dram (148ml)	AV40G	180
60 Dram (220ml)	AV60G	115



* 6 Dram are without Veterinary logo on cap.

Plastic Prescription Bottles

Amber plastic prescription bottles with child resistant closures with graduations in ml & oz.

SPECS	PART #	PIECES
1oz.	PPGOVAL1	200
2oz.	PPGOVAL2	200
4oz.	PPGOVAL4	200
8oz.	PPGOVAL8	100
16oz.	PPGOVAL16	50



Works great with our oral syringes & adapters. Use CDA24 for 2oz. - 8oz. and CDA28 for 16oz.

CUSTOM PROJECTS

CONNECTING SUPPLIERS AND CLIENTS, WE BRING THE WORLD TO YOU.

Due to our extensive network of global manufacturers, Air-Tite Products manages custom projects for nearly every industry. Small scale, large scale, high quality or on a budget, if your solution is not available off-the-shelf, contact us today to see what options may be available. When you contact us for a custom project, a team member will be assigned to your project with the sole intention of sourcing your solution.

Sample previous projects include:

- Branded printing and custom coloring on syringes and custom material-specific color ID caps for large chemical company
- Custom length needles developed for government agency for disease testing in large hogs
- Blank syringe barrels with separate plungers for paint sampling.
- Custom needle hub coloring for branded vaccine delivery for international animal health company.
- Custom bevel cuts on hypodermic needles for veterinary application.
- Special bulk packaging for wide variety of customers.
- 3rd party sterilization processing to create “bulk” sterile solutions for mass production in compound pharmaceutical industry.
- Development of unique length / gauge combinations in blunt cannulas for industrial manufacturing operation
- Amber tinted syringes for double blind pharmaceutical clinical trials
- Sourced sharp needles with no hub for specialty ophthalmic application.
- Unique syringe and needle combination for large scale food industry producer.
- Insulin-type fixed needle syringe custom combination for biopharmaceutical application.
- Detachable device for injection aspiration in aesthetics application.
- And many, many others



Whatever your need may be, contact us today and work with a team member to identify the appropriate solution for your specific application.

Welcome to the **Air-Tite** family!

We supply the world with high quality disposable + reusable hypodermic products for Veterinary, Medical, and Laboratory uses; offering personalized service with a 100% satisfaction guarantee. Our sales and service teams are always available to answer your questions and guide your selection with products that improve your business.

In addition to our own line of proven products, we represent and distribute other manufacturers such as:

Henke-Sass, Wolf GmbH (Germany), B. Braun (Germany),
Becton-Dickinson (USA), EXEL (USA), TSK laboratory (Japan),
Poulten & Graf (Germany), plus many more.

**TO REQUEST A CATALOG OR SAMPLES, OR TO PLACE
AN ORDER VISIT**

www.Air-Tite.com
Email atinfo@Air-Tite.com
or call 1.800.231.7762
office: 757.340.2501 fax:757.340.2912
We guarantee our quality 100%

HOURS: MON - FRI 8:30 AM - 6:00 PM Eastern Time

Orders placed by 2:00PM Eastern Time are packaged and shipped the same day

Orders for delivery in continental US are shipped UPS Ground, prepaid and added to your invoice. Expedited delivery available upon request.

Orders outside the continental US shipped Air Freight prepaid and added to your invoice.

Export sales welcomed. Please call for details.

Animal Health Manager:
Ashley Jones
Ashley@Air-Tite.com

Air-Tite Veterinary Distribution Partners


Animal Health International
Armor Animal Health
Clearview Enterprises
Greenwood Products
K Supply
Lambert Veterinary Supply
MWI
Patterson Veterinary Supply
Revival Animal Health
VSI
& More!

2-PART SYRINGES

2-Part Syringes - ideal for Veterinary applications requiring absence of rubber or silicone oil



- **Sterile**
on Page 3
- **Unsterile**
on Page 11



Scan For Website



ANIMAL HEALTH 