


Air-Tite Aesthetics



EXELINT
MEDICAL PRODUCTS

EXEL insulin syringes with fixed needle
Zero Dead Space – ideal for saving expensive
injectables

0.3cc x 30g x 5/16" #EU03308

0.3cc x 31g x 5/16" #EU03318

.5cc x 30g x 5/16" #EU05308

.5cc x 31g x 5/16" #EU05318

*Other configurations available including 1/2"
in 28g / 29g*

Air-Tite Products has a broad portfolio
for all subspecialties of aesthetics.



Premium quality, German engineered syringes
from Henke-Sass, Wolf.

Low Dead Space 1mL Syringes

- NORM-JECT® All-Plastic Luer Slip #A1
- SOFT-JECT® Luer Lock #RL1
- SOFT-JECT® Luer Slip #R1

*Larger sizes from 2.5mL to 60mL available and in
stock*



TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com 🌐 www.air-titeaesthetics.com

Air-Tite Aesthetics



Specialized needles for enhanced precision, efficiency, and comfort.



- Premium Quality Needles
- Hard-to-Find Sizes, for the Specialist
- Available as Small as 33g, which are 22% Thinner than Regular 30g Needles
- Regular Wall, Thin Wall, and Ultra-Thin Wall
- Compatible with all syringes
- ISO Certified
- Made in Japan

STERiJECT™ Needles

GAUGE	LENGTH	CATALOG #
22g	70mm (2-3/4")	TSK2276
	100mm (4")	TSK22101
25g	50mm (2")	TSK2551
27g	6mm (1/4") cap	TSK276
	13mm (1/2") blister	TSK2713U*
30g	4mm (3/16") cap	TSK304
	6mm (1/4") cap	TSK306
	13mm (1/2") cap	TSK3013U*
31g	25mm (1") cap	TSK3025
	13mm (1/2") cap	TSK3113
	13mm (1/2") blister	TSK3113B
32g	4mm (3/16") cap	TSK324
	9mm (3/8") blister	TSK3209**
	13mm (1/2") cap	TSK3213
33g	13mm (1/2") blister	TSK3213B
	13mm (1/2") cap	TSK3313

* Ultra Thin Wall, ideal for dermal fillers

** 32g x 9mm popular for retinal procedures

TO ORDER:

1 (800) 231-7762

1 (757) 340-2912

atinfo@air-tite.com www.air-titeaesthetics.com

Air-Tite Aesthetics



STERiGLIDE™ Microcannulas

Code #	Dimensions
TSK2250SG	22g x 50mm (2")
TSK2270SG	22g x 70mm (2-3/4")
TSK2538SG	25g x 38mm (1-1/2")
TSK2550SG	25g x 50mm (2")
TSK2725SG	27g x 25mm (1")
TSK2738SG	27g x 38mm (1-1/2")
TSK2750SG	27g x 50mm (2")
TSK3025SG	30g x 25mm (1")



Air-Tite also carries the TSK line of STERiJECT cannulas – a high-quality, standard cannula economically priced to fit your budget.

PREMIUM CANNULA

NEAR TIP PRODUCT DELIVERY
ULTRA THIN WALL

DOME SHAPED TIP DESIGN

- Easier Cannula Introduction
- Reduced Risk of Bruising





PROPRIETARY SURFACE TREATMENT

- Up to 50% Better Gliding
- Reduced Patient Discomfort

NEAREST-TO-TIP DELIVERY

- More Accurate Placement
- Reduced Product Loss



CLEAR SIDE-PORT MARKING

- More Accurate Placement
- Better Handling

TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com 🌐 www.air-titeaesthetics.com

Air-Tite Aesthetics


Hypodermic Needles, Syringes + Accessories
for Aesthetic Procedures



TSK MICROCANNULA
STERIGLIDE™

The cannula that glides

TSK 33G x 1/2"
NEEDLE



TSK
NEEDLE HOLDER

Safely secures the hard-shell plastic capsule of cannulas and needles for safer loading and removal from the syringe. # **TSKNH1**

TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com 🌐 www.air-titeaesthetics.com