


YOUR SYRINGE AND NEEDLE SPECIALISTS SINCE 1926

Air-Tite
— Products Co., Inc.

HEALTHCARE AND MEDICAL
CATALOG



Air-Tite Products Co., Inc.

565 Central Drive, Virginia Beach, VA 23454

T: 800.231.7762 F: 757.340.2912 L: 757.340.2501 | Air-Tite.com




WE SUPPLY SPECIALTY HYPODERMIC NEEDLES,
SYRINGES + ACCESSORIES FOR MANY
AESTHETIC PROCEDURES

**AWARD WINNING
MICROCANNULA FOR
BEST INJECTION TOOL**

STERiGLIDE™

The cannula that glides.



See page 19 for more information

Request samples or place orders at:

www.air-tite.com

atinfo@air-tite.com

800.231.7762



**ANTI - AGING &
BEAUTY AWARD**

Each box of **STERiGLIDE™**
is packaged with 20
microcannulas and pilot needles
for entry hole preparation.

TABLE OF CONTENTS

SYRINGES

2-Part Syringes	1
Henke-Ject Syringes	3
Exel Syringes	5
BD Syringes	7

SYRINGES WITH NEEDLES

Henke-Ject Syringes with Needles	9
EXEL Syringes with Needles	11
BD Syringes with Needles	13

NEEDLES

EZ Flow Needles	15
TSK Needles	17
TSK Microcannulas	19
Exel International Needles	21
HSW Needles	23
JBP Needles	25
BD Needles	27

IV AND BLOOD COLLECTION

Intravenous Sets	29
Blood Collection	31

SUPPLIES

Sharps Collection, Sharps Mail-Back System	33
Gloves	35
General Supplies	37

2-PART SYRINGES NORM-JECT[®] / HENKE-JECT[®]

ALL PLASTIC SYRINGES

Norm-Ject[®] syringes are latex free, contain no rubber or silicone oil. These syringes are the choice for any situation needing an inert, hypo-allergenic syringe. Because of their purity applications by users have included IVF, embryo-transfer, chromatography and many laboratory procedures. Manufactured only from laboratory grade polypropylene and polyethylene, these unique plastic syringes have a positive safety stop to help prevent accidental spills. Norm-Ject[®] syringes are individually sterile strip packed and are also available in bulk, unsterile packaging. 1mL - 20mL syringes are Norm-Ject from B Braun Melsungen AG. 30mL and 50mL syringes are Henke-Ject from Henke-Sass, Wolf GmbH. FDA registered, CE Mark, ISO 13485 certified.



CONVENTIONAL RUBBER TIP SYRINGES

Although in general use for medical injections, rubber tipped disposable syringes have several drawbacks compared to the Norm-Ject[®] all plastic syringes. Since they contain natural or synthetic rubber and silicone oil, they may contaminate any media exposed to the syringe. Adverse effects from rubber curing agents, which contain zinc, sulfur and other elements, have been observed on everything from chromatographs to in-vitro procedures. Synthetic rubber plungers are also much more reactive to organic solvents and the silicone oil can cause separation of emulsions. In short, your best choice is the Norm-Ject[®] syringe if you are interested in purity and accuracy.






STERILE 2-PART SYRINGES

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
A1	1mL Tuberculin	Luer Slip	Sterile	100	1800
A3	3mL Centric	Luer Slip	Sterile	100	2500
A5	5mL Centric	Luer Slip	Sterile	100	2000
A10	10mL Eccentric	Luer Slip	Sterile	100	1200
A20	20mL Eccentric	Luer Slip	Sterile	100	800
A30	30mL Eccentric	Luer Slip	Sterile	50	500
A50	50mL Eccentric	Luer Slip	Sterile	30	300

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
AL3	3mL	Luer Lock	Sterile	100	2500
AL5	5mL	Luer Lock	Sterile	100	2000
AL10	10mL	Luer Lock	Sterile	100	1200
AL20	20mL	Luer Lock	Sterile	100	800
AL30	30mL	Luer Lock	Sterile	50	500
AL50	50mL	Luer Lock	Sterile	30	300

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
AC50	50mL	Catheter	Sterile	30	300

** A1 Syringe has Low Dead Space Plunger Tip



STERILE HENKE-JECT® Conventional Syringes

STERILE HENKE-JECT®

Henke-Sass, Wolf GmbH Henke-Ject® (formerly Soft-Ject®) syringes are a high quality disposable rubber piston syringe offering substantial savings. A lubricated synthetic rubber tipped piston allows for smooth draw while a positive stop helps to prevent accidental pullout. These syringes are ETO gas sterilized and are latex and pyrogen free. All Henke-Ject® syringes are packaged in individual sterile strip packs. FDA registered, CE Mark, ISO 13485 certified.





STERILE HENKE-JECT®

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
R1**	1mL Tuberculin	Luer Slip	Sterile	100	2000
R2	2.5mL Centric	Luer Slip	Sterile	100	2500
R5	5mL Centric	Luer Slip	Sterile	100	2000
R10	10mL Eccentric	Luer Slip	Sterile	100	1200
R20	20mL Eccentric	Luer Slip	Sterile	100	800
R30HJ	30mL Eccentric	Luer Slip	Sterile	100	600
R50	50mL Eccentric	Luer Slip	Sterile	50	300

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
RL1 **	1mL	Luer Lock	Sterile	100	2000
RL3	3mL	Luer Lock	Sterile	100	2400
RL5	5mL	Luer Lock	Sterile	100	2000
RL10	10mL	Luer Lock	Sterile	100	1200
RL20	20mL	Luer Lock	Sterile	100	800
RL30HJ	30mL	Luer Lock	Sterile	100	600
RL50	50mL	Luer Lock	Sterile	50	300

**R1 and RL1 have low dead space plunger tip





EXEL INTERNATIONAL® DISPOSABLE SYRINGES

EXEL is known for unsurpassed quality at a low price. These syringes are latex free, sterile, non-pyrogenic, and non-toxic. The design includes bold print, clear barrels, and siliconized gaskets for smooth, consistent plunger motion and a positive stop.


1mL syringes utilize a low dead space plunger

5mL syringes graduated to 6mL

10mL syringes graduated to 12mL

20mL syringes graduated to 25mL

30mL syringes graduated to 35mL






EXEL INTERNATIONAL

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
ES1	1mL Tuberculin	Luer Slip	Sterile	100	1000
ES3	3mL Centric	Luer Slip	Sterile	100	1000
ES5	5mL Centric	Luer Slip	Sterile	100	800
ES10	10mL Centric	Luer Slip	Sterile	100	800
ES20	20mL Eccentric	Luer Slip	Sterile	50	500
ES30	30mL Eccentric	Luer Slip	Sterile	50	250
ES60	60mL Eccentric	Luer Slip	Sterile	25	150

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
EL1	1mL	Luer Lock	Sterile	100	1000
EL3	3mL	Luer Lock	Sterile	100	1000
EL5	5mL	Luer Lock	Sterile	100	800
EL10	10mL	Luer Lock	Sterile	100	800
EL20	20mL	Luer Lock	Sterile	50	500
EL30	30mL	Luer Lock	Sterile	50	250
EL60	60mL	Luer Lock	Sterile	25	150


Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
EC30	30mL	Catheter Tip	Sterile	50	250
EC60	60mL	Catheter Tip	Sterile	25	150



BD[®] SYRINGES



Premium quality from the world's largest syringe manufacturer, Becton Dickinson. BD syringes feature a clear barrel with bold scale markings, tapered plunger rod for ease of aspiration, and positive plunger stop. Available in both Luer-Lok and Luer Slip tips, BD syringes are made in America.






BD - BECTON DICKENSON

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
BD309659	1mL Tuberculin	Luer Slip	Sterile	200	1600
BD309656	3mL Centric	Luer Slip	Sterile	200	800
BD309647	5mL Centric	Luer Slip	Sterile	125	500
BD303134	10mL Centric	Luer Slip	Sterile	200	400
BD302831	20mL Centric	Luer Slip	Sterile	48	192
BD302833	30mL Centric	Luer Slip	Sterile	56	224
BD309654	50mL Centric	Luer Slip	Sterile	40	160

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
BD309628	1mL Polycarbonate	Luer Lock	Sterile	100	800
BD309657	3mL	Luer Lock	Sterile	200	800
BD309646	5mL	Luer Lock	Sterile	125	500
BD302995	10mL	Luer Lock	Sterile	200	400
BD302830	20mL	Luer Lock	Sterile	48	192
BD302832	30mL	Luer Lock	Sterile	56	224
BD309653	50mL	Luer Lock	Sterile	40	160

Catalog Number	Specifications	Tip Type	Sterility	Quantity per Box	Quantity per Case
BD309620	50mL	Catheter Tip	Sterile	40	160





HSW

Henke-Ject® (Formerly "Soft-Ject®") fixed needle syringes offer a crystal clear cylinder with easy to read scale. The piston, piston-rod, and cylinder consist of polypropylene. The fixed needle syringes are free of pyrogens, PVC, and latex, and are non-toxic. The 0.5cc syringes are in a sterile polybag with 10 pieces each, the 1cc syringes are sterile individually blister packaged.

The By-packs from Henke-Sass, Wolf offer a convenient syringe with needle combination packaged side-by-side in individual strip packs.





HENKE-JECT FIXED NEEDLE SYRINGES

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
RU13012	1cc Tuberculin	30g x 1/2"	Sterile	100	2000

HENKE-JECT SYRINGES / NEEDLE BY-PACKS

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
R22234	2.5mL LS	22g x 3/4"	Sterile	100	2500
R2221	2.5mL LS	22g x 1"	Sterile	100	2500
R22558	2.5mL LS	25g x 5/8"	Sterile	100	2500
R2018	20mL LS	18g x 1-1/2"	Sterile	100	800



EXEL

EXEL's Comfort Point fixed needle syringes offer a permanently attached needle, reduced dead space and smooth plunger action for greater filling accuracy. The auto-ground needle tip offers greater comfort. These syringes are sterile, non-toxic, and non-pyrogenic. All EXEL fixed needle syringes are packaged 10 per bag and 100 per box.

EXEL syringes with mounted needle offer convenience at an economic price. A wide selection of luer slip and luer lock syringes with standard size needles pre-attached. These syringes allow the user to reduce the number of inventory items and there is no need to open multiple packs to assemble a ready-to-use syringe.

EXEL FIXED NEEDLE SYRINGES

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
EU032913	0.3cc Insulin	29g x 1/2"	Sterile	100	500
EU03308	0.3cc Insulin	30g x 5/16"	Sterile	100	500
EU03318	0.3cc Insulin	31g x 5/16"	Sterile	100	500
EU052813	0.5cc Insulin	28g x 1/2"	Sterile	100	500
EU052913	0.5cc Insulin	29g x 1/2"	Sterile	100	500
EU05308	0.5cc Insulin	30g x 5/16"	Sterile	100	500
EU05318	0.5cc Insulin	31g x 5/16"	Sterile	100	500
EU12558	1cc tuberculin	25g x 5/8"	Sterile	100	500
EU12813	1cc insulin	28g x 1/2"	Sterile	100	500
EU12913	1cc insulin	29g x 1/2"	Sterile	100	500
EU1308	1cc insulin	30g x 5/16"	Sterile	100	500
EU1318	1cc insulin	31g x 5/16"	Sterile	100	500

EXEL LUER SLIP SYRINGES WITH MOUNTED NEEDLE

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
ES12558	1cc	25g x 5/8"	Sterile	100	1000
ES1251	1cc	25g x 1"	Sterile	100	1000
ES12612	1cc	26g x 1/2"	Sterile	100	1000
ES12712	1cc	27g x 1/2"	Sterile	100	1000
ES3201	3cc	20g x 1"	Sterile	100	1000
ES3211	3cc	21g x 1"	Sterile	100	1000
ES32234	3cc	22g x 3/4"	Sterile	100	1000
ES3221	3cc	22g x 1"	Sterile	100	1000
ES3231	3cc	23g x 1"	Sterile	100	1000
ES32558	3cc	25g x 5/8"	Sterile	100	1000

EXEL LUER LOCK SYRINGES WITH MOUNTED NEEDLE

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Quantity per Case
EL318112	3cc	18g x 1-1/2"	Sterile	100	1000
EL3201	3cc	20g x 1"	Sterile	100	1000
EL320112	3cc	20g x 1-1/2"	Sterile	100	1000
EL3211	3cc	21g x 1"	Sterile	100	1000
EL321114	3cc	21g x 1-1/4"	Sterile	100	1000
EL321112	3cc	21g x 1-1/2"	Sterile	100	1000
EL32234	3cc	22g x 3/4"	Sterile	100	1000
EL3221	3cc	22g x 1"	Sterile	100	1000
EL322114	3cc	22g x 1-1/4"	Sterile	100	1000
EL322112	3cc	22g x 1-1/2"	Sterile	100	1000
EL3231	3cc	23g x 1"	Sterile	100	1000
EL323112	3cc	23g x 1-1/2"	Sterile	100	1000
EL32558	3cc	25g x 5/8"	Sterile	100	1000
EL3251	3cc	25g x 1"	Sterile	100	1000
EL325112	3cc	25g x 1-1/2"	Sterile	100	1000
EL327114	3cc	27g x 1-1/4"	Sterile	100	1000
EL5201	5cc	20g x 1"	Sterile	100	800
EL5211	5cc	21g x 1"	Sterile	100	800
EL521112	5cc	21g x 1-1/2"	Sterile	100	800
EL5221	5cc	22g x 1"	Sterile	100	800
EL522112	5cc	22g x 1-1/2"	Sterile	100	800
EL10181	10cc	18g x 1"	Sterile	100	800
EL10201	10cc	20g x 1"	Sterile	100	800
EL10211	10cc	21g x 1"	Sterile	100	800
EL1021112	10cc	21g x 1-1/2"	Sterile	100	800
EL10221	10cc	22g x 1"	Sterile	100	800
EL1022112	10cc	22g x 1-1/2"	Sterile	100	800





BD

Becton Dickinson fixed needle syringes offer Ultra-Fine™ permanently attached needles. BD insulin syringes are the #1 selling brand in the US.

BD Mounted Needle syringes offer a convenient BD syringe and Precision Glide™ needle - no assembly required. The clear barrels on these syringes features bold scales and the needle helps improve dosage accuracy.

Air-Tite carries a mixture of BD Safety syringes in the most popular sizes. The Eclipse line allows clinicians to activate the safety mechanism after the injection while the pivoting shield offers low angle injections and efficient safety. The Safety-Lok™ line offers the only syringe of its type and affords clear visual confirmation of the sliding sleeve in the lock position. SafetyGlide™ has a single-hand activated shielding mechanism.

BD FIXED NEEDLE SYRINGES

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box
BD328438	3/10 mL Insulin	31g x 5/16"	Sterile	100
BD328468	1/2 mL Insulin	31g x 5/16"	Sterile	100
BD328418	1 mL Insulin	31g x 5/16"	Sterile	100

BD SYRINGES WITH MOUNTED NEEDLES

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box
BD309623	1mL TB Luer Slip	27g x 1/2"	Sterile	100
BD609637	1mL Luer Lock	20g x 1"	Sterile	100
BD309580	3mL Luer Lock	18g x 1-1/2"	Sterile	100
BD309578	3mL Luer Lock	20g x 1"	Sterile	100
BD309579	3mL Luer Lock	20g x 1-1/2"	Sterile	100
BD309572	3mL Luer Lock	22g x 1"	Sterile	100
BD309574	3mL Luer Lock	22g x 1-1/2"	Sterile	100
BD309570	3mL Luer Lock	25g x 5/8"	Sterile	100

BD SAFETY SYRINGES WITH MOUNTED NEEDLES

Catalog Number	Syringe Size	Needle	Sterility	Quantity per Box	Safety Brand
BD305781	3mL	25g x 5/8"	Sterile	50	Eclipse
BD305782	3mL	23g x 1"	Sterile	50	Eclipse
BD305946	1mL	26g x 3/8"	Sterile	100	Safety Glide
BD305906	3mL	22g x 1-1/2"	Sterile	50	Safety Glide





Smooth Entry • Smooth Delivery

Air-Tite Product's newest addition, the EZ FLO line of hypodermic needles, offer high quality and a broad selection at an affordable price. Ranging from 18 gauge to 30 gauge, these needles offer a standard bevel and regular-wall to combine comfort and durability. Each needle is sterile, individually packaged in a box of 100. Made of surgical stainless steel and a polypropylene hub. 18g to 25g are thin wall. 27g and 30g are regular wall thickness.



EZ FLO BY AIR-TITE

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NEZ18100	18	1"	Sterile	100	2000
NEZ18112	18	1-1/2"	Sterile	100	2000
NEZ18200	18	2"	Sterile	100	2000
NEZ19100	19	1"	Sterile	100	2000
NEZ19112	19	1-1/2"	Sterile	100	2000
NEZ20100	20	1"	Sterile	100	2000
NEZ20112	20	1-1/2"	Sterile	100	2000
NEZ21100	21	1"	Sterile	100	2000
NEZ21112	21	1-1/2"	Sterile	100	2000
NEZ21200	21	2"	Sterile	100	2000
NEZ22034	22	3/4"	Sterile	100	2000
NEZ22100	22	1"	Sterile	100	2000
NEZ22112	22	1-1/2"	Sterile	100	2000
NEZ22200	22	2"	Sterile	100	2000
NEZ22300	22	3"	Sterile	100	2000
NEZ22400	22	4"	Sterile	100	2000
NEZ23034	23	3/4"	Sterile	100	2000
NEZ23100	23	1"	Sterile	100	2000
NEZ23112	23	1-1/2"	Sterile	100	2000
NEZ25058	25	5/8"	Sterile	100	2000
NEZ25100	25	1"	Sterile	100	2000
NEZ25112	25	1-1/2"	Sterile	100	2000
NEZ25200	25	2"	Sterile	100	2000
NEZ27100	27	1"	Sterile	100	2000
NEZ27112	27	1-1/2"	Sterile	100	2000
NEZ27114	27	1-1/4"	Sterile	100	2000
NEZ27200	27	2"	Sterile	100	2000
NEZ30012	30	1/2"	Sterile	100	2000
NEZ30100	30	1"	Sterile	100	2000



Needles





TSK STERiJECT™ hypodermic needles are produced with a focus on needle sharpness and quality. All needles are thin-wall or ultra-thin wall. Made in Japan from one of the global leaders for ultra-fine tip needles, TSK, these needles are commonly used for aesthetic procedures, mesotherapy, and ophthalmology among other uses. Available in hard-to-find sizes for the specialist. 100 needles per box, surgical stainless steel, and poly hub.

TSK STERIJECT™ NEEDLES

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box	Quantity per Case
TSK2276	22	70mm	Thin	Sterile	100	3000
TSK22101	22	100mm	Thin	Sterile	100	3000
TSK2551	25	50mm	Thin	Sterile	100	3000
TSK2713U	27	13mm	Ultra-Thin	Sterile	100	6000
TSK304	30	4mm	Thin	Sterile	100	6000
TSK306	30	6mm	Thin	Sterile	100	6000
TSK3013U	30	13mm	Ultra-Thin	Sterile	100	6000
TSK3113	31	13mm	Thin	Sterile	100	6000
TSK324	32	4mm	Thin	Sterile	100	6000
TSK3206	32	6mm	Thin	Sterile	100	6000
TSK3209C	32	9mm	Thin	Sterile	100	6000
TSK3213	32	13mm	Thin	Sterile	100	6000
TSK3304	33	4mm	Thin	Sterile	100	6000
TSK3313	33	13mm	Thin	Sterile	100	6000



TSK Microcannulas



TSK STERiJECT™ CSH CANNULA



ROUNDED, CLOSED-TIP DESIGN

- Reduced Risk of Trauma
- Easier Handling



Conventional needle

VALUE PRICED CANNULA

- Traditional Product Delivery
- Ultra Thin Wall

ULTRA THIN WALL

- Less force required for injection
- Smoother and More Controlled Injection



Regular Wall

CODE #	DESCRIPTION	DIMENSIONS
TSK2250SP25	25-PACK STERiJECT CSH Microcannulas	22g x 50mm
TSK2538SP25	25-PACK STERiJECT CSH Microcannulas	25g X 38mm
TSK2550SP25	25-PACK STERiJECT CSH Microcannulas	25g x 50mm
TSK2738SP25	25-PACK STERiJECT CSH Microcannulas	27g x 38mm
TSK2750SP25	25-PACK STERiJECT CSH Microcannulas	27g x 50mm

TSK STERiGLIDE™ CANNULA

AWARD WINNING BEST INJECTOR



CODE #	DIMENSIONS
TSK 2238SG	22g x 38mm (1½")
TSK 2250SG	22g x 50mm (2.0")
TSK 2270SG	22g x 70mm (2¾")
TSK 2538SG	25g x 38mm (1½")
TSK 2550SG	25g x 50mm (2.0")
TSK 2725SG	27g x 25mm (1.0")
TSK 2738SG	27g x 38mm (1½")
TSK 2750SG	27g x 50mm (2.0")
TSK 3025SG	30g x 25mm (1.0")

Each box contains 20 microcannulas and pilot needles.

Needles

PREMIUM CANNULA

Smooth Gliding, Near Tip Product Delivery, Ultra-Thin Wall

Flexible enough to navigate anatomy, rigid enough to place accurately, with improved tactile sensitivity

DOME SHAPED TIP DESIGN

- Improved cannula introduction
- Reduced risk of trauma
- Easier handling




NEAREST-TO-TIP DELIVERY

- More Accurate Placement
- Reduced Product Loss



PROPRIETARY SURFACE TREATMENT

- Less force to advance
- Reduced Patient Discomfort



CLEAR SIDE-PORT MARKING

- Directional accuracy
- Better Handling





EXEL International[®] offers a variety of hypodermic needles with regular wall thickness, made of 304 surgical stainless steel, non-toxic and pyrogen-free. Each box is color coded to help the clinician easily identify the appropriate gauge. Each needle is individually packaged, sterilized, and ready to use.

EXEL INTERNATIONAL[®] HYPODERMIC NEEDLES

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
EN181	18	1"	Sterile	100	2,000
EN18112	18	1-1/2"	Sterile	100	2,000
EN191	19	1"	Sterile	100	2,000
EN19112	19	1-1/2"	Sterile	100	2,000
EN201	20	1"	Sterile	100	2,000
EN20112	20	1-1/2"	Sterile	100	2,000
EN211	21	1"	Sterile	100	2,000
EN21112	21	1-1/2"	Sterile	100	2,000
EN2234	22	3/4"	Sterile	100	2,000
EN221	22	1"	Sterile	100	2,000
EN22112	22	1-1/2"	Sterile	100	2,000
EN2334	23	3/4"	Sterile	100	2,000
EN231	23	1"	Sterile	100	2,000
EN23112	23	1-1/2"	Sterile	100	2,000
EN2558	25	5/8"	Sterile	100	2,000
EN251	25	1"	Sterile	100	2,000
EN25112	25	1-1/2"	Sterile	100	2,000
EN2612	26	1/2"	Sterile	100	2,000
EN2712	27	1/2"	Sterile	100	2,000
EN27114	27	1-1/4"	Sterile	100	2,000
EN27112	27	1-1/2"	Sterile	100	2,000
EN2834	28	3/4"	Sterile	100	1,000
EN3012	30	1/2"	Sterile	100	2,000
EN301	30	1"	Sterile	100	2,000
EN3212	32	1/2"	Sterile	100	1,000



SecureTouch needles feature one-handed, finger-stroke activation. No solid surface required. An audible click as well as tactile feedback and visual indicator confirm safety mechanism engagement. Needle remains locked inside activated cover. Fits any standard brand syringe.

EXEL® SECURETOUCH SAFETY NEEDLES

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
ENS18112	18	1-1/2"	Sterile	100	1.000
ENS201	20	1"	Sterile	100	1.000
ENS231	23	1"	Sterile	100	1.000
ENS2558	25	5/8"	Sterile	100	1.000
ENS251	25	1"	Sterile	100	1.000
ENS25112	25	1-1/2"	Sterile	100	1.000
ENS27114	27	1-1/4"	Sterile	100	1.000
ENS27112	27	1-1/2"	Sterile	100	1.000
ENS3012	30	1/2"	Sterile	100	1.000



Needles







Henke-Ject® (formerly Fine-Ject®) hypodermic needles by Henke-Sass, Wolf GmbH are high quality yet economical. Made of surgical stainless steel, these needles are sterile, disposable, individually wrapped, and available in thin wall up to 25 gauge.



HENKE-JECT® HYPODERMIC NEEDLES BY HENKE-SASS, WOLF GMBH

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box	Quantity per Case
NH15112	15	1-1/2"	Thin	Sterile	100	4,000
NH152	15	2"	Thin	Sterile	100	3,000
NH161	16	1"	Thin	Sterile	100	4,000
NH16112	16	1-1/2"	Thin	Sterile	100	4,000
NH172	17	2"	Thin	Sterile	100	3,000
NH18112	18	1-1/2"	Thin	Sterile	100	4,000
NH182	18	2"	Thin	Sterile	100	3,000
NH201	20	1"	Thin	Sterile	100	4,000
NH20112	20	1-1/2"	Thin	Sterile	100	4,000
NH211	21	1"	Thin	Sterile	100	4,000
NH212	21	2"	Thin	Sterile	100	3,000
NH221	22	1"	Thin	Sterile	100	4,000
NH22112	22	1-1/2"	Thin	Sterile	100	4,000
NH222	22	2"	Thin	Sterile	100	3,000
NH231	23	1"	Thin	Sterile	100	4,000
NH2558	25	5/8"	Thin	Sterile	100	4,000
NH251	25	1"	Thin	Sterile	100	4,000
NH252	25	2"	Thin	Sterile	100	3,000
NH2712	27	1/2"	Regular	Sterile	100	4,000
NH272	27	2"	Regular	Sterile	100	3,000
NH3012	30	1/2"	Regular	Sterile	100	4,000



Needles




JBP (S. Korea) Nanoneedles™ are the integration of craftsmanship and technology which has led to an ultra-thin wall needle making injections of high viscosity materials through extra fine needles possible. JBP uses a special process called Micron Sharpness to reduce frictional forces. With gauges ranging from 25 to 34, these needles are ideally suited for those seeking ultra-fine tips and the smallest of injection needles.



JBP NANONEEDLES™

Catalog Number	Gauge	Length	Wall Type	Package Type	Sterility	Quantity per Box	Quantity per Case
JBP2713B	27	13mm	Ultra-thin	Blister	Sterile	100	800
JBP3004	30	4mm	Ultra-thin	Capsule	Sterile	100	800
JBP3004B	30	4mm	Ultra-thin	Blister	Sterile	100	800
JBP3008	30	8mm	Ultra-thin	Capsule	Sterile	100	800
JBP3013	30	13mm	Ultra-thin	Capsule	Sterile	100	800
JBP3013B	30	13mm	Ultra-thin	Blister	Sterile	100	800
JBP3304	33	4mm	Ultra-thin	Capsule	Sterile	100	800
JBP3304B	33	4mm	Ultra-thin	Blister	Sterile	100	800
JBP3308	33	8mm	Ultra-thin	Capsule	Sterile	100	800
JBP3308B	33	8mm	Ultra-thin	Blister	Sterile	100	800
JBP3313	33	13mm	Ultra-thin	Capsule	Sterile	100	800
JBP3313B	33	13mm	Ultra-thin	Blister	Sterile	100	800
JBP3404	34	4mm	Ultra-thin	Capsule	Sterile	100	800
JBP3406	34	6mm	Ultra-thin	Capsule	Sterile	100	800
JBP3408	34	8mm	Ultra-thin	Capsule	Sterile	100	800



CONVENTIONAL NEEDLE

VS



JBP NANONEEDLE (UTW)



BD Precision Glide™ needles are available in a wide range of gauges, lengths and bevel designs.

BD Safety Glide™ Needles feature BD Activation Assist™ technology making a safe injection easier and more efficient.





BD® PRECISION GLIDE™ NEEDLES

Catalog Number	Gauge	Length	Wall Type	Sterility	Quantity per Box
BD305198	16	1-1/2"	Regular	Sterile	100
BD305216	16	1"	Thin	Sterile	100 NOKOR ADMIX
BD305195	18	1"	Regular	Sterile	100
BD305196	18	1-1/2"	Regular	Sterile	100
BD305215	18	1-1/2"	Thin	Sterile	100 NOKOR ADMIX
BD305180	18	1-1/2"	Regular	Sterile	100 (Blunt tip)
BD305200	19	1-1/2"	Thin	Sterile	100 Filter Needle
BD305175	20	1"	Regular	Sterile	100
BD305176	20	1-1/2"	Regular	Sterile	100
BD305167	21	1-1/2"	Regular	Sterile	100
BD305129	21	2"	Regular	Sterile	100
BD305155	22	1"	Regular	Sterile	100
BD305156	22	1-1/2"	Regular	Sterile	100
BD305145	23	1"	Regular	Sterile	100
BD305120	23	1-1/4"	Regular	Sterile	100
BD305122	25	5/8"	Regular	Sterile	100
BD305125	25	1"	Regular	Sterile	100
BD305127	25	1-1/2"	Regular	Sterile	100
BD305110	26	3/8"	Regular	Sterile	100
BD305115	26	5/8"	Regular	Sterile	100 Sub-Q
BD305109	27	1/2"	Regular	Sterile	100
BD305136	27	1-1/4"	Regular	Sterile	100
BD305106	30	1/2"	Regular	Sterile	100
BD305128	30	1"	Regular	Sterile	100

BD® SAFETY GLIDE NEEDLES

Catalog Number	Gauge	Length	Sterility	Quantity per Case
BD305915	21	1"	Sterile	500
BD305917	21	1-1/2"	Sterile	500
BD305900	22	1-1/2"	Sterile	500
BD305902	23	1"	Sterile	500
BD305901	25	5/8"	Sterile	500
BD305916	25	1"	Sterile	500



INTRAVENOUS SETS

EXEL IV ADMINISTRATION

Catalog Number	Specs	Sterility	Quantity per Case
EIV7815	78" Tubing, 15 drops/mL, vented / non-vented, Luer Slip	Sterile	50
EIV7815LL	78" Tubing, 15 drops/mL, vented / non-vented, Luer Lock	Sterile	50
EIV7860	78" Tubing, 60 drops/mL, vented / non-vented, Luer Slip	Sterile	50
EIV7860LL	78" Tubing, 60 drops/mL, vented / non-vented, Luer Lock	Sterile	50
EIV8420F	0.2 Micron IV filter, 4" tubing, roller clamp and occlusion clamp	Sterile	50
EIV10515	Mullti-site, 2 Y-injection, 105" tubing, 28mL priming volume	Sterile	50

EXEL IV CATHETERS

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EIV14114	14g x 1-1/4"	Sterile	50	200
EIV142	14g x 2"	Sterile	50	200
EIV16114	16g x 1-1/4"	Sterile	50	200
EIV162	16g x 2"	Sterile	50	200
EIV18114	18g x 1-1/4"	Sterile	50	200
EIV18112	18g x 1-1/2"	Sterile	50	200
EIV182	18g x 2"	Sterile	50	200
EIV201	20g x 1"	Sterile	50	200
EIV20114	20g x 1-1/4"	Sterile	50	200
EIV202	20g x 2"	Sterile	50	200
EIV221	22g x 1"	Sterile	50	200
EIV2434	24g x 3/4"	Sterile	50	200

EXEL WINGED INFUSION SETS


Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EWI1934	19g x 3/4"	Sterile	50	500
EWI2034	20g x 3/4"	Sterile	50	500
EWI2134	21g x 3/4"	Sterile	50	500
EWI2334	23g x 3/4"	Sterile	50	500
EWI2534	25g x 3/4"	Sterile	50	500
EWI2734	27g x 3/4"	Sterile	50	500

EXEL SAFETY WINGED INFUSION SETS

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EWI2134S	21g x 3/4" Safety	Sterile	50	500
EWI2334S	23g x 3/4" Safety	Sterile	50	500
EWI2534S	25g x 3/4" Safety	Sterile	50	500

EXEL INJECTION PLUGS

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EIP2	Injection Plug with Cap Luer Lock	Sterile	50	200



BLOOD COLLECTION

EXEL BLOOD COLLECTION SETS WITH 12" TUBING

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EBS2134	21g x 3/4"	Sterile	50	200
EBS2334	23g x 3/4"	Sterile	50	200
EBS2534	25g x 3/4"	Sterile	50	200



EXEL SAFETY BLOOD COLLECTION SETS WITH 12" TUBING

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EBS2134S	21g x 3/4"	Sterile	50	200
EBS2334S	23g x 3/4"	Sterile	50	200
EBS2534S	25g x 3/4"	Sterile	50	200




EXEL MULTI-TYPE BLOOD COLLECTION NEEDLES

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EBC201	20g x 1"	Sterile	100	1000
EBC20112	20g x 1-1/2"	Sterile	100	1000
EBC211	21g x 1"	Sterile	100	1000
EBC21112	21g x 1-1/2"	Sterile	100	1000
EBC221	22g x 1"	Sterile	100	1000
EBC22112	22g x 1-1/2"	Sterile	100	1000



BD VACUTAINER® BLOOD COLLECTION TUBES

Catalog Number	Specs	Sterility	Quantity per Case
BD366668	3mL Red, No Additive, w/ Silicone	Sterile	100
BD367835	3mL Lavender, EDTA no Silicone	Sterile	100
BD367812	4mL Red, No Additive, w/ Silicone	Sterile	100
BD367815	6mL Red, No Additive, w/ Silicone	Sterile	100
BD366430	10mL Red, No Additive, w/ Silicone	Sterile	100
BD367985	10mL Red/Grey, Plus SST Brand	Sterile	100




EXEL BLOOD COLLECTION ADAPTOR

For use with winged infusion sets

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EWI2BC	Adaptor	Sterile	100	1000

EXEL BLOOD COLLECTION HOLDER WITH MULTI-SAMPLE ADAPTOR

Catalog Number	Specs	Sterility	Quantity per Box	Quantity per Case
EBCA	Holder with Adaptor	Sterile	50	200



Sharps Collection Sharps Mail-Back System



Bemis® SharpSentinel® Containers

Part Number	Specifications	Quantity	Quantity Per case
SC101R	1 Quart (4.25"L x 4.25"W x 6.5"H - fits phlebotomy trays)	Each	100
SC105R	8 Quart(11.75"L x 7.75"W x 8.75"H - 2 gallons - table top or wall mount)	Each	32
SC108R	8 Quart (11.75"L x 7.75"W x 8.75"H - 2 gallons - table top or wall mount)	Each	30
SC112R	12 Quart (15"L x 6.5"W x 15"H - 3 gallons)	Each	12
SC212R	12 Quart Rotating Lid (15"L x 6.5"W x 15"H - 3 gallons)	Each	12
SC132R	8 Gallon (15.25"L x 11"W x 15.5"H)	Each	10
SC144R	11 Gallon (15.25"L x 11"W x 22"H)	Each	6
SC101H	1 Qt Holder - fits SC101R	Each	20
SC100B	SharpSentinel® Mounting Bracket (fit 5 quart (SC105R) & 8 quart (SC108R))	Each	10
SC100K	SharpSentinel® Bracket key - fits SC100B	Each	10
SC100W	3 Gallon Wire Bracket (Fits 12 quart SC112R & SC212R)	Each	5
SC5QWC	5 Quart Wall Cabinet	Each	
SC3GWC	3 Gallon Wall Cabinet	Each	
SCGH1	Glove Box Holder	Each	

Non-Sharps Medical Waste Mail-Back Systems **May not contain any LOOSE sharps

Part Number	Specifications	Dimension
MWM1G	1 Gallon System	8.375" diameter x 8" high
SCM5G	5 Gallon System	10.45" diameter x 14.51" high

Sharps Mail-Back System

Each Kit includes all components required for collection and transportation of regulated medical waste.

Each Kit includes:

Step-by-step instructions, waste shipping containers, prepaid and pre-addressed shipping labels, prepaid disposal & tracking documents.



5-1 quart

1 gallon

2 gallon

3 gallon

Part Number	Specifications	Dimension
SCM1Q	1 Quart System	Container Dimensions: 3.75" x 3.5" x 8" Each container holds approximately 42 3cc syringes.
SCM1Q-5	5-1 Quart System	Container Dimensions: 3.75" x 3.5" x 8" Each container holds approximately 42 3cc syringes.
SCM1G	1 Gallon System	Container Dimensions: 9" x 5.5" x 8.5" Each container holds approximately 98 3cc syringes.
SCM1G-2	2-1 Gallon System	Container Dimensions: 9" x 5.5" x 8.5" Each container holds approximately 98 3cc syringes.
SCM2G	2 Gallon System	Container Dimensions: 9" x 5.5" x 14.25" Each container holds approximately 245 3cc syringes.
SCM2G-2	2-2 Gallon System	Container Dimensions: 9" x 5.5" x 14.25" Each container holds approximately 245 3cc syringes.
SCM3G	3 Gallon System	Container Dimensions: 9" x 5.5" x 19.75" Each container holds approximately 358 3cc syringes.

Universal Brackets

Universal mounting brackets fit 1 quart, 1 gallon, 2 gallon, and 3 gallon Mail-Back systems (see above). Locking brackets include lock & key.

Part Number	Specifications
SCMSB123	Universal bracket
SCMSB123L	Universal bracket - locking

GLOVES



Latex Exam Gloves *100 gloves/box.*

Part Number	Specifications	Quantity	Quantity Per case
LS100	Powder FREE S, M, L, XL	100	1000
LS100S	Powder FREE Small	100	1000
LS100M	Powder FREE Medium	100	1000
LS100L	Powder FREE Large	100	1000
LS100X	Powder FREE XL	100	1000

Vinyl Exam Gloves *100 gloves per box.*

Part Number	Specifications	Quantity	Quantity Per case
VS100	Powder FREE S, M, L, XL	100	1000



Latex High Risk Gloves

50 gloves per box. Extra long, textured.

Part Number	Specifications	Quantity	Quantity Per case
LD300	Powder FREE S, M, L, XL	50	500

Nitrile Textured Gloves

100 gloves/box

Part Number	Specifications	Quantity	Quantity Per case
KT100	Powder FREE S, M, L, XL	100	1000

Black Nitrile Gloves

100 gloves/box

Part Number	Specifications	Quantity	Quantity Per case
KB100	Powder FREE S, M, L, XL	100	1000

SUPPLIES

Air-Tite Products Co. offers a wide range of supplies for the clinicians practice. From Sharps Containers to Personal Protective Equipment, from gauze to Instruments - Air-Tite carries these products for your convenience. Please see the current Price List or visit our website at www.air-tite.com for part numbers, pricing, ordering and more information.

HIGHLIGHT PRODUCTS

TSKNH1 is a high quality, reusable needle stand designed to allow a safe procedure for engaging the syringe to the needle, while also providing a convenient holder during the procedure.

IM57-400 syringe connector female-to-female with low dead space, Made in the USA.



GENERAL SUPPLIES

Access to nearly the entire range of Miltex surgical instruments including scissors, needle holders, forceps, clamps and many more.



Oral dose syringes are available.
Contact us for more details



Gauze, available in Sterile, Unsterile, Woven, and Non-woven



Supplies

In addition to the supplies mentioned above, Air-Tite Products offers Alcohol and Iodine Prep Pads, Cotton Products (Rolls, Balls, Tip Applicators), Disinfectant wipes, Hand Creams, Tongue depressors, Cohesive bandage wraps, Thermometers, Stethoscopes, Instrument Trays, Lubricating Jelly, and more. Please see the price list for our full offerings, or call our Customer Service team today!



TSK

Compared to conventional needles, the LDS hub design minimizes product loss which may lead to significant cost savings when injecting high cost drugs.



LDS hub

Regular hub

The LOW DEAD SPACE needle hub is designed to reduce the effective material left behind as close to zero as possible.



www.air-tite.com



atinfo@air-tite.com

800.231.7762

Welcome to the **Air-Tite** family!

We supply the world with high quality disposable + reusable hypodermic product for Laboratory. Medical and Veterinary use offering personalized service with a 100% satisfaction guarantee. Our sales and service team are always available to answer your questions and guide your selection with product that improve your business.

In addition to our own line of proven product, we represent and distribute other manufacturers such as:
Henke-Sass, Wolf GmbH (Germany), B. Braun (Germany), Becton-Dickinson (USA), EXEL (USA), TSK laboratory (Japan), Poulten Graf (Germany) plus many more.

TO REQUEST A CATALOG OR SAMPLES, OR TO PLACE AN ORDER VISIT

www.Air-Tite.com
Email atinfo@Air-Tite.com
or call 1.800.231.7762
office: 757.340.2501 fax:757.340.2912
We guarantee our quality 100%

HOURS: MON - FRI 8:30 am - 6:00 PM Eastern Time

Orders placed by 2:00PM Eastern Time are packaged and shipped the same day

Orders for delivery in continental US are shipped UPS Ground, prepaid and added to your invoice. Expedited delivery available upon request.
Orders outside the continental US shipped Air Freight prepaid and added to your invoice.

Export sales welcomed. Please call for details.

CUSTOMER SERVICE:	
Carolyn Ang	Carolyn@Air-Tite.com
Danna Williams.....	Danna@Air-Tite.com
Donna Rilee	DonnaR@air-tite.com
Cindy Turner	Cindy@Air-Tite.com
Jan Wahl	JanW@Air-Tite.com
DIRECTOR, SALES + MARKETING	CRAIG JOHNSON CraigJ@Air-Tite.com
HEALTHCARE SALES + MARKETING.....	WILL FOLEY WillF@Air-Tite.com
VETERINARY SALES + MARKETING.....	SAM GASSEL SamG@Air-Tite.com
LOGISTICS MANAGER	DAVID HARP DavidH@Air-Tite.com

ICE-STICK

AESTHETIC COOLING MASSAGER



This cooling device from Aesthetic Group in France provides a great non-drip solution in pre- and post-procedure for your patients. The ICE-STICK core can be removed from the plastic casing and placed in your freezer until needed. Used frequently in aesthetic medicine practices, this ICE-STICK can be cleaned and reused.

The ICE-STICK is composed of an ice-cold cartridge with refreshing aluminum head to cool areas of the face and body. Helpful prior to injection for pain reduction and to reduce the risk of hematomas.



