

Air-Tite Aesthetics

ICE-STICK

AESTHETIC COOLING MASSAGER

Aesthetic
Group

Ice-Stick

This cooling device from Aesthetic Group in France provides a great non-drip solution in pre- and post-procedure for your patients. The ICE-STICK core can be removed from the plastic casing and placed in your freezer until needed. Used frequently in aesthetic medicine practices, this ICE-STICK can be cleaned and reused.

The ICE-STICK is composed of an ice-cold cartridge with refreshing aluminum head to cool areas of the face and body. Helpful prior to injection for pain reduction and to reduce the risk of hematomas.

Warning – this product becomes very cool and shouldn't be applied directly to skin for a minimum of 5 minutes after removal from freezer. Clean and disinfect the product before use. Use only on the skin. Not a medical device. Do not boil or leave near a heated surface.



Air-Tite
Aesthetics

☎ **1 (800) 231-7762**
1 (757) 340-2501

✉ atinfo@air-tite.com

🌐 www.air-tite.com