

# NORM-JECT®



## Norm-Ject® All Plastic Syringes

Norm-Ject® syringes are latex free, contain no rubber, silicone oil, styrene, BPA, or DEHP and are certified DNA free. These syringes are the choice for any situation needing an inert, hypo-allergenic syringe. Because of their purity they are indicated for IVF, embryo-transfer, chromatography and many laboratory procedures. Manufactured only from laboratory grade polypropylene and polyethylene, these unique plastic syringes have a positive safety stop to prevent accidental spills. Norm-Ject® syringes are individually sterile strip packed and are also available in bulk, unsterile packaging.

## Conventional Rubber Tip Syringes

Although in general use for medical injections, rubber tipped disposable syringes have several drawbacks compared to the Norm-Ject® all plastic syringes. Since they contain natural or synthetic rubber and silicone oil, they may contaminate any media exposed to the syringe. Adverse effects from rubber curing agents, which contain zinc, sulfur and other elements, have been observed on everything from chromatographs to in-vitro procedures. Synthetic rubber plungers are also much more reactive to organic solvents and the silicone oil can cause separation of emulsions. In short, your best choice is the Norm-Ject® syringe if you are interested in purity and accuracy.

**Air-Tite**  
—Products Co., Inc.

Syringe and Needle Specialists Since 1926

### TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ [atinfo@air-tite.com](mailto:atinfo@air-tite.com)

🌐 [www.air-tite.com](http://www.air-tite.com)

# NORM-JECT<sup>®</sup>

## Syringes



### Sterile Syringes

Description	Air-Tite Part#	Syringes per Pac	Packs per Case
1ml Tuberculin <sup>1</sup>	A1	100	18
3ml Luer Slip cen	A3	100	25
5ml Luer Slip cen <sup>2</sup>	A5	100	20
10ml Luer Slip ecc <sup>2</sup>	A10	100	12
20ml Luer Slip ecc <sup>2</sup>	A20	100	8
30ml Luer Slip ecc	A30	50	10
50ml Luer Slip ecc <sup>2</sup>	A50	30	100
50ml Catheter Tip <sup>2</sup>	AC50	30	10
3ml Luer Lock	AL3	100	25
5ml Luer Lock <sup>2</sup>	AL5	100	20
10ml Luer Lock <sup>2</sup>	AL10	100	12
20ml Luer Lock <sup>2</sup>	AL20	100	8
30ml Luer Lock	AL30	50	10
50ml Luer Lock <sup>2</sup>	AL50	30	10

### Non-Sterile Syringes

Description	Air-Tite Part#	Syringes per Pac	Packs per Case
1ml Luer Slip mini bulk	ABC1	100	70
3ml Luer Slip Mini bulk	ABC3LS	100	60
5ml Luer Slip cen <sup>2</sup> mni bulk	ABC5LS	100	32
10ml Luer Slip cen <sup>2</sup> mni bulk	ABC10LS	100	19
20ml Luer Slip cen <sup>2</sup> mni bulk	ABC20LS	100	8
30ml Luer Slip ecc mni bulk	ABC30LS	50	10
50ml Luer Slip ecc <sup>2</sup> mni bulk	ABC50LS	30	16
50ml Catheter Tip <sup>2</sup> min bulk	ABC50C	30	10
3ml Luer Lock mini bulk	ABC3LL	100	60
5ml Luer Lock <sup>2</sup> mini bulk	ABC5LL	100	32
10ml Luer Lock <sup>2</sup> mini bulk	ABC10LL	100	19
20ml Luer Slip <sup>2</sup> mini bulk	ABC20LL	100	10
30ml Luer Slip mini bulk	ABC30LL	50	16
50ml Luer Slip <sup>2</sup> mini bulk	ABC50LL	30	16

<sup>1</sup> Dose saver design with low dead space plug on the piston to minimize waste - recommended for Botox<sup>®</sup>

<sup>2</sup> 5 ml syringes are graduated to 6ml, 10ml syringes are graduated to 12ml, 20ml syringes are graduated to 24ml, and 50ml syringes are graduated to 60ml.

**Air-Tite**  
—Products Co., Inc.

Syringe and Needle Specialists Since 1926

#### TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com