

B | BRAUN OMNIFIX®

- > German Engineered
- > Fine 0.01 ml graduation for precise dosage of the smallest volumes
- > Highly Transparent barrel with high contrast black graduation for easy reading
- > Rubber Piston with double sealing ring allows slow aspiration and injection
- > LDS plunger allows for minimal residual volume
- > PVC, latex, and DEHP free
- > Individual sterile strip packs
- > FDA registered. ISO 9001 Certified.



Air-Tite
— Products Co., Inc.

Syringe and Needle Specialists Since 1926

TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com

OMNIFIX[®] Low Dead Space Syringes

1 ML (LOW DEAD SPACE)

Part #: BBL1

100 Pieces per Box

B. Braun Omnifix syringes are high quality German engineered disposable rubber piston luer lock syringes. These syringes have a highly transparent polypropylene barrel, high contrast black graduations for easy reading, polystyrene plungers with an integrated backstop to help prevent unintentional plunger withdrawal, and a rubber piston with a double sealing ring made of polyisoprene allowing for slow aspiration and injection.

All OMNIFIX syringes are packaged in individual sterile strip packs. FDA registered. ISO 9001 certified. Manufacturer #: 9167006V-02
100 Pieces per Box



TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com

Air-Tite

—Products Co., Inc.

Syringe and Needle Specialists Since 1926