

# B | BRAUN OMNIFIX®

- > German Engineered
- > Fine 0.01 ml graduation for precise dosage of the smallest volumes
- > Highly Transparent barrel with high contrast black graduation for easy reading
- > Rubber Piston with double sealing ring allows slow aspiration and injection
- > LDS plunger allows for minimal residual volume
- > PVC, latex, and DEHP free
- > Individual sterile strip packs
- > FDA registered. ISO 9001 Certified.



**Air-Tite**  
— Products Co., Inc.

Syringe and Needle Specialists Since 1926

**TO ORDER:**

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ [atinfo@air-tite.com](mailto:atinfo@air-tite.com)

🌐 [www.air-tite.com](http://www.air-tite.com)

# OMNIFIX®

## Low Dead Space Syringes

1 ML ( LOW DEAD SPACE )

Part # : BBL1

100 Pieces per Box

B. Braun Omnifix syringes are high quality German engineered disposable rubber piston luer lock syringes. These syringes have a highly transparent polypropylene barrel, high contrast black graduations for easy reading, polystyrene plungers with an integrated backstop to prevent unintentional plunger withdrawal, and a rubber piston with a double sealing ring made of polyisoprene allowing for slow aspiration and injection.

All OMNIFIX syringes are packaged in individual sterile strip packs. FDA registered. ISO 9001 certified. Manufacturer # : 9167006V-02  
100 Pieces per Box



**TO ORDER:**

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ [atinfo@air-tite.com](mailto:atinfo@air-tite.com)

🌐 [www.air-tite.com](http://www.air-tite.com)

**Air-Tite**  
Products Co., Inc.

Syringe and Needle Specialists Since 1926