


SOFT-JECT®

HSW SOFT-JECT® – 3-part disposable syringes HSW SOFT-JECT® 1 ml – 50mL in Luer Slip, Luer Lock and Catheter version

- Premium syringe
- Smooth flow, tight, high transparency barrel
- Latex-free
- Pyrogen-free, PVC-free, non-toxic
- Double seal for increased safety
- Lubricated synthetic rubber piston allows fine dosing and smooth draw
- Available in Luer, Luer Lock and Catheter
- 10 ml syringes with 0.2 ml increments
- 50 ml and 100 ml syringes with extended graduation for wider range of application
- Also available with by-packed or mounted needle
- HSW SOFT-JECT® Luer Lock syringes are compatible with most commercially available syringe pumps



Luer Lock



Luer Slip

Air-Tite

— Products Co., Inc.

Syringe and Needle Specialists Since 1926

TO ORDER:


☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com

SOFT-JECT® Syringes



LUER SLIP SYRINGES	PART #	# PER BOX
1mL Tuberculin*	R1	100
2.5mL Centric	R2	100
5mL Centric	R5	100
10mL Eccentric**	R10	100
20mL Eccentric	R20	100
30mL Eccentric	R30	60

120 Catheter Tip	RC100	30
------------------	-------	----

LUER LOCK SYRINGES	PART #	# PER BOX
1mL Tuberculin*	RL1	100
3mL	RL3	100
5mL	RL5	100
10mL**	RL10	100
20mL	RL20	100
30mL	RL30	60
50mL**	RL50	50

*1mL has a low dead space plunger

**10mL graduated to 12mL, 50mL graduated to 60mL

LUER SLIP SYRINGES WITH BY-PACKED NEEDLES

SPECS	PART #	# PER BOX
2.5mL x 20g x 1"	R2201	100
2.5mL x 20g x 1-1/2"	R220112	100
2.5mL x 22g x 3/4"	R22234	100
2.5mL x 22g x 1"	R2221	100
2.5mL x 25g x 5/8"	R22558	100
20mL x 18g x 1-1/2"	R2018	100

TO ORDER:

Air-Tite
Products Co., Inc.

Syringe and Needle Specialists Since 1926

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com