

TSK

TSK STERJECT® Premium Needles

Specialized needles for enhanced precision, efficiency, and comfort.

- Premium Quality Needles
- Hard-to-Find Sizes, for the Specialist
- Available as Small as 33g, which are 22% Thinner than Regular 30g Needles
- Regular Wall, Thin Wall, and Ultra-Thin Wall
- FDA Approved and CE Marked
- ISO Certified
- Made in Japan

Needle Walls

Thin Wall



Regular Wall

A thin wall needle has a larger inner diameter than a corresponding regular wall needle. Outer diameters of the two needles are identical. The advantage of a thin wall needle is increased flow with flow equal to a needle one gauge size larger. For example, a 20 gauge thin wall has the same flow volume as a 19 gauge regular wall, but is thinner for more patient comfort.



Air-Tite

— Products Co., Inc.

Syringe and Needle Specialists Since 1926

TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com

TSK

GAUGE	LENGTH	CATALOG #
22g	70mm (2-3/4")	TSK2276
	100mm (4")	TSK22101
25g	50mm (2")	TSK2551
27g	6mm (1/4") cap	TSK276
	13mm (1/2") blister	TSK2713U*
30g	4mm (3/16") cap	TSK304
	6mm (1/4") cap	TSK306
	13mm (1/2") cap	TSK3013U*
	25mm (1") cap	TSK3025
31g	13mm (1/2") cap	TSK3113
	13mm (1/2") blister	TSK3113B
32g	4mm (3/16") cap	TSK324
	9mm (3/8") blister	TSK3209**
	13mm (1/2") cap	TSK3213
	13mm (1/2") blister	TSK3213B
33g	13mm (1/2") cap	TSK3313

* Ultra Thin Wall, ideal for dermal fillers

** 32g x 9mm popular for retinal procedures



Shown on HSW Syringes, available through Air-Tite Products

TSK needle in blister pack



TSK needle in capsule pack



Air-Tite

Products Co., Inc.

Syringe and Needle Specialists Since 1926

TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com