


HENKE-JECT®

3-PART SYRINGES

HSW HENKE-JECT® 3-part disposable syringes

1 mL - 30 mL in Luer Slip and Luer Lock Versions

- Premium syringe
- Smooth flow, tight, high transparency barrel
- Latex-free
- Pyrogen-free, PVC-free, non-toxic
- Double seal for increased safety
- Lubricated synthetic rubber piston allows fine dosing and smooth draw
- 10 ml syringes with 0.2 ml increments
- Also available with by-packed or mounted needle
- HSW HENKE-JECT® Luer Lock syringes are compatible with most commercially available syringe pumps



Air-Tite
Products Co., Inc.

Syringe and Needle Specialists Since 1926

TO ORDER:


☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com

HENKE-JECT® 3-Part Syringes



LUER SLIP SYRINGES	PART #	# PER BOX
1mL Tuberculin*	R1	100
2.5mL Centric	R2	100
5mL Centric	R5	100
10mL Eccentric**	R10	100
20mL Eccentric	R20	100
30mL Eccentric	R30HJ	100

LUER LOCK SYRINGES	PART #	# PER BOX
1mL Tuberculin*	RL1	100
2.5mL	RL2	100
5mL	RL5	100
10mL**	RL10	100
20mL	RL20	100
30mL	RL30HJ	100

*1mL has a low dead space plunger

**10mL graduated to 12mL

LUER SLIP SYRINGES WITH BY-PACKED NEEDLES

SPECS	PART #	# PER BOX
2.5mL x 22g x 3/4"	R22234	100
2.5mL x 22g x 1"	R2221	100
2.5mL x 25g x 5/8"	R22558	100
20mL x 18g x 1-1/2"	R2018	100

TO ORDER:

Air-Tite
Products Co., Inc.
Syringe and Needle Specialists Since 1926

☎ 1 (800) 231-7762
☎ 1 (757) 340-2912

✉ atinfo@air-tite.com
🌐 www.air-tite.com