

# 2-PART SYRINGES



## 2-Part All Plastic Syringes

Norm-Ject® and Air-Tite Premium 2-Part syringes are latex free, contain no rubber silicone oil, styrene, BPA, or DEHP. These syringes are the choice for any situation needing an inert, non-reactive syringe. Because of their purity, these syringes have found a wide range of specialty applications including IVF, embryo-transfer, chromatography, and many laboratory procedures. Manufactured only from laboratory grade polypropylene and polyethylene, these unique plastic syringes have a positive safety stop to help prevent accidental spills. These 2-part syringes are individually sterile strip packed and are also available in bulk, unsterile packaging.



## Conventional Rubber Tip Syringes

Although in general use for medical injections, rubber tipped disposable syringes have several drawbacks compared to the Norm-Ject® and Air-Tite Premium all-plastic syringes. Since they contain natural or synthetic rubber and silicone oil, they may contaminate any media exposed to the syringe. Adverse effects from rubber curing agents, which contain zinc, sulfur and other elements, have been observed on everything from chromatographs to in-vitro procedures. Synthetic rubber plungers are also much more reactive to organic solvents and the silicone oil can cause separation of emulsions. In short, your best choices are the Norm-Ject® and Air-Tite Premium syringes if you are interested in purity and accuracy.



**Air-Tite**  
—Products Co., Inc.  
Syringe and Needle Specialists Since 1926

### TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ [atinfo@air-tite.com](mailto:atinfo@air-tite.com)

🌐 [www.air-tite.com](http://www.air-tite.com)

# 2-PART SYRINGES



## Sterile Syringes

Description	Air-Tite Part #	Syringes per Box	Boxes per Case
1mL Luer Slip Tuberculin <sup>1</sup>	A1	100	18
3mL Luer Slip (cen)	A3	100	25
5mL Luer Slip (cen) <sup>2</sup>	A5	100	20
10mL Luer Slip (ecc) <sup>2</sup>	A10	100	12
20mL Luer Slip (ecc) <sup>2</sup>	A20	100	8
50mL Luer Slip (ecc) <sup>2</sup>	AMS50*	30	12
1mL Luer Lock <sup>1</sup>	AML1*	100	20
3mL Luer Lock	AL3	100	25
5mL Luer Lock <sup>2</sup>	AL5	100	20
10mL Luer Lock <sup>2</sup>	AL10	100	12
20mL Luer Lock <sup>2</sup>	AL20	100	8
50mL Luer Lock <sup>2</sup>	AML50*	30	12


## Non-Sterile Syringes

Description	Air-Tite Part #	Syringes per Pack	Packs per Case
1mL Luer Slip Mini Bulk <sup>1</sup>	ABC1	100	70
3mL Luer Slip (cen) Mini Bulk	ABC3LS	100	60
5mL Luer Slip (cen) Mini Bulk <sup>2</sup>	ABC5LS	100	32
10mL Luer Slip (cen) Mini Bulk <sup>2</sup>	ABC10LS	100	19
20mL Luer Slip (cen) Mini Bulk <sup>2</sup>	ABC20LS	100	10
50mL Luer Slip (ecc) Mini Bulk <sup>2</sup>	AMB50LS*	30	12
1mL Luer Lock Mini Bulk <sup>1</sup>	AMB1LL*	100	42
3mL Luer Lock Mini Bulk	ABC3LL	100	60
5mL Luer Lock Mini Bulk <sup>2</sup>	ABC5LL	100	32
10mL Luer Lock Mini Bulk <sup>2</sup>	ABC10LL	100	19
20mL Luer Lock Mini Bulk <sup>2</sup>	ABC20LL	100	10
50mL Luer Lock Mini Bulk <sup>2</sup>	AMB50LL*	30	12

\* Air-Tite Premium Lab/Vet 2-Part brand. Not for human use.

<sup>1</sup> Dose saver design with low dead space plunger to minimize waste.

<sup>2</sup> 5mL syringes are graduated to 6mL, 10mL syringes are graduated to 12mL, 20mL syringes are graduated to 24mL, and 50mL syringes are graduated to 60mL.



**Air-Tite**  
Products Co., Inc.

Syringe and Needle Specialists Since 1926

### TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com