

# TSK



## TSK STERJECT® Premium Needles

Specialized needles for enhanced precision, efficiency, and comfort.

- Premium Quality Needles
- Hard-to-Find Sizes, for the Specialist
- Available as Small as 33g, which are 22% Thinner than Regular 30g Needles
- Regular Wall, Thin Wall, and Ultra-Thin Wall
- FDA Approved and CE Marked
- ISO Certified
- Made in Japan



### Needle Walls

Thin Wall      Ultra-Thin Wall



Regular Wall

A thin wall needle has a larger inner diameter than a corresponding regular wall needle. Outer diameters of the two needles are identical. The advantage of a thin wall needle is increased flow with flow equal to a needle one gauge size larger. For example, a 20 gauge thin wall has the same flow volume as a 19 gauge regular wall, but is thinner for more patient comfort.

# Air-Tite

—Products Co., Inc.

Syringe and Needle Specialists Since 1926

**TO ORDER:**

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ [atinfo@air-tite.com](mailto:atinfo@air-tite.com)

🌐 [www.air-tite.com](http://www.air-tite.com)



# STERiJECT™ Needles

See Low Dead Space Needles for additional TSK Needle sizes



GAUGE	LENGTH	CATALOG #
25g	13mm (1/2")	TSK2513U*
27g	13mm (1/2") 19mm (3/4") 38mm (1-1/2")	TSK2713U* TSK2719U* TSK2738U*
30g	4mm (3/16") 6mm (1/4") 13mm (1/2")	TSK304U* TSK306U* TSK3013U*
31g	13mm (1/2")	TSK3113
32g	4mm (3/16") 6mm (1/4") 9mm (3/8") 13mm (1/2")	TSK324 TSK3206 TSK3209C TSK3213
33g	4mm (3/16") 13mm (1/2")	TSK3304 TSK3313

\*Ultra-Thin Wall, ideal for dermal fillers


\*\*32g x 9mm popular for retinal procedures



Shown on HSW Syringes,  
available through Air-Tite Products



TSK Needle in Capsule Pack



**Air-Tite**  
Products Co., Inc.  
Syringe and Needle Specialists Since 1926

TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com