

Since 1926


Air-Tite Products Co., Inc.

AIR-TITE BRAND NEEDLES

FOR LABORATORY AND VETERINARY USE

AIR-TITE PRODUCTS OFFERS THE BROADEST RANGE OF LAB AND VET USE DISPOSABLE HYPODERMIC NEEDLES.

- From ¼" to 5"
- 14g to 27g
- Short bevels on some needles for repetitive dosing durability and sharpness
- Specialty needles ideally designed for unique purposes
- #304 Surgical Stainless Steel Cannula
- Polypropylene Hubs
- Sterile blister packs
- Custom projects welcome



Air-Tite

Products Co., Inc.

Syringe and Needle Specialists Since 1926



STANDARD SIZE NEEDLES:

	1/2"	3/4"	1"	1-1/2"
14g				N14112
16g		N1634	N161	N16112
18g	N1812	N1834	N181	N18112
19g			N191	N19112
20g	N2012	N2034	N201	N20112
21g			N211	N21112
22g		N2234	N221	N22112
23g		N2334	N231	
25g			N251	N25112
26g	N2612			
27g	N2712			

TO ORDER:

☎ 1 (800) 231-7762

☎ 1 (757) 340-2912

✉ atinfo@air-tite.com

🌐 www.air-tite.com

SHORT BEVEL NEEDLES

ideal for repetitive dosing:

Catalog	Gauge	Length
N1814	18g	1/4"
N2014	20g	1/4"
N2038	20g	3/8"
N2238	22g	3/8"

LONG NEEDLES

	2"	3"	3-1/4"	3-1/2"	4"	5"
14g	N142		N14314			
16g		N163			N164	N165
18g	N182	N183				
19g		N193				
22g		N223		N22312	N224	
27g	N272					

BLOOD COLLECTION NEEDLES

	1"	1-1/2"
16g Single-Type	NBC161	NBC16112
18g Single-Type	NBC181	NBC18112
18g Multi-Type	NBC181M	NBC18112M
20g Multi-Type	NBC201	NBC20112

**WE ALSO CARRY HSW ECO METAL
HUB FROM 18G – 22G
12 PER PACK, NON-STERILE
ACCESS TO ALL HSW SIZES UPON REQUEST**



IDEAL FOR SPECIALTY PURPOSES

Large Animal Blood Draw:	Sonography:
16g x 3"	22g x 3"
16g x 4"	22g x 3-1/2"
16g x 5"	22g x 4"
18g x 3"	

Companion Animal

Vaccinations:

N2558 – 25g x 5/8"

N2234 – 22g x 3/4"

Air-Tite

Products Co., Inc.

Syringe and Needle Specialists Since 1926

TO ORDER:

☎ **1 (800) 231-7762**

☎ **1 (757) 340-2912**

✉ atinfo@air-tite.com

🌐 www.air-tite.com