



Smooth Entry • Smooth Delivery



Air-Tite Product's newest addition, the EZ FLO line of hypodermic needles, offer high quality and a broad selection at an affordable price. Ranging from 18 gauge to 30 gauge, these needles offer a standard bevel, thin wall from 18g to 25g and regular-wall in 27g and 30g to combine comfort and durability. Each needle is sterile, individually packaged in a box of 100. Made of surgical stainless steel and a polypropylene hub.

Since 1926

Air-Tite
Products Co., Inc.

TO ORDER:

1 (800) 231-7762

atinfo@air-tite.com

www.air-tite.com

EZ FLO By AIR-TITE

Air-Tite
 Products Co., Inc.
 Syringe and Needle Specialists Since 1926

Catalog Number	Gauge	Length	Sterility	Quantity per Box	Quantity per Case
NEZ18034	18	3/4"	Sterile	100	2000
NEZ18100	18	1"	Sterile	100	2000
NEZ18112	18	1-1/2"	Sterile	100	2000
NEZ18200	18	2"	Sterile	100	2000
NEZ19100	19	1"	Sterile	100	2000
NEZ19112	19	1-1/2"	Sterile	100	2000
NEZ20100	20	1"	Sterile	100	2000
NEZ20112	20	1-1/2"	Sterile	100	2000
NEZ21100	21	1"	Sterile	100	2000
NEZ21112	21	1-1/2"	Sterile	100	2000
NEZ21200	21	2"	Sterile	100	2000
NEZ22034	22	3/4"	Sterile	100	2000
NEZ22100	22	1"	Sterile	100	2000
NEZ22112	22	1-1/2"	Sterile	100	2000
NEZ22200	22	2"	Sterile	100	2000
NEZ22300	22	3"	Sterile	100	2000
NEZ22400	22	4"	Sterile	100	2000
NEZ23034	23	3/4"	Sterile	100	2000
NEZ23100	23	1"	Sterile	100	2000
NEZ23112	23	1-1/2"	Sterile	100	2000
NEZ25058	25	5/8"	Sterile	100	2000
NEZ25100	25	1"	Sterile	100	2000
NEZ25112	25	1-1/2"	Sterile	100	2000
NEZ25200	25	2"	Sterile	100	2000
NEZ27100	27	1"	Sterile	100	2000
NEZ27112	27	1-1/2"	Sterile	100	2000
NEZ27114	27	1-1/4"	Sterile	100	2000
NEZ27200	27	2"	Sterile	100	2000
NEZ30012	30	1/2"	Sterile	100	2000
NEZ30100	30	1"	Sterile	100	2000

